

# Vivo 3 Clinician's Manual



# Table of Contents

1	Introduction.....	9
1.1	Manufacturer Information.....	10
1.2	What is the Vivo 3? .....	10
1.2.1	Non Invasive Interfaces .....	11
1.2.2	Mobility and Usage Environment .....	11
1.2.3	Continuous Operation.....	11
1.2.4	Service Life .....	11
1.3	Intended Use.....	11
1.4	Intended Users.....	11
1.4.1	Respiratory health care specialists.....	11
1.4.2	Lay Operators .....	12
1.4.3	Service Personnel.....	12
1.5	Contraindications .....	13
1.6	About this Manual .....	14
1.6.1	Audience .....	14
1.6.2	Icons in this manual .....	14
2	Safety Information .....	16
2.1	General Use — Warnings and Precautions.....	16
2.2	Electricity — Warnings and Precautions.....	19
2.2.1	Electromagnetic Compatibility and Electrostatic Discharge (EMC and ESD).....	20
2.3	Environment — Warnings and Precautions.....	20
2.3.1	Disposal .....	21
2.4	Patient Circuit — Warnings and Precautions.....	22
2.5	Filter Usage — Warnings and Precautions.....	24
2.6	Invasive Use — Warnings and Precautions .....	25
2.7	Humidification and Heating — Warnings and Precautions.....	26
2.8	Cleaning and Maintenance — Warning and Precautions.....	27
2.9	Oxygen Usage — Warning and Precautions .....	28
2.10	Mobile Use — Warning and Precautions .....	30
3	Product Description.....	32
3.1	Main Components.....	32
3.2	Ventilator Front .....	33
3.3	Ventilator Back .....	35
3.4	Power Management.....	36
3.5	Menus.....	37
3.5.1	Use the Menu.....	37
3.5.2	The Monitor Page .....	40
3.5.3	The Setup Pages .....	41
3.5.4	The Alarm Pages .....	42
3.5.5	The Others Pages.....	42

3.6	Symbols on the Vivo 3 .....	46
3.7	Symbols on the Display.....	50
4	Treatment Functions and Settings .....	51
4.1	Treatment Modes .....	51
4.1.1	PCV+A — Pressure Controlled Ventilation (Assisted) .....	51
4.1.2	PSV — Pressure Support Ventilation .....	51
4.1.3	S — Spontaneous.....	52
4.1.4	S/T — Spontaneous/Timed .....	52
4.1.5	T — Timed.....	52
4.1.6	CPAP — Continuous Positive Airway Pressure .....	53
4.2	Treatment Settings.....	53
4.2.1	Ventilation mode .....	53
4.2.2	Inspiratory Pressure (IPAP) .....	53
4.2.3	Expiratory Positive Airway Pressure (EPAP).....	54
4.2.4	Continuous Positive Airway Pressure (CPAP).....	54
4.2.5	Breath Rate .....	55
4.2.6	Inspiratory Time (Insp. Time) .....	55
4.2.7	Minimum Inspiratory Time (Min Insp. Time) .....	55
4.2.8	Maximum Inspiratory Time (Max. Insp. Time) .....	56
4.2.9	Inspiratory Trigger (Insp.Trigger) .....	56
4.2.10	Expiratory Trigger (Exp.Trigger) .....	57
4.2.11	Rise Time .....	57
4.2.12	Target Volume .....	58
4.2.13	Max Pressure .....	58
4.2.14	Min Pressure .....	59
4.2.15	Ramp Up .....	59
4.2.16	Ramp Down .....	60
4.2.17	Ramp Pressure .....	60
4.2.18	Humidifier .....	61
4.2.19	Humidifier Level .....	61
4.2.20	Heated Circuit .....	62
4.2.21	Heated Circuit Level.....	62
5	Prepare the Vivo 3 for Use .....	63
5.1	Checking the Vivo 3 before First Use .....	63
5.2	Placing the Vivo 3.....	63
5.3	Connecting the Vivo 3 to Power Supply.....	64
5.4	Connecting the Patient Circuit .....	65
5.4.1	Connect the Heated Patient Circuit .....	66
5.4.2	Disconnect the Heated Patient Circuit from the Ventilator .....	66
5.5	Power up the Vivo 3 .....	66
5.6	Adjusting the Settings .....	67
6	How to Use the Vivo 3 .....	69

6.1	Switch On the Vivo 3 .....	69
6.2	Switch Between Clinical Mode and Home Mode .....	69
6.3	Start the Treatment .....	69
6.4	Stop the Treatment .....	70
6.5	Switch Off the Vivo 3 .....	71
6.6	Monitor Treatment .....	71
	6.6.1 Treatment Values Monitored by the Vivo 3 .....	72
	6.6.2 Treatment Values Monitored by External Equipment .....	73
6.7	Using Accessories .....	73
	6.7.1 Using the Attachable Humidifier .....	74
	6.7.2 Using the Heated Circuit .....	79
	6.7.3 Using the SpO <sub>2</sub> Sensor .....	80
	6.7.4 Using the Oxygen Supply Adapter .....	82
	6.7.5 Using the Nurse Call Connection .....	83
	6.7.6 Using the Remote Alarm Unit .....	84
	6.7.7 Using the Protective Cover .....	85
	6.7.8 Using the Lightweight Mobility Bag .....	85
	6.7.9 Using the Y-Cable .....	86
	6.7.10 Using the Vivo 3 with the Trolley .....	86
6.8	Transfer Data between the Vivo 3 and a PC .....	87
	6.8.1 Log Data .....	88
	6.8.2 Transfer Data using a Memory Card .....	88
	6.8.3 Transfer Data using an USB Cable .....	89
6.9	Basic Troubleshooting .....	89
7	Alarms .....	90
	7.1 Operator's Position .....	91
	7.1.1 Checking the Operator's Position .....	91
	7.2 Handle Alarms .....	91
	7.2.1 Identify an Alarm Condition .....	91
	7.2.2 Pause the Alarm Sound .....	93
	7.2.3 Reset an Alarm .....	93
	7.2.4 View Historical Alarms .....	93
	7.2.5 Adjust the Alarm Sound Level .....	94
	7.3 Physiological Alarms .....	95
	7.3.1 High Pressure Alarm .....	95
	7.3.2 Low Pressure Alarm .....	96
	7.3.3 High Breath Rate Alarm .....	97
	7.3.4 Low Breath Rate Alarm .....	97
	7.3.5 High Minute Volume Alarm .....	98
	7.3.6 Low Minute Volume Alarm .....	98
	7.3.7 Rebreathing Alarm .....	99

	7.3.8	Apnea Alarm .....	99
	7.3.9	Disconnection Alarm .....	100
	7.3.10	High EPAP Alarm .....	101
	7.3.11	Low EPAP Alarm .....	101
	7.3.12	High SpO <sub>2</sub> Alarm.....	102
	7.3.13	Low SpO <sub>2</sub> Alarm.....	102
	7.3.14	High Pulse Rate Alarm.....	103
	7.3.15	Low Pulse Rate Alarm .....	103
7.4		Technical Alarms .....	104
	7.4.1	High Pressure Limitation Alarm .....	104
	7.4.2	Power Fail Alarm.....	104
	7.4.3	SpO <sub>2</sub> Sensor Failure / Disconnected Alarm .....	104
	7.4.4	SpO <sub>2</sub> Artifact .....	105
	7.4.5	Ambient Pressure Compensation Lost Alarm .....	106
	7.4.6	Low Battery Alarm.....	106
	7.4.7	Critically Low Battery Alarm .....	106
	7.4.8	High Patient Air Temp. (High Patient Air Temperature) .....	107
	7.4.9	Flow Sensor Failure .....	107
	7.4.10	Internal Function Failure .....	108
8		Cleaning and Maintenance .....	109
	8.1	Cleaning the Vivo 3 .....	109
		8.1.1 Clean the Main Unit Externally .....	109
		8.1.2 Air Pathway Disinfection .....	110
		8.1.3 Clean the Patient Circuit .....	110
	8.2	Clean and Replace Patient Air Inlet Filters .....	110
	8.3	Change of Patients.....	112
	8.4	Alarm Tests .....	112
		8.4.1 Mandatory Test .....	112
		8.4.2 Optional Tests .....	113
	8.5	Repair.....	115
9		Technical Specifications .....	116
	9.1	Ventilator Size and Weight.....	116
	9.2	Power Supply .....	116
	9.3	Environmental Conditions.....	117
	9.4	Pneumatic Diagram .....	119
	9.5	Technical Data.....	120
	9.6	Data Accuracy.....	122
	9.7	Emission and Immunity Declaration .....	123
		9.7.1 Guidance and Manufacturer's Declaration – Electromagnetic Immunity.....	124
		9.7.2 Guidance and Manufacturer's Declaration – Electromagnetic Emission.....	126

9.7.3	Recommended separation distances between portable and mobile RF communications equipment and the Vivo 3 .....	127
9.7.4	Recommended separation distances between external power conductors and the Vivo 3 .....	128
9.8	Delivery Settings .....	128
10	Accessories .....	130
10.1	Accessories/Consumables .....	130
10.2	Accessories .....	132
<b>Appendices</b> .....		137
A	Patient Settings Record .....	138
B	FAA Compliance.....	139





# 1 Introduction

## WARNING!



### Risk of Personal Injury

The Vivo 3 must only be used:

- For the intended treatment in accordance with this operating manual and with the instructions given by the responsible clinical personnel.
- In accordance with the operating conditions specified in this operating manual.
- In original and unmodified shape and only with accessories specified or approved by Breas Medical.

Every other use may lead to risk of personal injury!

## CAUTION!



Read this Clinician's manual thoroughly so that you completely understand how the Vivo 3 is operated and maintained before taking it into use, to ensure correct usage, maximum performance and serviceability.

Non-professional caregivers (e.g. family members) should consult the medical equipment provider's respiratory therapist if they have any questions about the function, proper use, operation, service or maintenance of the Vivo 3.

## WARNING!



The Vivo 3 is not designed for life support treatment:

- The Vivo 3 should not be used for life support treatment.
- The Vivo 3 shall only be used by patients with spontaneous breathing.
- The Vivo 3 should not be used for ventilator dependent patients.

## 1.1 Manufacturer Information



### Address

Breas Medical AB  
Företagsvägen 1  
SE-435 33 Mölnlycke Sweden

### Web Address

[www.breas.com](http://www.breas.com)

### E-mail address

[breas@breas.com](mailto:breas@breas.com)

### Phone

+46 31 868800, Order: +46 31 868820, Technical support: +46 31 868860

### Fax

+46 31 868810

## Local Representatives Information

[www.breas.com/contact-us/](http://www.breas.com/contact-us/)

## 1.2 What is the Vivo 3?

The Vivo 3 is an internally powered pressure ventilator capable of delivering invasive or non-invasive ventilatory support for the care of individuals who require long-term support from mechanical ventilation for respiratory insufficiency or respiratory failure, with or without obstructive sleep apnea.

The Vivo 3 can be operated in the following modes:

- PCV+A (Assisted Pressure Controlled Ventilation). See page 51.
- PSV (Pressure Support Ventilation). See page 51.
- S (Spontaneous). See page 52.
- S/T (Spontaneous/Timed). See page 52.
- T (Timed). See page 52.
- CPAP (Continuous Positive Airway Pressure). See page 53.

The following modes can be combined with the Target Volume setting:

- PCV+A
- PSV
- S/T
- T

### **1.2.1 Non Invasive Interfaces**

The ventilator system can be used non-invasively with nasal mask, full/total face mask, and nasal pillow interfaces.

### **1.2.2 Mobility and Usage Environment**

The ventilator system is classified as transit-operable and is intended to be used in homes, public spaces, institutions and hospitals.

The ventilator system is intended to be used together with portable applications such as wheelchairs, personal family vehicles, ground ambulances and civil aircraft (not helicopter). It is not intended for use during emergency transports.

### **1.2.3 Continuous Operation**

The ventilator can be used for continuous operation up to 24 hours/day at least for 90 days without restarting.

### **1.2.4 Service Life**

The expected service life of the Vivo 3 is 5 years or 20,000 hours.

## **1.3 Intended Use**

Vivo 3 is intended to provide non-invasive or invasive ventilation for patients weighing over 10 kg (22 lbs) who require long-term support or mechanical ventilation for respiratory insufficiency or respiratory failure, with or without obstructive sleep apnea.

Vivo 3 is intended for spontaneously breathing patients.

## **1.4 Intended Users**

This section describes the intended users of the Vivo 3, their qualifications and their related documents.

### **1.4.1 Respiratory health care specialists**

Health care professionals such as physicians and respiratory therapists, assigned to form the clinical authority when it comes to operating mechanical ventilators. They have a good understanding of the human respiratory system and a general understanding of mechanical ventilators.

They are allowed to change the clinical settings of a ventilator and prescribe new settings. They may also operate software applications for follow-up on patient's ventilator treatment.

## **Training**

The respiratory health care specialists shall be trained to a good knowledge of the Vivo 3, its capabilities and the settings that can be made. This training consists of reading the Clinician's manual in complete and it shall be conducted before operating the Vivo 3.

## **Related Documents**

The Clinician's manual is intended for the respiratory health care specialists. It shall be available for training of new personnel and as reference when operating the Vivo 3.

### **1.4.2 Lay Operators**

Day-to-day caregivers, patients, relatives and other non-professional users that operate the Vivo 3 within the prescribed settings.

They are allowed to operate the Vivo 3 with the Home mode activated. The lay operator may also perform basic maintenance that doesn't require special equipment or a service environment.

In Home mode, the device is locked in order to limit settings accessibility and hide features/controls.

The User Manual contains the information intended for patients and lay operators.

## **Training**

The lay operator shall be trained to basic knowledge of the Vivo 3 and in the specific operations they are assigned to perform. The training shall be based on the user manual and the responsible clinical personnel shall assess the level of training required for each lay operator.

## **Related Documents**

The User manual is intended for lay operator. It shall be available for the training and as reference when operating the Vivo 3.

### **1.4.3 Service Personnel**

Certified service personnel with responsibility to maintain the equipment in proper working order. They have a technical education and/or relevant experience of technical work on electrical equipment. If local or national regulations requests additional authorization or competence, these shall be complied to.

Certified service personnel may perform any repairs, upgrades or service operations that they have been certified to perform, as long as they have the required equipment and the operation is performed in an appropriate environment. They may also operate software applications for follow-up on ventilators usage and for troubleshooting.

## Training and Certification

Service personnel shall be trained on the Vivo 3 and certified by Breas for being allowed to perform any service, repairs or other operations on the Vivo 3. The training consists of reading the services manual in complete. After completed training the certification test may be performed.

## Related Documents

- The Service Manual.
- The User Manual.
- Service bulletins, available for certified service personnel on the Breas extranet.

## 1.5 Contraindications

The Vivo 3 is not a life-support ventilator and is contraindicated in patients who are unable to tolerate more than brief interruptions in ventilation.

If a patient has any of the following conditions, therapy with positive airway pressure may be contraindicated and the prescribing clinician shall decide if the benefit of ventilatory assistance outweighs the risks:

- Untreated pneumothorax
- Pneumomediastinum
- Inability to maintain a patent airway or adequately clear excessive respiratory secretions
- Severe acute systemic complications (shock, unstable arrhythmias, myocardial ischemia)
- Severe bullous lung disease
- Risk of vomiting
- Pathologically low blood pressure, especially if associated with intravascular volume depletion
- Cerebrospinal fluid leak, recent cranial surgery or trauma

The use of the Vivo 3 is contraindicated in an MRI environment.

## Adverse Effects

Patients should report any unusual chest pain, severe headache or increased breathlessness.

The following side effects may occur during the course of therapy with the Vivo 3, patients are advised to report any new or changing adverse effects to their physician:

- Nasal, mouth or throat dryness
- Nosebleeds
- Abdominal bloating
- Ear or sinus discomfort
- Eye irritation
- Skin rashes

## 1.6 About this Manual

### 1.6.1 Audience

This manual is primarily intended for care providers, clinical personnel, physicians and others who require a working knowledge of the Vivo 3 system. The manual comprises detailed information of the settings and functions of the Vivo 3, to be handled by trained health care personnel only.








- Patients and lay caregivers may read the user manual for reference purposes, after appropriate guidance from the responsible care provider.
- Service personnel may order the Vivo 3 Service Manual that contains detailed technical information for maintenance, service and repair.

### 1.6.2 Icons in this manual

In this manual, icons are used to highlight specific information. The meaning of each icon is explained in the table below.

Icon	Explanation
	<b>Warning!</b> Risk of death or personal injury.
	<b>Warning!</b> Risk of Cross-contamination.
	<b>Warning!</b> Risk of electric shock.

Icon	Explanation
	<p><b>Warning!</b> Hot surface, risk of burns.</p>
	<p><b>Warning!</b> Flammable material, risk of fire.</p>
	<p><b>Caution!</b> Risk of equipment damage, loss of data, extra work, or unexpected results.</p>
	<p><b>Note</b> Information that may be valuable but is not of critical importance, tips.</p>
	<p><b>Reference</b> Reference to other manuals with additional information on a specific topic.</p>

## 2 Safety Information

### 2.1 General Use — Warnings and Precautions

#### **WARNING!**



##### **Risk of Personal Injury**

The Vivo 3 is not designed for life support treatment:

- The Vivo 3 shall only be used by patients with spontaneous breathing.
- The Vivo 3 should not be used for ventilator dependent patients.



##### **Risk of Faulty Treatment**

If the patient is admitted to a hospital or is prescribed any other form of medical treatment, always inform the medical staff that the patient is on mechanical ventilation treatment.



##### **Risk of Insufficient Ventilation**

Usage outside the specified operating conditions may cause reduced performance. The Vivo 3 must only be used in accordance with the operating conditions specified in this operating manual



##### **Risk of Faulty Treatment**

Do not use the Vivo 3 in the event of

- Suspected damage to the device, including the occurrence of Internal Functional Failure alarms.
- Unexpected patient symptoms during treatment.
- Unexplainable or sudden changes of pressure, performance or sound during operation.
- Delivered air being abnormally hot or emitting an odor.

Contact your responsible care provider for an inspection.



##### **Risk of Suffocation**

Do not remove the protection for the heated patient circuit connector, the connector can be fitted with the protection in place.

If removed, children can choke or suffocate if swallowing it.



##### **Risk of Asphyxia or Personal Injury**

Using incorrect settings may cause personal injury or severe medical conditions, such as hypercarbia, producing arterial acidemia.

The therapy settings shall be based on a physician's description. Changes to settings must be made by authorized clinical personnel only.



**Risk of Faulty Treatment**

The responsible organization should periodically reassess the settings of the therapy for effectiveness.

**Risk of Faulty Treatment**

Always prepare the Vivo 3 as described in this manual before use.

**Risk of Unnoticed Critical Conditions**

- All the physiological alarms of the Vivo 3 must be set at safe levels that will effectively warn the user of any risk.

The alarm levels should be assessed considering the patient's treatment settings.

- Any change of treatment settings or change of components in the ventilation system may require readjustment of the alarm levels.
- The alarm sound level should be set to a clearly audible level. Setting the alarm sound level below that of the ambient sound level can impede recognition of alarm conditions.

**Risk of Unnoticed Critical Conditions**

When using the Vivo 3 invasively, physiological alarms, such as the low volume alarm and the low breath rate alarm, must be carefully set, to ensure safe use.

**Risk of Reduced Safety and Performance**

Accessories that have not been tested with the Vivo 3 might affect safety features and performance negatively. Only use the Vivo 3 with accessories approved by Breas Medical.

Incompatible parts can result in degraded performance and change of pressure gradient.

If unapproved accessories are used, Breas Medical has no responsibility for the safe and effective use of the Vivo 3.

The responsible organization is responsible for ensuring the compatibility of the ventilator and all parts used to connect to the patient before use.

**WARNING!****Risk of Electric Shock**

Modifying or using the ventilator with accessories that are not specified or approved by Breas may cause cardiac arrhythmia.

The Vivo 3 must only be used in original and unmodified shape and only with accessories specified or approved by Breas Medical.

Inadequate use of device or accessories may cause loss of treatment or decreased performance.

## **WARNING!**



### **Risk of Burns**

Covering breathing tubes with a blanket or heating them with an overhead heater can affect the quality of the therapy or injure the patient.

## **CAUTION!**



Read the Clinician's manual thoroughly and understand the operation of the Vivo 3 before operating or using the machine.



Always use the therapeutic pressure setting, as individually determined with the configuration of the equipment and accessories.



Proper placement and positioning of the patient interface is critical to the consistent operation of this equipment.



Handle the Vivo 3 with care.



Make sure to place and pack the device in a way that prevents unintentional start of the machine.



Due to the internal battery, the Vivo 3 may start if the Start/Stop button is pressed even without the mains being connected.



Do not use the Vivo 3 with nitric oxide, helium or helium mixtures. This may affect patient air flow and volume measurements.

## 2.2 Electricity — Warnings and Precautions



### WARNING!

#### Risk of Electric Shock

High voltage contact may cause cardiac arrhythmia.

- Do not operate the Vivo 3 if it has a damaged power cord, power supply or casing.
- To avoid electrical shock, only clean the Vivo 3 according to instructions in this manual. Do not soak or immerse the Vivo 3 into any fluids.
- Use the approved power supply units only.

Use of unapproved power supply units may compromise the electrical isolation and lead to risk of electric shock.

- Do not use more than one multiple portable socket-outlet or extension cord.  
If a multiple portable socket-outlet is used, it must not be placed on the floor.
- The operator must not touch accessible contacts of connectors and the patient simultaneously.
- Nurse Call must only be connected to a safety extra low voltage system with an isolation from mains voltage which complies with the requirements of IEC 60601-1.



### WARNING!

#### Risk of Faulty Treatment

Electromagnetic Interference may cause electrical equipment to malfunction.

- The aspects of electromagnetic compatibility must be considered.
  - The Vivo 3 should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, the Vivo 3 should be observed to verify normal operation in that configuration.
  - Mobile or transportable radio transmitters may interfere with the Vivo 3.
  - Further guidance for safe installation of the ventilator can be found in the chapter about emission and immunity declaration.
- If a portable AC power supply is used, make sure that the voltage variations are within the operating limits of the Vivo 3.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Vivo 3, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## 2.2.1 Electromagnetic Compatibility and Electrostatic Discharge (EMC and ESD)

Use of accessories, transducers and cables other than those specified or provided by Breas could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Electromagnetic disturbance may impair the safety and performance of the Vivo 3. The electromagnetic field levels at the Vivo 3 should not exceed 20 V/m.

Unsteady indicated values for delivered volumes or pressures and the occurrence of alarm conditions without apparent cause may be an indication of loss of performance due to electromagnetic disturbances.

Measures for keeping electromagnetic field levels low should include but are not be limited to:

- Normal precautions with regard to relative humidity and conductive characteristics of clothing in order to minimize the build-up of electrostatic charges.
- Avoiding the use of radio emitting devices (e.g. cellular or cordless telephones, microwave ovens and high-frequency surgery apparatus) closer than 1 meter to the Vivo 3.
- Avoiding the use of known sources of Electromagnetic Interference, (e.g. RFID, diathermy equipment), in the presence of the Vivo 3.

Please note some of these RF emitters may not be visible and the Vivo 3 can potentially be exposed to fields from these RF emitters without the user's awareness. If abnormal performance of the Vivo 3 is observed, and the RF emitters cannot be identified and removed, the Vivo 3 may need to be reoriented or relocated.



See the section *Emission and Immunity Declaration* for detailed information and further guidance for mitigating electromagnetic disturbance.

## 2.3 Environment — Warnings and Precautions

### WARNING!



#### Risk of Intoxication

Do not use the Vivo 3 in a toxic environment.



#### Risk of Faulty Treatment

If a room humidifier is used, place it at least 2 meters away from the Vivo 3.



### **Risk of Faulty Treatment**



Do not use or store the Vivo 3 in a magnetic resonance (MR) environment. Use of the Vivo 3 in an MR environment may result in malfunction of the Vivo 3 and pose unacceptable risk to the patient, medical staff or other persons.



### **Risk of faulty Treatment**

The performance of the Vivo 3 may deteriorate at altitudes or ambient temperatures outside the operation conditions specified in the section *Environmental Conditions*.

- Do not use the ventilator while positioned in a warm place, such as direct sunlight or close to a radiator as this might lead to temperature outside the specifications.
- Do not use the ventilator in an hyperbaric chamber, as this would cause an ambient pressure outside the specifications.
- Do not use the ventilator immediately after storage or transport outside the recommended operating conditions.

### **WARNING!**



#### **Risk of Fire**

Do not use the Vivo 3 in environments where explosive gases or flammable anesthetic agents present.

### **WARNING!**



#### **Risk of Electric Shock**

Water on and in the device may cause an electric conductive path. Do not expose the Vivo 3 to rain or snowfall.

## **2.3.1 Disposal**

The Vivo 3, any accessories and any replaced parts must be disposed of in accordance with the local environmental regulations regarding the disposal of used equipment and waste.

## 2.4 Patient Circuit — Warnings and Precautions

### WARNING!



#### **Risk of Insufficient Ventilation**

Insufficient ventilation may cause transient hypoxia.

The Vivo 3 ventilator is intended to be used with patient circuits with intentional leakage and compliant to ISO 17510. Recommended leak rate: 20 to 50 liters per minute at 10 cmH<sub>2</sub>O.

Failure to use a mask or accessory that minimizes rebreathing of carbon dioxide or permits spontaneous breathing can cause asphyxiation.



#### **Risk of Abnormal Exhalation Volume Measurement**

Non-invasive ventilation: The exhaled volume of the patient can differ from the measured exhaled volume due to unintentional leaks around the mask.

Invasive ventilation: The exhaled volume of the patient can differ from the measured exhaled volume due to unintentional leaks around the cannula.

For correct measurements, minimize unintentional leaks.



#### **Risk of Reduced Safety and Performance**

Accessories that have not been tested with the Vivo 3 might affect safety features and performance negatively. Only use the Vivo 3 with accessories approved by Breas Medical.

Incompatible parts can result in degraded performance and change of pressure gradient.

If unapproved accessories are used, Breas Medical has no responsibility for the safe and effective use of the Vivo 3.

The responsible organization is responsible for ensuring the compatibility of the ventilator and all parts used to connect to the patient before use.



#### **Risk of Reduced Performance**

Filters and patient connected parts must be replaced regularly to ensure correct function of the Vivo 3.



#### **Risk of Insufficient Ventilation**

Insufficient ventilation may cause transient hypoxia.

Before use:

- Make sure that the patient circuit and joined parts are undamaged and correctly connected, in order to avoid unwanted leakage.
- Make sure that the leakage port of the circuit or mask is not blocked or obstructed. This port prevents rebreathing by flushing the exhaled air.
- For active exhalation valve circuits: Check the function of the exhalation valve and that is not blocked or obstructed.
- The Vivo 3 should be turned on and the function of the leakage port should be checked before use: The pressurized air from the Vivo 3 causes a continuous flow of air through the leakage port, enabling flushing of exhaled air.

**Risk of Insufficient Ventilation**

Unapproved patient circuits may come loose.

To prevent disconnection of the patient circuit or patient circuit system during use, especially during ambulatory use, only patient circuits in compliance with ISO 5367 or ISO 80601-2-74.

**Risk of Suffocation**

Periodically check for moisture in the patient circuit.

When present, remove the moisture. Before attempting to dry the circuit, disconnect it from the Vivo 3 to ensure no water flows back into the Vivo 3.

The frequency at which these checks must be performed will depend on the patient's condition and the device used. The responsible caregiver should assess this on an individual basis in accordance with the patient's needs.

**Risk of Suffocation**

If the patient needs assistance to take off the patient interface, the patient shall not be left alone. This is to avoid the risk of re-breathing of CO<sub>2</sub> in case of accidental ventilator failure.

Do not breathe through the connected patient circuit unless the ventilator is turned on and operating properly.

**Risk of Suffocation**

If the patient is using a full face mask (covering mouth and nose), the mask must be equipped with a safety entrainment valve.

**Risk of Constriction**

Entanglement with cables or tubing constricting airways may cause asphyxiation.

Do not leave long lengths of air tubing or cables around the top of the bed. It could twist around the patient's head or neck while sleeping.

**Risk of Insufficient Ventilation**

Insufficient ventilation may cause transient hypoxia.

The use of equipment such as endotracheal tubes, oral/nasal tubes, adaptors etc. with small inner diameters or high resistance filters, humidifiers etc. increase the resistance in the patient circuit which may interfere with the operation of the patient disconnect function. It may also interfere with the device trigger function.

This impact can be reduced by conducting a pre-use test.

**Risk of Excessive Carbon Dioxide**

Insufficient carbon dioxide removal may cause arterial acidemia.

For reducing the risk of rebreathing CO<sub>2</sub>, make sure that the leakage port is located as near the patient interface as possible. This is extra important for treatments with low pressure, as this reduces the flow through the leakage port.

## WARNING!



### **Risk of Electric Shock**

Do not use antistatic or electrically conductive hoses or tubing with the ventilator breathing system. This could result in electrical shock.

## WARNING!



### **Risk of Cross-Contamination**

Patient circuits might get contaminated by exhaled gases. To avoid cross-contamination, always use a properly cleaned or a new patient circuit when the Vivo 3 is to be used by a new patient.

## NOTE



For masks and accessories, always follow the manufacturer's instructions.

## 2.5 Filter Usage — Warnings and Precautions

## WARNING!



### **Risk of Insufficient Ventilation**

Insufficient ventilation may cause transient hypoxia.

Replace or clean the inlet filters as specified in the *Maintenance* chapter.

Using old or clogged filters may cause the Vivo 3 to operate at higher temperatures than intended.

When operating the Vivo 3, make sure that the air inlet and filters are not obstructed or occluded.



### **Risk of Insufficient Ventilation**

Insufficient ventilation may cause transient hypoxia.

Do not use high resistance bacteria filter at the air outlet of the Vivo 3. High resistance bacteria filters placed between the air outlet and the patient interface may interfere with the operation of the patient disconnect function. It may also interfere with the device trigger function.



### **Risk of Insufficient Ventilation**

Nebulization or humidification can increase the resistance of breathing system filters and the operator must monitor the breathing system filter frequently for increased resistance and blockage to ensure the delivery of the therapeutic pressure.

Increased resistance may interfere with the operation of the patient disconnected function.



## **WARNING!**



### **Risk of Cross-Contamination**

Deep tissue or mucosal contact with infectious agents may cause infections:

- If assigning the Vivo 3 to a new patient, always replace the filters.
- Always use the Vivo 3 with patient air inlet filters installed.



### **Risk of Cross-Contamination**

Deep tissue or mucosal contact with infectious agents may cause infections.

If the Vivo 3 is used by several patients, a low resistance bacteria filter shall be used between the air outlet and the patient circuit, for preventing patient cross-contamination. Reuse of bacteria filter, patient circuit or mask may expose patients to contagious agents.

## **2.6 Invasive Use — Warnings and Precautions**

## **WARNING!**



### **Risk of Therapeutic System Performance or Failure**

Invasive use of the attachable humidifier may cause reduced performance. The attachable humidifier is not intended for invasive use.



### **Risk of Faulty Treatment**

If using the Vivo 3 for invasive ventilation, humidification is recommended. Make sure to use an external heated humidifier that complies to invasive use according to ISO 80601-2-74.



### **Risk of Unnoticed Critical Conditions**

When using the Vivo 3 invasively, physiological alarms, such as the low volume alarm and the low breath rate alarm, must be carefully set, to ensure safe use.



### **Risk of Rebreathing**

The Vivo 3 is equipped with a low leakage alarm. The low leakage alarm is not a substitute for operator vigilance in ensuring that the leakage ports remains clear at all times. Periodically check the leakage ports during therapy.



### **Risk of Rebreathing**

Be aware of risks for rebreathing exhaled gases:

- A blocked, clogged or in other ways occluded leakage port causes exhaled gases to be rebreathed.
- With low pressures the risk of rebreathing increases.
- I:E (inspiration time : expiration time) ratios close to 1:1 increase the risk of rebreathing as the expiratory time allowing exhaled gases to be cleared from the circuit before the next breath is decreased.

The Vivo 3 is equipped with a low leakage alarm. However, the alarm is not a substitute for operator vigilance for ensuring that the leakage ports remains clear at all times.

Periodically check the leakage ports during therapy.

See the section 2.4 *Patient Circuit — Warnings and Precautions*, page 22 for mitigating the risks of rebreathing.

## **2.7 Humidification and Heating — Warnings and Precautions**

### **WARNING!**



#### **Risk of Insufficient Ventilation**

The attachable humidifier is not intended for invasive use.



#### **Risk of Suffocation**

When the attachable humidifier is installed, the Vivo 3 must be located below the patient and on a flat surface. This is to prevent personal injury from accidental spillage or from excess water or condensation flowing down the patient tube and into the patient's mask. Extra cautions should be taken for patients who are unable to guard their airways or cannot pull off the mask.



#### **Risk of Suffocation or Loss of Ventilation**

Incorrect placing of the ventilator may cause transient hypoxia.

If using an external humidifier, it shall be placed below both the patient and the Vivo 3. This is to prevent personal injury from accidental spillage or from excess water or condensation flowing down to the ventilator or down the patient tube and into the patient's mask. Extra cautions should be taken for patients who are unable to guard their airways or cannot pull off the mask.



#### **Risk of Inflammation**

Incorrect connection of the ventilator may cause inflammation.

The attachable humidifier is only enabled when the Vivo 3 is connected to the AC Power supply.

**Risk of Suffocation**

Installation of a water trap may be required if the condensation is extensive in the patient circuit when using a heated humidifier.  
The water trap prevents condensed water in the patient circuit from reaching the patient airways and causing personal injury.

**Risk of Suffocation**

Do not use the attachable humidifier during mobile use.  
Due to movements, water spillage from the humidifier or condensed water may flow to the patient and cause suffocation.

**WARNING!****Risk of Electric Chock**

Do not use the attachable humidifier during mobile use.  
Internal water spillage may cause electric chocks and may damage the device.

**Risk of electric shock**

If using the protective cover or the carry bag, first remove the attachable water chamber. Water spillage may cause electric shocks.

**Risk of Electric Shock**

Electrostatic discharge or contact with leakage current may cause an electric shock.  
The attachable humidifier must be detached before filling. Do not fill above the *Maximum Water Level* indication on the water chamber.

**WARNING!****Risk of Burns**

After using the ventilator, wait one minute before you open the water chamber since it can get hot under certain conditions (for example if the humidifier runs out of water).

**CAUTION!**

The use of external humidifier may require readjustment of the low-pressure alarm.

## 2.8 Cleaning and Maintenance — Warning and Precautions

This manual contains instructions for cleaning and maintenance that can be carried out by the care provider or users with physical ability and working knowledge of the system.

## **WARNING!**



### **Risk of Faulty Treatment**

Service and Maintenance of the shall not be performed when the is in use.

## **WARNING!**



### **Risk of Electric Shock**

Cleaning with excessive water or opening the device's casing without certified training may cause electric shocks.

The Vivo 3 should be regularly cleaned and maintained in accordance with this operating manual.



### **Risk of Electric Shock**

High voltage contact may cause cardiac arrhythmia.

Repairs and modifications must be carried out by authorized technicians only and in accordance with instructions from Breas Medical

- The Vivo 3 must not be opened, repaired or modified by unauthorized personnel. If subjected to unauthorized operations, Breas Medical is no longer responsible for the device's performance and safety and all warranties will become invalid.
- The Vivo 3 must not be modified or interconnected to unapproved equipment.

## **CAUTION!**



Do not attempt to autoclave or sterilize the Vivo 3.

## **2.9 Oxygen Usage — Warning and Precautions**

When using the Vivo 3 with oxygen, always follow the oxygen provider's instructions and use only medical grade oxygen complying with local regulations.

## **WARNING!**



Do not use a humidifier between the oxygen source and the ventilator, in order to humidify the oxygen flow.

If humidification is required, use the attachable humidifier or an external humidifier after the patient air outlet.



### **Risk of faulty Treatment**

At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure delivered, the patient's breathing pattern, mask selection and leak rate.

To monitor the oxygen concentration, use an external patient air oxygen saturation monitoring unit which complies with ISO 80601-2-55 and is equipped with a high oxygen level alarm.



### **Risk of fire**

The presence of oxygen can speed up combustion of inflammable materials.



### **Risk of Fire**

When oxygen is used with the Vivo 3, the oxygen flow must be turned off when the Vivo 3 is not operating. Oxygen delivered into the patient tubing may accumulate within the machine enclosure. Oxygen accumulated in the machine enclosure increases the risk of fire.



### **Risk of Fire**

Ventilate the room adequately. Do not smoke in a room where oxygen is being used.



### **Risk of Fire**

Naked light bulbs and other sources of ignition must be kept a minimum of 2 meters (6 feet) away from the oxygen cylinder or any part of the patient circuit.



### **Risk of Fire**

Do not use aerosols or solvents close to the oxygen supply, even when the oxygen supply is shut off.



### **CAUTION!**

Supplemental oxygen flow must not exceed 30 l/min or 100 kPa.



## CAUTION!

Supplemental oxygen is added before the volume measurement sensor and thereby included in the measurements. However, the oxygen concentration still has influence on the volume measurement of the delivered air.

This measurement is based on a normal oxygen concentration of 21%. If the oxygen concentration is higher, the actual inspired volume will deviate from the monitored volume as follows:

- 40% oxygen concentration: -2.5% deviation
- 60% oxygen concentration: -5% deviation
- 80% oxygen concentration: -7.5% deviation

## 2.10 Mobile Use — Warning and Precautions

This section applies if using the Vivo 3 during transit, for example on a wheel chair or in a car.

### No Attachable Humidifier During Mobile Use



#### WARNING!

##### Risk of Suffocation

Do not use the attachable humidifier during mobile use.

Due to movements, water spillage from the humidifier or condensed water may flow to the patient and cause suffocation.



#### WARNING!

##### Risk of Electric Chock

Do not use the attachable humidifier during mobile use.

Internal water spillage may cause electric chocks and may damage the device.

### Only Approved Power Supplies



#### WARNING!

##### Risk of Electric Chock

If connecting the Vivo 3 to an external DC power source, always use the DC DC Power Supply accessory. Connecting directly to an external DC power source may compromise the electrical isolation and cause an electric shock.

If connecting the Vivo 3 to a portable AC power generator, make sure its voltage variations are within the operating limits of the Vivo 3.

## Chock and Weather Protection



### **CAUTION!**

During mobile use, protect the device with either of the Protective Cover accessory or the Lightweight Mobility Bag accessory.  
Do not use the Vivo 3 while in a carry bag.


### 3 Product Description

This section describes the main Vivo 3 medical electric equipment.

For information about accessories and user replaceable spare parts., see .10 *Accessories*, page 130.


#### 3.1 Main Components

Part	Function	Breas Part no.
Ventilator unit	Main unit.	229000




---

RRC PS90M AC/DC Power supply	Power supply adapter for the Vivo 3.	006994
------------------------------	--------------------------------------	--------






---

Mains power cord	Power cord between the mains socket and the power supply adapter.	003520
------------------	---	--------






Part	Function	Breas Part no.
<b>User Manual</b> 	Instructions for use	User's manual: 007231 Clinician's manual: 007232
<b>Carry bag</b> 	For transportation, when not in use.	007013
<b>Patient circuit, 15 mm</b> 	Delivers air to the patient	006712


### 3.2 Ventilator Front

This section describes the front panel.



No.	Item	Description
1	Screen	Displays pages with information, settings and commands.
2	Directional buttons	Moves between and selects objects on the current page.
3	On/Off button	Starts /stops treatment. Turn the ventilator On/Off.
4	Navigation buttons	Select pages according to the corresponding label on the display. The navigation buttons can temporarily be designated functions for replying to questions or requests from event windows.
5	Audio pause button	Pauses the alarm sound.
6	Audio Pause LED	Lights yellow when the alarm sound is paused.
7	Alarm LED	Flashes during active alarms.

### 3.3 Ventilator Back



No.	Item	Description
1	Power port	Port for connecting the power supply. See 5.3 <i>Connecting the Vivo 3 to Power Supply</i> , page 64
2	Oxygen port	Connection for low pressure/bleed-in oxygen
3	Patient air outlet	Connection for patient circuit
4	Patient air inlet	Patient air inlet, with filter holder
5	Air bypass module	Directs the patient air flow. Removed if the attachable humidifier is used.
6	USB port	For data transfer to a PC.
7	SD Card port	For copying records and logs to a PC.
8	Communication port	For connection to SpO <sub>2</sub> sensor or Accessory box (for connection of nurse call or remote alarm)

The ventilator can be carried by hand using the handle.

## 3.4 Power Management



### WARNING!

Do not connect the ventilator directly to the battery of a wheelchair or any similar external power source, as this can compromise the ventilator performance and result in degradation of the health of the patient.  
For DC power, always use the DC/DC power supply adapter (accessory).

The Vivo 3 has a power management system that automatically selects the best available power source, according to the priority list below:

1. Mains power via Breas Mains power supply
2. External DC via Breas DC/DC power supply adapter
3. Internal battery

The power source is indicated by a symbol at the top of the display.



Power Source	Icon
Mains	
External DC	
Internal battery	

For connecting the Vivo 3 to both mains and external DC at the same time, the accessory *Power supply Y—cable* is required.

If a power source fails, the Vivo 3 will switch to the next source in priority and show a message on the display. If all available sources fail, the *Power Failure* alarm is given and the Vivo 3 shuts down.

## Charging the Internal Battery

The internal battery is automatically charged when the Vivo 3 is connected to external power. The state of charge is indicated by a symbol on the display. The green part of the battery shows the state of charge.

Battery State of Charge	Icon
High	
Charging	

## Preparing the Battery for Storage


If storing the Vivo 3, longer than 1 month, make sure that the internal battery is about half charged before storing, in order to maintain maximum battery capacity. Optimal storage temperature is 5 to 30°C (41 to 86°F).

## 3.5 Menus

### 3.5.1 Use the Menu

The menu consists of four function sections:


- Monitor
- Setup
- Alarm
- Others




- a. Navigation buttons
- b. Active page indicator
- c. Directional buttons
- d. Page number indicator (for functions with several pages)

### Selecting the Section to Display

- 1 Press the navigation button for the requested function page.  
⇒The page is now displayed.



- 2 For functions with several pages grouped together, press the navigation button again to browse the pages.



- 3 For pages with menus, settings or additional information, press the Up and Down arrow buttons to select an item on the page.

### Select an Item on a Page


This procedure describes how to navigate between selectable items. Read-only information cannot be selected.

#### Selectable Items

- Sub pages (indicated by an arrow to the right of the item text)
- Settings
- Commands

- 1 Press the Up or Down arrow button to select the first selectable item on the page.

⇒ The selected item is highlighted.



- 2 To select another item on the list, press the Up or Down arrow until it is highlighted.


Note that items with read only information cannot be selected.

### Enter a Sub Page

- 1 Select the sub page (indicated by an arrow to the right of the item text) using the Up or Down arrow button.

- 2 Enter the sub page by pressing the right arrow button.

⇒ The sub page is now displayed.




- 3 Leave the sub page using the left arrow button

### Change a Setting

- 1 Select the setting (indicated by an arrow to the right of the item text) using the Up or Down arrow button.


- 2 Change the setting using the right or left arrow buttons to select between the predefined values.  
⇒ The currently displayed value will remain selected when leaving the page.



- 3 Leave the sub page using the navigation buttons.

### Execute a Command


- 1 To execute a command, select it and press the right arrow button.  
⇒ The command execution is started. Additional actions related to the command might be requested in event windows.



## 3.5.2 The Monitor Page

The monitor page displays the treatment data. It consists of a bar graph illustrating the current pressure and text area displaying current monitored values in text.





### Bar graph legend

1. Current pressure
2. Low pressure alarm level
3. Set inspiration pressure level
4. High pressure alarm level

For information about monitored values, see 6.6.1 *Treatment Values Monitored by the Vivo 3*, page 72.

### 3.5.3 The Setup Pages

The setup pages contain settings related to the treatment.

#### **WARNING!**



#### **Risk of Asphyxia or Personal Injury**

Using incorrect settings may cause personal injury or severe medical conditions, such as hypercarbia, producing arterial acidemia.

The therapy settings shall be based on a physician's description. Changes to settings must be made by authorized clinical personnel only.

## Dynamic Content

The content on the setup pages are dynamic: Only applicable settings are displayed. For example, if having the setting *Heated Circuit* set to *Off*, you will not see the settings for heat level setting.

The most significant setting is *Ventilation Mode*. If changing the ventilation mode, always review and confirm the following settings as new settings may have been added.



For detailed information about each setting, see 4 *Treatment Functions and Settings*, page 51.

### 3.5.4 The Alarm Pages

The alarm pages contain alarm settings and an alarm history list.

For information about the alarms and their possible settings, see 7 *Alarms*, page 90.

### 3.5.5 The Others Pages

The Others pages contain non-clinical settings and information.

#### 3.5.5.1 Device Settings

This section describes device settings that don't affect the ventilating function of the Vivo 3.

Setting	Description
Language	<b>Value range:</b> The available languages are listed with their native names accompanied by a flag representing a country where the language is spoken. <b>Default value:</b> English (The language selection menu is displayed when the Vivo 3 is started for the first time.)
Pressure Unit	<b>Value range:</b> <ul style="list-style-type: none"><li>• cmH<sub>2</sub>O</li><li>• mbar</li><li>• hPa</li></ul> <b>Default value:</b> cmH <sub>2</sub> O

Setting	Description
Display Light	<p><b>Value range:</b></p> <ul style="list-style-type: none"> <li>On — Always lit at the selected intensity</li> <li>Auto — Always lit, with automatic adjustment of the intensity with regards to the ambient light.</li> <li>Delayed — The display is dimmed after about 30 seconds (time depends on usage mode and power source). If any button is pressed or any alarm occurs, the display is lit again.</li> </ul> <p><b>Default value:</b> On</p>
Light Intensity	<p><b>Value range:</b> 1 to 5</p> <p><b>Default value:</b> 5</p>

### 3.5.5.2 Pre-use Test

The pre-use test is performed using a guide in the graphical user interface.

### 3.5.5.3 Device Information

Information	Description
Product name	
Battery Percentage	Battery charge level
Battery Time Left	Estimated remaining usage time on battery power
Firmware version	Technical information for service personnel
Bootloader version	Technical information for service personnel
Board Revision	Technical information for service personnel
Interface Board Rev	Technical information for service personnel

### 3.5.5.4 Device Memory

Command	Description
Save Memory Data on Card	Copies the memory data from the internal memory to the memory card.
Erase Memory Data	Erases memory data from the internal memory

The following data is stored and logged:

- Running hours
- Technical alarms
- Settings
- Asset data
- Treatment hours
- Treatment settings
- Device serial number
- Physiological alarms
- Detailed log (containing at least 24 h data of clinical parameters)
- Breath log (containing at least 30 day data of clinical parameters)
- Usage log (containing at least 1 year data of non-clinical events, alarms and settings)

### 3.5.5.5 Clock

Setting	Description
Time	Sets the time for the Vivo 3. The time is used for logs and reports.
Date	Sets the date for the Vivo 3. The date is used for logs and reports.
Time Format	Select whether to use 12 hr clock or 24 hr clock.
Date Format	Select the date format to use.
Alarm Clock	Activates an alarm clock on the Vivo 3. One short signal every 5 seconds.
Alarm Clock Time	Sets the alarm time for the alarm clock.
Alarm Clock Volume	Sets the volume for the alarm clock.

#### NOTE



The alarm clock works only when the ventilator is running on mains.

### 3.5.5.6 Connectivity

This section contain commands for turning radio based connections on or off.

### 3.5.5.7 Compliance Data

Setting / information	Description
Min daily use	The minimum daily use, in hours, to reach compliance per day.
Reset Compliance Data	Reset of data.
Start Date	Start date for treatment.
Total Usage Hours	The total number of hours the ventilator has been running in operating mode during the download period.
Total Days	The total number of days in the download period.
Days with Usage	The number of days in the download period where the ventilator has been running in operating mode (all day, or part of the day).
Average Usage Hours	The average number of hours per day the ventilator has been running in operating mode. Only days where the ventilator has been running in operating mode are part of the value (days without treatment are excluded in the calculation).
Days Compliant	Number and percentage of days when minimum daily use has been reached.
AHI	Apnea Hypopnea Index

### 3.5.5.8 Patient Operating Time








Setting / information	Description
Patient Operating Time	Total hours of usage
Reset Patient Operating Time	Reset the counter to 0.



### 3.5.5.9 About

Setting / information	Description
Regulatory FCC	Contains FCC Compliance statement
Regulatory IC	Contains IC Statment






## 3.6 Symbols on the Vivo 3

### Symbols on the Product Information Label

Symbol	Description
	Manufacturer information
	Nemko certification mark (NRTL/SCC accredited)
<b>IP22</b>	Degree of protection provided by enclosure. See 9.3 <i>Environmental Conditions</i> , page 117 for explanation.
<b>Rx Only</b>	(Symbol only applicable in U.S.)  Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Read user instructions
	RTCA/DO-160 G categorization. Category M This category is defined for equipment and interconnected wiring located in areas where apertures are electro-magnetically significant and not directly in view of radio receiver's antenna. This category may be suitable for equipment and associated interconnecting wiring located in the passenger cabin or in the cockpit of a transport aircraft.
<b>CE</b> 2797	CE marking applies in accordance with directive MDD 93/42/EEC, as amended by 2007/47/ EC.
	Internal Battery Information
	Read the “Disposal” section for information about recycling and disposal, see 21.
	Product number










Symbol	Description
	Serial number
	Date of manufacture

### Symbols on the Vivo 3 Back





Symbol	Description
	Oxygen connection port. Max 30 l/min and 100 kPa.
	Power connection port. Use approved power supplies only.
	USB port
	SD memory card port
	I/O port for accessory box/SpO <sub>2</sub>

## Additional Symbols on Parts and Accessories


This section describes additional symbols for Vivo 3 detachable parts and accessories from Breas Medical. Each item, or its package, have the symbols that applies for the specific accessory.

Symbol	Description
	Caution, hot surface
	Caution symbol, read the Accessory's instructions for more information.
	Keep away from rain.
	Single patient use
	CE marking applies in accordance with directive MDD 93/42/EEC, as amended by 2007/47/ EC.
	Applied part, type BF Electrically connected to the patient but not directly to the heart.
	IEC protection Class II: Double insulated equipment.
	Remote alarm port
	Alarm nurse call port



Symbol	Description
	Effort belt port (not used)
	AC power port
	DC power port
	SPO <sub>2</sub> port

### 3.7 Symbols on the Display



1. SD card inserted and working
2. Nurse call connected
3. Remote alarm connected
4. SpO2 connected
5. Humidifier activated
6. Heated circuit activated
7. “Clock radio wake up” activated
8. All alarms shut off
9. Powered by mains.
10. Powered by external battery (no icon = powered by internal battery)
11. Battery status (red — yellow — green). A flash indicates charging.

## 4 Treatment Functions and Settings

This chapter describes the modes, settings and parameters that controls the ventilation of the Vivo 3.

Lay operators must only use the User Manual, not the Clinician's Manual.

### 4.1 Treatment Modes

This section describes the ventilating modes of the Vivo 3.

#### 4.1.1 PCV+A — Pressure Controlled Ventilation (Assisted)

In PCV+A mode, the ventilation is controlled by set values for pressure, inspiratory time, and rise time.

The inspiration is started by either of the following triggers:

- The patient actively takes a breath.
- The interval set by the breath rate, if the patient does not actively trigger a breath.

The inspiration stops and the expiration starts by either of the following triggers:

- The inspiration time expires.
- The High Pressure Alarm limit is reached.

#### 4.1.2 PSV — Pressure Support Ventilation

In the PSV mode, the patient's spontaneous breathing is assisted by the ventilator. The patient controls the start of inspiration through the inspiratory trigger and the start of exhalation by the expiratory trigger.

The set pressure is used as a target pressure, if the flow is decreased to the expiratory trigger level before the set pressure is reached, the expiration starts.

When an inspiration is started either when the patient triggers a breath, or when the breath rate setting initiates an inspiration in case of a prolonged apnea, the ventilator delivers a flow up to a certain preset pressure limit. In case of a patient-initiated breath, the patient continues the breath for as long as they wish and cycles off when a percentage of drop in peak inspiratory flow (expiratory trigger) has been reached.

Spontaneous breaths stop and an exhalation starts in three cases:

- The inspiration flow has dropped to the value set for expiratory trigger.
- The inspiration time is longer than the limit for maximal inspiration time or when inspiration time reaches 3 seconds.
- The limit for the high-pressure alarm is reached.

### 4.1.3 S — Spontaneous

In the S mode, the patient's spontaneous breathing is assisted by the ventilator. The patient controls the start of inspiration through the inspiratory trigger and the start of exhalation by the expiratory trigger.

The set pressure is used as a target pressure, if the flow is decreased to the expiratory trigger level before the set pressure is reached, the expiration starts.

When an inspiration is started by the patient triggering a breath, the ventilator delivers a flow up to a certain preset pressure limit. The patient continues the breath for as long as he/she wishes and cycles off when a percentage of drop in peak inspiratory flow (expiratory trigger) has been reached.

Spontaneous breaths stop and an exhalation starts in three cases:

- The inspiration flow has dropped to the value set for expiratory trigger.
- The inspiration time is longer than the limit for maximal inspiration time or when inspiration time reaches 3 seconds.
- The High Pressure Alarm limit is reached.

### 4.1.4 S/T — Spontaneous/Timed

In the S/T mode, the patient's spontaneous breathing is assisted by the ventilator. The patient controls the start of inspiration through the inspiratory trigger and the start of exhalation by the expiratory trigger.

The set pressure is used as a target pressure, if the flow is decreased to the expiratory trigger level before the set pressure is reached, the expiration starts.

When an inspiration is started either when the patient triggers a breath, or when the breath rate setting initiates an inspiration in case of a prolonged apnea, the ventilator delivers a flow up to a certain preset pressure limit. In case of a patient initiated breath, the patient continues the breath for as long as they wish and cycles off when a percentage of drop in peak inspiratory flow (expiratory trigger) has been reached.

Spontaneous breaths stop and an exhalation starts in three cases:

- The inspiration flow has dropped to the value set for expiratory trigger.
- The inspiration time is longer than the limit for maximal inspiration time or when inspiration time reaches 3 seconds.
- The limit for the high-pressure alarm is reached.

### 4.1.5 T — Timed

In the T mode the ventilation is controlled by set values for pressure, breath rate, inspiratory time, and rise time.

The inspiration starts by the interval set by the breath rate.

The inspiration stops and the expiration starts by either of the following triggers:

- The inspiration time expires.
- The High Pressure Alarm limit is reached.

### 4.1.6 CPAP — Continuous Positive Airway Pressure

In CPAP mode the Vivo 3 is applying a continuous positive pressure to the airways. The flow will automatically be adjusted to maintain the set CPAP level.

## 4.2 Treatment Settings

This section describes settings and parameters that affects the ventilating function of the Vivo 3. On the machine, these settings are in the **Settings** menu.

### 4.2.1 Ventilation mode

This setting is for selecting the ventilation mode for the Vivo 3.

Property	Description
Available values	The Vivo 3 can be operated in the following modes: <ul style="list-style-type: none"><li>• PCV+A (Assisted Pressure Controlled Ventilation). See page 51.</li><li>• PSV (Pressure Support Ventilation). See page 51.</li><li>• S (Spontaneous). See page 52</li><li>• S/T (Spontaneous/Timed). See page 52</li><li>• T (Timed). See page 52</li><li>• CPAP (Continuous Positive Airway Pressure). See page 53.</li></ul>
Default value	PSV

### 4.2.2 Inspiratory Pressure (IPAP)

The Inspiratory Pressure setting defines the airway pressure during the inspiratory phase.

Inspiratory pressure is also known as Inspiratory Positive Airway Pressure (IPAP).

The air pressure is measured by a sensor, and software regulates the blower that generates the pressure.

Property	Description
Value range	4 cmH <sub>2</sub> O to 40 cmH <sub>2</sub> O.
Default value	10 cmH <sub>2</sub> O

Property	Description
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	<ul style="list-style-type: none"> <li>• PCV+A</li> <li>• PSV</li> <li>• S</li> <li>• S/T</li> <li>• T</li> </ul>

### 4.2.3 Expiratory Positive Airway Pressure (EPAP)

The EPAP setting defines the airway pressure during the expiration phase.

Property	Description
Value range	2 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O. For inspiratory pressure settings below 22 cmH <sub>2</sub> O, the maximum value for EPAP is limited to the inspiratory pressure minus 2 cmH <sub>2</sub> O.
Default value	5 cmH <sub>2</sub> O
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	<ul style="list-style-type: none"> <li>• PCV+A</li> <li>• PSV</li> <li>• S</li> <li>• S/T</li> <li>• T</li> </ul>

### 4.2.4 Continuous Positive Airway Pressure (CPAP)

The CPAP setting defines the pressure that will be applied to the airways.

Property	Description
Value range	4 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O.
Default value	10 cmH <sub>2</sub> O
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	CPAP

#### 4.2.5 Breath Rate

The Breath Rate setting defines the minimum number of breaths the Vivo 3 will deliver as long as no inspiratory trigger effort from the patient is detected. The cycles will be ventilator-initiated breaths.

Property	Description
Value range	0 BPM to 40 BPM The available range for Breath Rate might also be limited by the set values for <i>Inspiratory Time</i> and <i>Target Volume</i> .
Default value	12 BPM
Setting resolution	1 BPM
Available in modes	PCV+A PSV S S/T T

#### 4.2.6 Inspiratory Time (Insp. Time)

The Inspiratory Time setting defines the length of each inspiration from start of inspiration to cycling off to expiration.

Property	Description
Value range	0.3 s to 3 s The available range for Inspiratory Time might also be limited by the set values for <i>Breath Rate</i> .
Default value	1.5 s
Setting resolution	0.1 s
Available in modes	PCV+A PSV S/T T

#### 4.2.7 Minimum Inspiratory Time (Min Insp. Time)

The Minimum Inspiratory Time setting defines a minimum length for each inspiration. If the Minimum Inspiratory Time is set to Off, the length of the inspiration and/or maximum inspiratory time is dependent on the set Expiratory Trigger.

Property	Description
Value range	<ul style="list-style-type: none"><li>Off</li><li>0.3 s to 3 s</li></ul>
Default value	Off

Property	Description
Setting resolution	0.1 s
Available in modes	<ul style="list-style-type: none"> <li>• PSV</li> <li>• S</li> <li>• S/T</li> </ul>

#### 4.2.8 Maximum Inspiratory Time (Max. Insp. Time)

The Maximum Inspiratory Time setting defines a maximum length for each inspiration. If the Maximum Inspiratory Time is set to Off, the length of the inspiration and/or minimum inspiratory time is dependent on the set Expiratory Trigger.

Property	Description
Value range	<ul style="list-style-type: none"> <li>• Off</li> <li>• 0.3 s to 3 s</li> </ul>
Default value	Off
Setting resolution	0.1 s
Available in modes	<ul style="list-style-type: none"> <li>• PSV</li> <li>• S</li> <li>• S/T</li> </ul>

#### 4.2.9 Inspiratory Trigger (Insp.Trigger)

The inspiratory trigger setting defines the patient's effort required to initiate a ventilator assisted breath.

When the patient starts a breath, an increasing flow is created in the patient circuit. If the patient's effort reaches the set inspiratory trigger level an inspiration is initiated. If the patient cannot trigger a breath, the ventilator will deliver breaths according to the set Breath Rate.

The Assisted breath modes in PCV are turned off if the inspiratory trigger is set to Off.

Property	Description
Value range	Auto, 1 to 9
Default value	Auto



Property	Description
Setting resolution	1
Available in modes	<ul style="list-style-type: none"> <li>• PCV+A</li> <li>• PSV</li> <li>• S</li> <li>• S/T</li> </ul>

#### 4.2.10 Expiratory Trigger (Exp.Trigger)

The Expiratory Trigger setting defines the moment when the ventilator will cycle from the inspiratory to the expiratory phase.

Property	Description
Value range	Auto, 1 to 9, corresponding to 10% to 90% of peak flow
Default value	Auto
Setting resolution	1
Available in modes	<ul style="list-style-type: none"> <li>• PSV</li> <li>• S</li> <li>• S/T</li> </ul>

#### 4.2.11 Rise Time

The Rise Time setting controls the speed of the pressure increase from start of inspiration to set inspiratory pressure.

A low setting will give a faster pressure increase and therefore a longer plateau at the set inspiratory pressure. A high setting will give a slow increase and therefore a shorter plateau.

Property	Description
Value range	Auto, 1 to 9
Default value	Auto
Setting resolution	1
Available in modes	<ul style="list-style-type: none"> <li>• PCV+A</li> <li>• PSV</li> <li>• S</li> <li>• S/T</li> <li>• T</li> </ul>

### 4.2.12 Target Volume

The Target Volume setting defines the tidal volume that the Vivo 3 will aim for while ventilating the patient in a pressure mode. This is done by adapting the Inspiratory Pressure between two adjustable pressure limits: Min IPAP and Max IPAP.

Property	Description
Value range	100 ml to 1500 ml. The available range for Target Volume might also be limited by the set values for <i>Breath Rate</i> .
Default value	100 ml
Setting resolution	Below 500 ml: 20 ml Above 500 ml: 50 ml
Available in modes	<ul style="list-style-type: none"><li>• PCV+A</li><li>• PSV</li><li>• S</li><li>• S/T</li><li>• T</li></ul>

### 4.2.13 Max Pressure

The Max pressure setting is only used when Target Volume is activated.

Max pressure defines the upper pressure limit up to where the Vivo 3 can increase the pressure to reach the set Target Volume. If Target Volume is not reached at Max pressure, the Vivo 3 will continue to ventilate at this Max pressure setting.

Property	Description
Value range	4 cmH <sub>2</sub> O to 40 cmH <sub>2</sub> O.
Default value	Same as the value for <i>IPAP</i>
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	<ul style="list-style-type: none"><li>• PCV+A</li><li>• PSV</li><li>• S</li><li>• S/T</li><li>• T</li></ul>

The air pressure is measured by a sensor, and software regulates the blower that generates the pressure.

#### 4.2.14 Min Pressure

The Min pressure setting is only used when Target Volume is activated.

Min pressure defines the lower pressure limit down to where the Vivo 3 can decrease the pressure to maintain the set Target Volume. If the actual volume is above Target Volume at Min pressure, the Vivo 3 will continue to ventilate at this Min pressure setting.

Property	Description
Value range	4 cmH <sub>2</sub> O to 20 cmH <sub>2</sub> O.
Default value	Same as the value for <i>IPAP</i>
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	<ul style="list-style-type: none"><li>• PCV+A</li><li>• PSV</li><li>• S</li><li>• S/T</li><li>• T</li></ul>

#### 4.2.15 Ramp Up

The ramp up setting defines a ramp time for increasing the airway pressures.

At the beginning of Ramp Up time, IPAP will start at 2 cmH<sub>2</sub>O above the initially set EPAP.

Property	Description
Value range	<ul style="list-style-type: none"><li>• Off</li><li>• 10 min to 60 min</li></ul>
Default value	Off
Setting resolution	10 min
Available in modes	<ul style="list-style-type: none"><li>• PCV+A</li><li>• PSV</li><li>• S</li><li>• S/T</li><li>• T</li><li>• CPAP</li></ul>

#### 4.2.16 Ramp Down

The ramp Down setting defines a ramp time for decreasing the airway pressures.

At the end of the Ramp Down time, IPAP will end 2 cmH<sub>2</sub>O above the initially set EPAP.

Property	Description
Value range	<ul style="list-style-type: none"><li>• Off</li><li>• 10 min to 60 min</li></ul>
Default value	Off
Setting resolution	10 min
Available in modes	<ul style="list-style-type: none"><li>• PCV+A</li><li>• PSV</li><li>• S</li><li>• S/T</li><li>• T</li><li>• CPAP</li></ul>

#### 4.2.17 Ramp Pressure

The ramp pressure setting is used for defining the start pressure when Ramp Up is used:

IPAP shall start 2 cmH<sub>2</sub>O above Ramp Pressure, or 2 cmH<sub>2</sub>O above EPAP if Ramp Pressure is lower than EPAP

CPAP mode:

Property	Description
Value range	<ul style="list-style-type: none"><li>• 2 cmH<sub>2</sub>O to the set value for EPAP (PSV and PCV+A modes)</li><li>• 2 cmH<sub>2</sub>O to the set value for CPAP (CPAP mode)</li></ul>
Default value	<ul style="list-style-type: none"><li>• The set value for EPAP (PSV and PCV+A modes)</li><li>• The set value for CPAP (CPAP mode)</li></ul>

Property	Description
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Available in modes	<ul style="list-style-type: none"> <li>• PCV+A</li> <li>• PSV</li> <li>• S</li> <li>• S/T</li> <li>• T</li> <li>• CPAP</li> </ul>

#### 4.2.18 Humidifier

This setting defines whether a humidifier shall be used.

By default, this setting can be modified in Clinical mode only, but the Vivo 3 may be configured to allow modification in home mode as well.



#### WARNING!

Read the section 2.7 *Humidification and Heating — Warnings and Precautions*, page 26 before activating the humidifier, to make sure all conditions are met and considered.

Property	Description
Value range	<ul style="list-style-type: none"> <li>• On</li> <li>• Off</li> </ul>
Default value	Off
Available in modes	All modes

#### 4.2.19 Humidifier Level

This setting defines the level of humidification, if *Humidifier* is set to *On*.

By default, this setting can be modified in Clinical mode only, but the Vivo 3 may be configured to allow modification in home mode as well.

Property	Description
Value range	1 to 5
Default value	1
Setting resolution	1
Available in modes	All modes

## 4.2.20 Heated Circuit

This setting defines whether a heated patient circuit shall be used as patient circuit.

By default, this setting can be modified in Clinical mode only, but the Vivo 3 may be configured to allow modification in home mode as well.



### WARNING!

Read the section 2.7 *Humidification and Heating — Warnings and Precautions*, page 26 before activating the heated circuit, to make sure all conditions are met and considered.

Property	Description
Value range	<ul style="list-style-type: none"><li>• On</li><li>• Off</li></ul>
Default value	Off
Available in modes	All modes

## 4.2.21 Heated Circuit Level

This setting defines the heating level, if *Heated Circuit* is set to *On*.

By default, this setting can be modified in Clinical mode only, but the Vivo 3 may be configured to allow modification also in home mode.


Property	Description
Value range	1 to 5
Default value	1
Setting resolution	1
Available in modes	All modes

## 5 Prepare the Vivo 3 for Use

### 5.1 Checking the Vivo 3 before First Use

Before using the Vivo 3, perform the following checks.

- 1 Ensure that you have the equipment mentioned in 3.1 *Main Components*, page 32
- 2 Ensure that the equipment is in good condition.
- 3 If stored more than 1 month, connect the Vivo 3 to the power supply to recharge the internal battery.
- 4 Check that the air inlet filters are installed.




### 5.2 Placing the Vivo 3




#### WARNING!

Read the chapter 2.3 *Environment — Warnings and Precautions*, page 20 carefully to make sure all conditions are met and considered.

- 1 Place the Vivo 3 on a solid, flat surface. The Vivo 3 should be placed lower than the patient in order to prevent the device from falling on the patient, as well as preventing condensed water from reaching the patient.



- 2 Make sure that nothing can block the patient air inlet at the back of the Vivo 3, such as a curtain etc.



### 5.3 Connecting the Vivo 3 to Power Supply

This chapter describes how to connect a Breas power supply.



#### WARNING!


Read the chapter 2.2 *Electricity — Warnings and Precautions*, page 19 carefully to make sure all conditions are considered and met.



#### CAUTION!


Mains isolation is provided by the approved AC/DC and DC/DC adapters. The ventilator must only be powered by the approved AC/DC or DC/DC adapter (accessory) listed in 10 *Accessories*, page 130.

- 1 Plug the DC plug by turning it into the power port at the back of the Vivo 3.




- 2 Make sure that the mains cord is connected to the power adapter (a).
- 3 Turn the DC plug clockwise 90 degrees (b).
- 4 a) If using mains power, connect the AC plug to the mains supply outlet:





Item 1 in the figure: AC/DC power supply

b) If using DC power, connect the DC plug to the DC power outlet:



Item 2 in the figure: DC/DC power supply



### CAUTION!

Make sure that the power outlet is not blocked, so that the cord can be unplugged without difficulties.

## 5.4 Connecting the Patient Circuit



### WARNING!

Read the chapter 2.4 *Patient Circuit — Warnings and Precautions*, page 22 Carefully to make sure all conditions are considered and met.

The Vivo 3 is intended to be used with leakage circuits only. Recommended leak rate: 20 to 50 liters per minute at 10 cmH<sub>2</sub>O.

- 1 Check that the circuit is clean and undamaged.
- 2 Connect the tube to the air outlet.

If having a heated circuit, make sure to connect the end with the heating plug to the air outlet. Refer to 5.4.1 *Connect the Heated Patient Circuit*, page 66 and 5.4.2 *Disconnect the Heated Patient Circuit from the Ventilator*, page 66.


- 3 Check whether the patient interface has an integrated leakage port.
  - If yes, connect the patient circuit to the patient interface.
  - If no, first connect a leakage valve to the patient circuit, then connect it to the patient interface.



A bacterial filter can be attached on the ventilator outlet connector.

#### 5.4.1 Connect the Heated Patient Circuit

- 1 Check the patient circuit for damages.
- 2 Connect the ventilator cuff to the air outlet of the ventilator. A clicking sound is heard when the latches are fitted correctly.



- 3 If applicable, perform a pre-use test on the ventilator. For specific information, see the User Manual for the ventilator.
- 4 Connect the patient interface cuff to the patient interface.
- 5 On the ventilator, activate the circuit heating. .

#### 5.4.2 Disconnect the Heated Patient Circuit from the Ventilator

- 1 Disconnect from the ventilator:  
Press the latches and pull the circuit off from the ventilator. Do not pull by the tube.
- 2 Disconnect from the patient interface:  
Hold by the cuff and pull the circuit off from the patient interface. Do not pull by the tube.

### 5.5 Power up the Vivo 3

This function check may be performed as a pre-use check or whenever the function of the Vivo 3 needs to be checked.

## Procedure

### Prerequisites

- The Vivo 3 shall be connected to the power supply.
  - The Vivo 3 shall be turned on and in *Standby* mode.
- 1 Press and hold the **On/Off** button until the *Starting Treatment* progress bar is filled.
  - 2 When starting treatment, an automatic alarm signal test is performed. Check that the test is performed successfully, this is indicated by:
    - A short beep indicating functional audio signaling.
    - The alarm LED first lights yellow, then red, indicating functional visual signaling.
    - The audio pause LED lights yellow.
    - In about a second, both LEDs are turned off.

If the test fails, do not use the Vivo 3. Contact your supplier of the Vivo 3 for a technical check.

- 3 Perform a pre-use test (select **Others** —> **Pre-use Test**).
- 4 Disconnect the power cord for more than 5 seconds.

Check that the device switches to the internal battery and that the information message *Switched to Internal Battery* is shown on the screen.

If not, contact your supplier of the Vivo 3.

- 5 Reconnect the power cord.

Check that the Vivo 3 switches to the mains supply.

If not, contact your supplier of the Vivo 3.

## 5.6 Adjusting the Settings

### WARNING!



#### Risk of Asphyxia or Personal Injury

Using incorrect settings may cause personal injury or severe medical conditions, such as hypercarbia, producing arterial acidemia.

The therapy settings shall be based on a physician's description. Changes to settings must be made by authorized clinical personnel only.




For detailed information about treatment modes and settings, see 4 *Treatment Functions and Settings*, page 51.

The Vivo 3 shall be turned on and in *Standby* mode.

Follow the instructions below when setting up the Vivo 3.

- 1 Turn on the Vivo 3.



- 2 Press the **Setup** button.
- 3 Within the limits prescribed by the physician, adjust the settings for the best possible breathing comfort for the patient.

If changing the treatment mode, always review and consider the other settings. The new treatment mode may have additional settings specific for the mode.

- 4 Check that alarms can be received at the operator's position. For detailed information, see 7.1 *Operator's Position*, page 91.
- 5 Make a copy of the *Patient Settings Record* on page 138 and document the settings.


The Vivo 3 is now ready for use. If restarting the Vivo 3, it will start with the settings used when it was powered off.

## 6 How to Use the Vivo 3

For a summary description of the means of initiating and terminating the inspiratory phase in each mode of the ventilator, please refer to 4.1 *Treatment Modes*.

### 6.1 Switch On the Vivo 3

- 1 Make sure that the power supply is connected.
- 2 Press the **On/Off** button.  
⇒The Vivo 3 is now in standby mode with the display turned on.



### 6.2 Switch Between Clinical Mode and Home Mode

The Clinical mode gives access to all functions and settings, while the Home mode gives access only to a limited set of functions and settings. Press and hold down the left and right arrow buttons simultaneously for 3 seconds to switch modes.

### 6.3 Start the Treatment

#### WARNING!




#### Risk of Asphyxia or Personal Injury

Using incorrect settings may cause personal injury or severe medical conditions, such as hypercarbia, producing arterial acidemia.

The therapy settings shall be based on a physician's description. Changes to settings must be made by authorized clinical personnel only.

- 1 Make sure that the Vivo 3 is in *Standby* mode.



- 2 Press and hold the **On/Off** button until the *Starting Treatment* progress bar is filled.

⇒The Vivo 3 now performs a function test indicated by a short beep and then starts the treatment. Check that the self test is performed successfully, this is indicated by:

- A short beep indicating successful alarm signalling
- The alarm LED first lights yellow, then red
- The audio pause LED lights yellow
- In about a second, both LEDs are turned off


If the function test beep absents, take the Vivo 3 out of use and contact your supplier of the Vivo 3 e.



## 6.4 Stop the Treatment


- 1 Press and hold the **On/Off** button until the *Stopping Treatment* progress bar is filled.

⇒A request for a complementary action is now displayed for stopping the treatment. If no complementary action is taken within 6 seconds, the Vivo 3 reverts to normal operation.



2 Do one of the following:

- Press the **Mute Alarm** button (a). (Always available.)
- Press **Ramp** (b) for a ramped stop. (Available if a *Ramp Down* time has been specified.)




## 6.5 Switch Off the Vivo 3

1 Make sure that the treatment is stopped and the Vivo 3 is in *Standby* mode.


2 Press the **On/Off** button.

⇒ When the message “Do you want to turn off the ventilator?” is displayed, confirm by pressing the **Mute Alarm** button within 6 seconds. Otherwise, the Vivo 3 will revert to standby mode.



3 Press the **Mute Alarm** button.

⇒ The Vivo 3 is now turned off.




## 6.6 Monitor Treatment

1 Press the **Monitor** navigation button.

⇒ The Monitor page is now displayed.

The monitor page displays treatment data monitored by the Vivo 3. It consists of a bar graph illustrating the current pressure and text area displaying monitored values in text.



#### Bar graph legend

1. Current pressure
2. Low pressure alarm level
3. Set inspiration pressure level
4. High pressure alarm level

### 6.6.1 Treatment Values Monitored by the Vivo 3

- **Ppeak**

*P<sub>peak</sub>* (*Peak pressure*) is the highest pressure that is recorded during the latest inspiratory phase.

- **EPAP**

*EPAP* (*Expiratory Positive Airway Pressure*) is the lowest pressure that is recorded during the latest expiratory phase.

- **Leakage**

*Leakage* is the average calculated leak (l/min) over the last breath, with a breath by breath update.

- **MVe**



*Mve (Minute Volume, expiratory)* is calculated as the Tidal Volume multiplied with the Total Breath Rate.

- **Vte**

*Vte (Tidal volume, expiratory)* is the expired tidal volume for each breath.

- **% in TgV**

Percentage of breaths within the target volume window.

- **Rise Time**

*Rise Time* displays the duration of the pressure increase during the inspiration phase.

- **SpO<sub>2</sub>**

*SpO<sub>2</sub> (Saturation of Peripheral Oxygen)* displays the patient's oxygen saturation, as measured by the SpO<sub>2</sub> module.

This value is only displayed if the SpO<sub>2</sub> module is connected.

Displayed ranges: According to the manufacturer's specifications.

- **Pulse Rate**

*Pulse Rate* displays the patient's pulse rate as measured by the SpO<sub>2</sub> module.

This value is only displayed if the SpO<sub>2</sub> module is connected.

Displayed ranges: According to the manufacturer's specifications.

## 6.6.2 Treatment Values Monitored by External Equipment

### Monitoring of Expiratory CO<sub>2</sub>

For monitoring expiratory CO<sub>2</sub>, an external monitoring device shall be connected to the patient circuit.

The device shall comply with ISO 80601-2-55.

### Monitoring of Oxygen Saturation

For monitoring the patient air oxygen saturation, an external monitoring device shall be connected to the patient circuit.

The device shall comply with ISO 80601-2-55 and have a high oxygen level alarm.

## 6.7 Using Accessories

This section describes how to use accessories provided by Breas Medical.

## 6.7.1 Using the Attachable Humidifier



### WARNING!

Read the section 2.7 *Humidification and Heating — Warnings and Precautions*, page 26 before using the Vivo 3 with the attachable humidifier.



### WARNING!

#### Risk of Burns

If you remove and open the humidifier just after using the ventilator, ensure that you don't touch the heater inside the humidifier since it can be very hot.



### WARNING!

The humidifier is for single patient use only.



### NOTE

The ventilator can also be used with active stand-alone and passive humidifiers

The humidifier is intended to humidify the patient air. It is intended for non-invasive use only.

### About the Attachable Humidifier

The attachable humidifier is intended for non-invasive use only.


The humidifier requires a connection to the Mains power supply to work.

The information in the below table is applicable to the recommended breathing system configuration, which is Attachable humidifier and Heated circuit.


Property	Value
Humidifier classification	Category 2 humidifier according to ISO 80601-2-74
Rated Flow	20-50 l/min
Operating Conditions	+15°C to +35°C. Humidity: RH from 15% to 95%.
Gas leakage	< 0.2 l/min at 50 cmH <sub>2</sub> O
Max humidification output	> 13 mg/l, tested according to ISO 80601-2-74

## Using the Attachable Humidifier for the First Time — Overview

- 1 Take out the air bypass unit by pressing the locking latch (a) and then pulling it out.




- 2 Fill the humidifier with water.
- 3 Insert the humidifier.



## Detach the Humidifier from the Vivo 3

- 1 If any treatment is running, stop it.
- 2 Push down the locking latch (a) and then pull the humidifier out.



- 3 If you will use the Vivo 3 without the humidifier, install the air bypass unit in place of the humidifier.

## Fill the Humidifier

### Duration of operation between humidifier refills


Humidifier level (5): 12 hours. At lower settings, the duration will be longer.

### CAUTION!



- The water chamber must be detached from the Vivo 3 when filling water into the water chamber.
- Use only distilled or sterilized water or boiled, chilled tap water in the humidifier water chamber. This is to reduce bacteria and mineral deposits.
- Do not fill the water chamber with hot water.
- Do not overfill the water chamber. The water chamber has a capacity of 350 ml and the maximum filling level is indicated on the chamber.
- After using the ventilator, wait one minute before you open the water chamber since it can get hot under certain conditions (for example if the humidifier runs out of water).


- 1 Detach the water chamber, see above.
- 2 Make sure the water meets the quality requirements. It shall be either:
  - Distilled
  - Sterilized
  - Boiled and chilled tap water.
- 3 Hold the humidifier with the air path openings (a) up and fill water into either of the air path openings. Make sure not to overfill (b).



- 4 If the outside of the humidifier is wet, dry it with a lint free cloth before attaching it to the Vivo 3.

### Attach the Humidifier to the Vivo 3

- 1 If any treatment is running, stop it.
- 2 If the air path bypass unit is installed to the Vivo 3, remove it by pressing the locking latch (a) and then pulling it out.




- 3 Make sure the humidifier is correctly filled and push it into the Vivo 3 so the locking latch is engaged.  
⇒A click indicates that the humidifier is correctly installed.

### Activate the Humidifier


- 1 Press the **Setup** button twice.  
⇒The second setup page is now displayed, with humidifier settings.
- 2 Press the **Up arrow** button until the setting **Humidifier** is selected.
- 3 Press the **Right arrow** button to turn the humidifier on.  
⇒The value is shifted to On and the humidifier symbol is lit on the display.  
⇒The setting for Humidifier Level is displayed.
- 4 Press the **Down arrow** to select the **Humidifier Level** setting.
- 5 Press the **Left arrow** button to decrease the humidification or the **Right arrow** button to increase the humidification.

## When not Using the Humidifier

- 1 Take out the humidifier by pushing down the locking latch (a) and then pull the humidifier out.



- 2 Empty the humidifier of water.
- 3 Insert the air path bypass unit.




## Cleaning the Humidifier



### **WARNING!**

#### **Risk of Burns**

If you remove and open the humidifier just after using the ventilator, ensure that you don't touch the heater inside the humidifier since it can be very hot.



- 1 Open the humidifier by pulling the lid up. Don't detach the silicone airpaths (a). Also, don't dismount the cartridge heater (b).
- 2 Wash the humidifier by hand or in a dishwasher.
- 3 After cleaning, ensure that the silicone gasket (c) is positioned correctly: When closing the humidifier; the grooved edge must be positioned downwards.

### 6.7.2 Using the Heated Circuit



#### WARNING!

Read section 2.7 *Humidification and Heating — Warnings and Precautions*, page 26 before using the Vivo 3 with the heated circuit.



#### NOTE

The heated circuit requires connection to the Mains power supply to work.



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for the heated circuit.



For information about connecting a patient circuit with heated circuit, see 5.4 *Connecting the Patient Circuit*, page 65.


### Activate the heated circuit

- 1 Press the **Setup** button twice.  
⇒The second setup page is now displayed, with heated circuit settings.
- 2 Press the **Up arrow** button until the setting **Heated Circuit** is selected.
- 3 Press the **Right arrow** button to turn the heated circuit on.  
⇒The value is shifted to On and the heated circuit symbol is lit on the display.  
⇒The setting for Heated Circuit Level is displayed.
- 4 Press the **Down arrow** to select the **Heated Circuit Level** setting.
- 5 Press the **Left arrow** button to decrease the heating or the **Right arrow** button to increase the heating.

#### 6.7.3 Using the SpO<sub>2</sub> Sensor

The SpO<sub>2</sub> module (consisting of a SpO<sub>2</sub> sensor, an electronic unit) is intended to measure functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate.





The SpO<sub>2</sub> module can be connected to the Vivo 3 (item 1 above) using the SpO<sub>2</sub> adapter cable (007079) or to the accessory box (007000) in order to monitor and store SpO<sub>2</sub> measurements. The SpO<sub>2</sub> measurements will be stored in the data memory which can be downloaded to a PC and viewed in the PC software.



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for SpO<sub>2</sub> module.



### CAUTION!

- When using the Vivo 3 with the SpO<sub>2</sub> sensor, the Vivo 3 displays functional oxygen saturation measured by the sensor.
- The following information concerns the light emitted by the SpO<sub>2</sub>
  - Peak Wavelength (red): 660 nm
  - Peak Wavelength (infrared): 905 nm
  - Maximum Optical Output Power: ≤ 15 mW
  - For more information regarding the oxygen probe's range of peak wavelengths, max optical power and usage, please refer to the respective probe manual.
- Environmental factors may influence the function or accuracy of the pulse oximeter, such as ambient light, physical movement, diagnostic testing, low perfusion, electromagnetic interference, dysfunctional haemoglobin, presence of certain dyes and inappropriate positioning of the pulse oximeter probe.
- A functional test cannot be used to assess the accuracy of a pulse oximeter probe or a pulse oximeter.

#### 6.7.4 Using the Oxygen Supply Adapter



### WARNING!

Read the section 2.9 *Oxygen Usage — Warning and Precautions*, page 28 before using the Vivo 3 with oxygen.

Supplemental oxygen with a flow up to 30 l/min can be added by an oxygen source with rotameter such as oxygen cylinder, central oxygen supply system or an oxygen concentrator.

## Supplemental Oxygen Supply Requirements

Property	Requirement
Maximum flow	30 l/min.
Maximum pressure	100 kPa
Supply source	Source equipped with rotameter. Examples of supply sources: <ul style="list-style-type: none"><li>• Oxygen cylinder</li><li>• Central oxygen supply</li><li>• Oxygen concentrator</li></ul>
Connector	The oxygen source shall be equipped with the Breas Low Pressure Oxygen Adapter, art. no. 005032


### Connect the Oxygen supply

- 1 Connect the oxygen adapter (article no. 005032) to the oxygen supply's tube.
- 2 Connect the oxygen adapter to the oxygen port at the back of the Vivo 3. See 3.3 *Ventilator Back*, page 35 for detailed information.
- 3 If using a device for monitoring the oxygen saturation, connect it according to the supplier instructions.

### 6.7.5 Using the Nurse Call Connection

The Vivo 3 can be connected to a nurse call system through a port on the accessory box. When connected, alarms from the Vivo 3 will be forwarded to the nurse call system.

- 1 Connect the Accessory box to the Communication port at the back of the Vivo 3.
- 2 Connect the nurse call cable to the Nurse call/ Remote alarm port on the accessory box.



⇒The Nurse call symbol is now lit on the display.

- 3 Trigger an alarm on the Vivo 3 and check that it activates the nurse call system. For detailed information about triggering alarms, see 8.4 *Alarm Tests*, page 112.


### 6.7.6 Using the Remote Alarm Unit



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for Remote Alarm.

The Remote Alarm enables care providers and clinical personnel to monitor the Vivo 3 alarms remotely. The Remote Alarm forwards alarms from the Vivo 3.

- 1 Connect the Accessory box to the Communication port at the back of the Vivo 3.
- 2 Connect the remote alarm cable to the Nurse call/ Remote alarm port on the accessory box.



- 3 Start the remote alarm unit.
- 4 Trigger an alarm on the Vivo 3 and check that it activates the remote alarm system. For detailed information about triggering alarms, see 8.4 *Alarm Tests*, page 112.

### 6.7.7 Using the Protective Cover

The protective cover is intended for mobile use of the Vivo 3 in hospitals, institutions and home care environments. It can be used while the Vivo 3 is operating, for example mounted on a wheelchair, in a personal vehicle, or carried by hand.

The protective cover protects the Vivo 3 from environmental impact such as shock, water spill, sunlight, dust and dirt, under normal handling.

**The protective cover does not protect the Vivo 3 from rain or snowfall.**

- The protective cover can be used together with external power supply units.
- The protective cover cannot be used together with the attachable humidifier. When using the protective cover, the air path bypass unit shall be installed in place of the attachable humidifier.

### 6.7.8 Using the Lightweight Mobility Bag

The lightweight Mobility Bag is intended for mobile use of the Vivo 3 in hospitals, institutions and home care environments. It can be used while the ventilator is operating, for example mounted on a wheelchair, in a personal vehicle, or carried by hand.

The lightweight Mobility Bag protects the ventilator from, water spill, sunlight, dust and dirt, under normal handling.

The lightweight mobility bag cannot be used together with the attachable humidifier. When using the lightweight mobility bag, the air path bypass unit shall be installed in place of the attachable humidifier.

It does not protect against environmental impact such as shock, from rain or snowfall.


### 6.7.9 Using the Y-Cable

The Y-cable is used for connecting the Vivo 3 to both mains and external DC at the same time. see 5.3 *Connecting the Vivo 3 to Power Supply*, page 64. When both power sources are available, the mains will be used.

### 6.7.10 Using the Vivo 3 with the Trolley


This section describes how to use the Vivo 3 and a trolley with mounting bracket.

Mount and dismount the Vivo 3 as shown in the picture:



The bottom plate is mounted to the trolley using four screws.

Before mounting the bottom plate, the tilt bar must be dismantled:



Be careful when handling the trolley with the ventilator mounted, in order to avoid any risk of the trolley falling. The trolley can be tipped 10° and return to vertical position, when loaded in accordance with the weight specifications below.

- The total maximum load of the trolley is 24 kg (52 lbs).
- The maximum load of the trolley basket is 0.9 kg (2 lbs).

## 6.8 Transfer Data between the Vivo 3 and a PC

### **WARNING!**



Read the section 2.2 *Electricity — Warnings and Precautions*, page 19.

### **CAUTION!**



Do not eject the memory card or disconnect the USB cable while the Vivo 3 is transferring data. Doing so may result in loss of data and/or damaged equipment.

### **NOTE**



In order to view and present patient data, the Vivo PC Software must be used.



Instructions on how to manage data in the Vivo PC Software can be found in the software's help.

## 6.8.1 Log Data

Data is stored in three different logs:


- Usage log, containing 1 year of data sampled and averaged every 24 hours. Data include non-clinical events, physiological and technical alarms, and settings.
- Breath log, containing 1 month of data sampled every breath. Only sampled during treatment. Data include clinical parameters.
- Detailed log, containing 24 hours of data sampled at 62.5 Hz. Only sampled during treatment. Data include flow and pressure data. This log is stored directly on SD card.

The log data is maintained when device is powered off or in case of power failure. When logs are filled up, the oldest data will be discarded.

## 6.8.2 Transfer Data using a Memory Card

The Vivo 3 has an SD memory card slot at the back of the device.

- 1 Insert an SD card to the memory card slot.



- 2 Open the **Device Memory** page.  
**Others > Device Memory**
- 3 Select **Save Memory Data on Card** and press the right arrow button to start the action.  
⇒ The data transfer starts. When finished, a confirmation message is displayed.  
The logs can be viewed in a PC with the Vivo PC Software installed.
- 4 Remove the memory card from the Vivo 3 and insert it to a memory card reader connected to a PC with the Vivo PC Software installed.
- 5 Start the Vivo PC Software for accessing the data.  
For detailed information about using the Vivo PC Software, see the software's help file.




### 6.8.3 Transfer Data using an USB Cable



#### WARNING!

The PC must be placed outside the patient area (i.e. more than 2 meters from the patient).

- 1 Connect the USB cable to the Vivo 3.  
Make sure it is fitted correctly.



- 2 Connect the other end of the cable to the PC.
- 3 Start the software for viewing or retrieving data.

## 6.9 Basic Troubleshooting

Problem	Action
The ventilator doesn't start.	Check the connection of the power cord between the ventilator and the power outlet.
The ventilator starts, but the patient circuit doesn't deliver any air.	Ensure that both ends of the patient circuit are correctly connected. Straighten the patient circuit or replace it.
The humidifier doesn't work properly.	If the humidifier is incorrectly assembled, disassemble it and then assemble it correctly. If the air is dry despite using the humidifier, increase the level of humidification. If the humidifier doesn't deliver any heat, check that the Vivo 3 is powered from the AC power supply; the humidifier doesn't work on battery.

## 7 Alarms



### **WARNING!**

#### **Risk of Unnoticed Critical Conditions**

- All the physiological alarms of the Vivo 3 must be set at safe levels that will effectively warn the user of any risk.  
The alarm levels should be assessed considering the patient's treatment settings.
- Any change of treatment settings or change of components in the ventilation system may require readjustment of the alarm levels.
- The alarm sound level should be set to a clearly audible level. Setting the alarm sound level below that of the ambient sound level can impede recognition of alarm conditions.



### **CAUTION!**

Never leave a patient unattended during an alarm condition.

### **Distributed Alarm System**

There are two distributed alarm systems for the ventilator:

- The cable connected remote alarm unit provided by Breas Medical
- The nurse call cable provided by Breas Medical

### **Alarm Function During Power Failure**

At power failure, the following will happen:

1. When the mains power is broken, the internal battery takes over the power supply and the following message is displayed: Switched to Internal Battery.
2. When 20 minutes (+-5 minutes) of the battery charge remains, the following alarm (Medium) is displayed: Last Power Source Low Capacity. An audio signal is heard.
3. When 5 – 7 minutes of the battery charge remains, the following alarm (High) is displayed: Critical Low Last Power Source. An audio signal is heard.
4. When 0 minutes of the battery charge remains, the following message is displayed: The internal battery is completely discharged. Connect mains or DC to charge.
5. An alarm is triggered, and the Vivo 3 is powered off.

The alarm settings are maintained during power failure.

## 7.1 Operator's Position

The alarm priority indications are designed to be recognized from a distance of 4 meters and by an angle of 50 ° from the normal of the Vivo 3 display.

### 7.1.1 Checking the Operator's Position

- 1 Activate an alarm. For detailed instructions, see 8.4 *Alarm Tests*, page 112.
- 2 From the operator's position, make sure that the audible alarm signal is heard and that it is possible to recognize the alarm priority level by either of the visual signals (flashing LED lights and display message on the screen). The sound pressure level range is 50-73 dBA for high and medium priority alarm signals.

For detailed information, see 7.2.1 *Identify an Alarm Condition*, page 91.

- 3 If the test fails, consider the following actions:
  - Find a better position for the operator.
  - Adjust the alarm sound level, see 7.2.5 *Adjust the Alarm Sound Level*, page 94.
  - Add a remote alarm unit to the system, see 10 *Accessories*, page 130.

## 7.2 Handle Alarms

### 7.2.1 Identify an Alarm Condition

If an alarm condition is detected, the Vivo 3 main unit and the remote alarm unit (if connected) will alarm without delay. The alarms will remain active until the alarm condition is resolved.

Active alarms are indicated by:

- Audible signal, see page 91.
- Alarm message on the screen, see page 92.
- The alarm LED, see page 92.

#### Alarm Audio signal

- **High priority alarms**  
3 short signals followed by 2 more after 0.5 s. The signal sequence repeats every 3rd second.



- **Medium priority alarms**  
3 signals. The signal sequence repeats every 6th second.



For information about adjusting the alarm sound level, see 7.2.5 *Adjust the Alarm Sound Level*, page 94.

## Information Message Audio Signal

Short signal every 5 seconds



## Alarm Sound Pressure Level

The alarm sound pressure level is adjustable within 55–80 dBA.

## Alarm Message on the Screen

The name of the active alarm is displayed on the screen.

- **High priority alarms**  
Red highlight color.
- **Medium priority alarms**  
Yellow highlight color.



For detailed information about specific alarms, see 7.3 *Physiological Alarms*, page 95 and 7.4 *Technical Alarms*, page 104.

## Display of Multiple Alarms


If several alarms are active, active high priority alarms have precedence over medium priority alarms: All high priority alarm conditions must be resolved before any medium alarms are displayed.

If several alarms of the same priority are active at the same time, the alarm descriptions are looped in the display.

A “>>” symbol indicates that more alarms are to be displayed in the loop.

## Alarm LED signal


- **High priority alarms**  
Red light, flashing quickly (0.5 s. interval).
- **Medium priority alarms**  
Yellow light, flashing slowly (2 s. interval).



## 7.2.2 Pause the Alarm Sound

The audible signal can be paused for 60 seconds by pressing the Audio Pause button. Pressing the button again reactivates the sound.

If a new alarm condition occurs during the audio pause period, the audible signal will be reactivated.



## 7.2.3 Reset an Alarm

To reset an alarm, correct the cause of the alarm condition.

⇒Once the cause is corrected, the alarm disappears from the display.



### WARNING!

If an alarm condition cannot be corrected, take the Vivo 3 out of use and contact your supplier of the Vivo 3 e.



## 7.2.4 View Historical Alarms

To view historical alarms, press the Alarm button three times.




## 7.2.5 Adjust the Alarm Sound Level


- 1 Press the Alarm navigation button twice.  
⇒The second alarm page is now displayed.




- 2 Press the Down or Up arrow button to select the Alarm Sound Level setting.



- 3 Press the Left or Right arrow buttons to adjust the sound level.



- Press the Up arrow button to finish the adjustment by deselecting the setting.



- When finished with the sound level adjustment, check that the alarm can be received at the operator's position, see 7.1 *Operator's Position*, page 91.

## 7.3 Physiological Alarms

### 7.3.1 High Pressure Alarm

Property	Description
Alarm text	<b>High Pressure</b>
Priority	High
Alarm condition	A High Pressure alarm will be given when the patient pressure reaches the set High Pressure alarm limit for three consecutive breaths.
Possible cause	<ul style="list-style-type: none"> <li>Mismatch between pressure setting and alarm setting.</li> <li>Coughing during inspiration.</li> <li>Changes in airway resistance and or compliance.</li> </ul>
Reset criteria	A full breath is performed with maximum pressure below the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings. The actual breath is however terminated if the High Pressure alarm limit is reached.
Setting range	<ul style="list-style-type: none"> <li>5 cmH<sub>2</sub>O to 45 cmH<sub>2</sub>O</li> <li>Off</li> </ul> <p>Note that the High pressure alarm cannot be set lower than the value set for the Low pressure alarm.</p>
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Setting Display	The High Pressure alarm setting is displayed by a red line in the pressure bar graph.



### 7.3.2 Low Pressure Alarm

Property	Description
Alarm text	<b>Low Pressure</b>
Priority	High
Alarm condition	A Low Pressure alarm will be given when the Vivo 3 pressure fails to reach the low pressure alarm limit for 15 seconds.
Possible cause	<ul style="list-style-type: none"> <li>• Disconnection of patient circuit.</li> <li>• Mismatch between pressure setting and alarm setting.</li> <li>• Leakage from the mask or other components of the patient circuit.</li> </ul>
Reset criteria	The pressure rises above the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"> <li>• 1 cmH<sub>2</sub>O to 40 cmH<sub>2</sub>O</li> <li>• Off</li> </ul> <p>Note that the Low pressure alarm cannot be set higher than the value set for IPAP or the High pressure alarm.</p>
Setting resolution	Below 10 cmH <sub>2</sub> O: 0.5 cmH <sub>2</sub> O Above 10 cmH <sub>2</sub> O: 1.0 cmH <sub>2</sub> O
Setting display	The Low Pressure alarm setting is displayed by a red line in the pressure bar graph.





### 7.3.3 High Breath Rate Alarm

Property	Description
Alarm text	<b>High Breath Rate</b>
Priority	Medium
Alarm condition	A High Breath Rate alarm will be given when the alarm limit has been exceeded for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Mismatch between the Breath Rate setting and the alarm setting.</li><li>• Increased Breath Rate.</li><li>• Too sensitive setting of the inspiratory trigger setting.</li></ul>
Reset criteria	The breath rate goes below the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 10 bpm to 50 bpm</li><li>• Off</li></ul>
Setting resolution	1 bpm.

### 7.3.4 Low Breath Rate Alarm

Property	Description
Alarm text	<b>Low Breath Rate</b>
Priority	High
Alarm condition	A Low Breath Rate alarm will be given when the delivered total breath rate is below the alarm limit for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Mismatch between the Breath Rate setting and the alarm setting.</li><li>• The patient cannot trigger breaths because the inspiratory trigger setting is too high.</li><li>• Decrease in the patient's spontaneous breathing.</li><li>• Circuit disconnection.</li></ul>
Reset criteria	The breath rate goes above the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 4 bpm to 30 bpm.</li><li>• Off</li></ul>
Setting resolution	1 bpm.

### 7.3.5 High Minute Volume Alarm

Property	Description
Alarm text	<b>High MV</b>
Priority	Medium
Alarm condition	A High Minute Volume alarm will be given when the minute volume is above the set alarm limit for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Disconnection of patient circuit.</li><li>• Mismatch between pressure setting and alarm setting.</li><li>• Leakage from the mask or other components of the patient circuit.</li><li>• Increased breath rate.</li></ul>
Reset criteria	The minute volume goes below the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 1 l. to 40 l.</li><li>• Off</li></ul>
Setting resolution	0.5 l.

### 7.3.6 Low Minute Volume Alarm

Property	Description
Alarm text	<b>Low MV</b>
Priority	High
Alarm condition	A Low Minute Volume alarm will be given when the minute volume is below the set alarm limit for 15 1 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Disconnection of patient circuit.</li><li>• Mismatch between pressure setting and alarm setting.</li><li>• Leakage from the mask or other components of the patient circuit.</li><li>• Decreased breath rate.</li></ul>
Reset criteria	The minute volume goes above the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 1 l. to 30 l.</li><li>• Off</li></ul>
Setting resolution	0.5 l.

### 7.3.7 Rebreathing Alarm

Property	Description
Alarm text	<b>Rebreathing</b>
Priority	High
Alarm condition	A Rebreathing alarm will be given if the intentional leakage is too low for more than 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Obstructed or occluded patient circuit.</li><li>• Incorrect patient circuit.</li></ul>
Reset criteria	The leakage is back within limits.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• On</li><li>• Off</li></ul>

### 7.3.8 Apnea Alarm

Property	Description
Alarm text	<b>Apnea</b>
Priority	High
Alarm condition	An Apnea alarm will be given when no patient-triggered breath is detected for the set period of time.
Possible cause	<ul style="list-style-type: none"><li>• Patient stopped breathing.</li><li>• Patient decreases spontaneous breathing.</li><li>• Circuit disconnection.</li><li>• Inspiratory Trigger is set too high.</li></ul>
Reset criteria	Inspiratory effort detected by the Vivo 3.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 5 to 60 s.</li><li>• Off</li></ul>
Setting resolution	5 s below 15 s. 15 s above 15 s. 60 s above 60 s.

### 7.3.9 Disconnection Alarm



#### CAUTION!

No single alarm can reliably detect all disconnections due to the number of possible combinations of therapy settings, circuit configurations and patient interfaces. To verify that patient disconnection can be detected, including if the patient interface becomes accidentally detached from the patient, it is advised to test the functionality of the Disconnection Alarm upfront with the complete set-up as used during treatment, including items such as filters, circuit, connectors, interface (mask, cannula etc.)

Property	Description
Alarm text	<b>Disconnection</b>
Priority	High
Alarm condition	A Disconnection alarm will be given when the measured flow exceeds the expected leakage flow at the set Pressure for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Too high leakage in the patient circuit.</li><li>• The patient has removed the mask.</li><li>• Circuit disconnection.</li><li>• Pilot pressure tube disconnection</li></ul>
Reset criteria	The leakage is back within limits.
Ventilator action	The Vivo 3 will continue treatment according to the current settings
Setting range	<ul style="list-style-type: none"><li>• On</li><li>• Off</li></ul>

### 7.3.10 High EPAP Alarm

Property	Description
Alarm text	<b>High EPAP</b>
Priority	Medium
Alarm condition	A High EPAP alarm will be given when EPAP has gone 30% above the set value for 3 breaths.
Possible cause	<ul style="list-style-type: none"><li>• Blocked leakage port.</li><li>• Too short expiratory time.</li><li>• Changes in airway resistance and or compliance.</li></ul>
Reset criteria	EPAP has gone below the alarm limit (lower than 30% above the set value).
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• On</li><li>• Off</li></ul>

### 7.3.11 Low EPAP Alarm

Property	Description
Alarm text	<b>Low EPAP</b>
Priority	Medium
Alarm condition	A Low EPAP alarm will be given when EPAP has gone 30% below the set value for 3 breaths.
Possible cause	<ul style="list-style-type: none"><li>• Excessive leakage.</li></ul>
Reset criteria	EPAP has gone above the alarm limit (higher than 30% below the set value).
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• On</li><li>• Off</li></ul>

### 7.3.12 High SpO<sub>2</sub> Alarm

Property	Description
Alarm text	<b>High SpO<sub>2</sub></b>
Priority	Medium
Alarm condition	A High SpO <sub>2</sub> alarm will be given when the measured SpO <sub>2</sub> exceeds the alarm limit for 30 seconds.
Possible cause	Too high flow of bleed-in oxygen.
Reset criteria	The SpO <sub>2</sub> value goes back below the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	<ul style="list-style-type: none"><li>• 80 % to 100 %</li><li>• 90 % to 100 %</li><li>• Off</li></ul>
Setting resolution	1 %

This alarm requires a connected SpO<sub>2</sub> sensor.

### 7.3.13 Low SpO<sub>2</sub> Alarm

Property	Description
Alarm text	<b>Low SpO<sub>2</sub></b>
Priority	High
Definition	A Low SpO <sub>2</sub> alarm will be given when the measured SpO <sub>2</sub> is below the alarm limit for 30 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Too low flow of bleed-in oxygen.</li><li>• Oxygen inlet is disconnected.</li><li>• Delivered tidal volumes are too small.</li></ul>
Setting range	70% to 100%
Setting resolution	1%
Ventilator action	The ventilator will continue treatment with the same settings.

This alarm requires a connected SpO<sub>2</sub> sensor.

### 7.3.14 High Pulse Rate Alarm

Property	Description
Alarm text	<b>High Pulse Rate</b>
Priority	Medium
Alarm condition	A High Pulse Rate alarm will be given when the measured pulse rate exceeds the alarm limit for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Insufficient ventilatory support.</li><li>• Too low flow of bleed-in oxygen.</li><li>• The EPAP value is set too high.</li><li>• Bad positioning of the finger probe.</li></ul>
Reset criteria	The pulse rate goes back below the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	20 to 250 bpm (beats per minute) Off
Setting resolution	5 bpm (beats per minute)

This alarm requires a connected SpO<sub>2</sub> sensor.

### 7.3.15 Low Pulse Rate Alarm

Property	Description
Alarm text	<b>Low Pulse Rate</b>
Priority	High
Alarm condition	A low pulse rate alarm will be given when the measured pulse rate goes below the alarm limit for 15 seconds.
Possible cause	<ul style="list-style-type: none"><li>• Bad positioning of the finger probe.</li><li>• Too low flow of bleed-in oxygen.</li><li>• Insufficient ventilatory support.</li></ul>
Reset criteria	The pulse rate goes back above the alarm limit.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Setting range	20 to 250 bpm (beats per minute) Off
Setting resolution	5 bpm (beats per minute)

This alarm requires a connected SpO<sub>2</sub> sensor.

## 7.4 Technical Alarms

### 7.4.1 High Pressure Limitation Alarm

Property	Description
Alarm text	<b>High Pressure Limitation</b>
Priority	High
Alarm condition	The High Pressure Limitation alarm is given if the high pressure alarm limitation limit (60 cmH <sub>2</sub> O) is reached.
Possible cause	<ul style="list-style-type: none"><li>• Mismatch between pressure setting and alarm setting.</li><li>• Coughing during inspiration.</li><li>• Changes in airway resistance and or compliance.</li></ul>
Reset criteria	A full breath with pressure below the set alarm limit.
Ventilator action	The current breath is terminated and then the Vivo 3 will continue treatment according to the current settings.

### 7.4.2 Power Fail Alarm

Property	Description
Alarm text	The display is turned off. The Power Fail alarm is indicated by the alarm LED and the audible signal only.
Priority	High The display is turned off. The Power Fail alarm is indicated by the alarm LED and the audible signal only.
Alarm condition	The Power Fail alarm is given if the last power source fails to provide enough power for running the ventilator.
Possible cause	The last available power source cannot deliver power to the ventilator. Battery discharged or battery failure.
Reset criteria	External power supply connected to ventilator.
Ventilator action	The Vivo 3 stops the treatment, turns off the display and gives the Power Fail alarm for at least 2 minutes. If power is restored within the alarm time, the ventilator will automatically resume treatment with current settings. When powered up again, the power failure will be logged.

### 7.4.3 SpO<sub>2</sub> Sensor Failure / Disconnected Alarm

This alarm requires a connected SpO<sub>2</sub> sensor.



Property	Description
Alarm text	<b>SpO2 Disconnected</b>
Priority	Medium
Alarm condition	An SpO <sub>2</sub> Sensor Failure/Disconnection alarm will be given if one of the conditions below appears: <ul style="list-style-type: none"> <li>• An error signal is received from the sensor</li> <li>• No signal at all from the sensor is received within 2 seconds.</li> </ul>
Possible cause	Faulty or disconnected sensor.
Reset criteria	Normal communication with the sensor is reestablished. An information message remains until acknowledged by the user.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.

#### 7.4.4 SpO<sub>2</sub> Artifact

This alarm requires a connected SpO<sub>2</sub> sensor.

Property	Description
Alarm text	<b>Poor SpO2 Signal</b>
Priority	Medium
Alarm condition	A poor SpO <sub>2</sub> Signal alarm will be given if perfusion is too low or artifacts are detected by the sensor.
Possible cause	Check the sensor and its placement on the patient.
Reset criteria	An OK signal is received from the sensor or the sensor is disconnected. An information message remains until acknowledged by the user.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.

#### 7.4.5 Ambient Pressure Compensation Lost Alarm

Property	Description
Alarm text	<b>Pressure Comp Lost</b>
Priority	Medium
Alarm condition	An Ambient Pressure Compensation Lost alarm will be given when the automatic ambient pressure compensation functionality is out of order.
Ventilator action	The Vivo 3 will continue treatment according to the current settings. Normal atmospheric pressure at sea level will be used as approximation for the temporary ambient pressure compensation. If used at other altitude, delivered and measured pressures may deviate.
Reset	Reset of ventilator.

#### 7.4.6 Low Battery Alarm

Property	Description
Alarm text	<b>Low Battery. Connect to Mains.</b>
Priority	Medium
Alarm condition	This alarm will be given when the last battery source (internal battery) has 15 minutes of operating time left with current settings.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Reset Criteria	Connection of another power source within 15 minutes or the treatment will shut down.

#### 7.4.7 Critically Low Battery Alarm

Property	Description
Alarm text	<b>Critically Low Battery. Connect to Mains</b>
Priority	High
Alarm condition	This alarm will be given when the last battery source (internal battery) has 5 minutes of operating time left with current settings.
Ventilator action	The Vivo 3 will continue treatment according to the current settings.
Reset Criteria	Connection of another power source within 5 minutes or the treatment will shut down.

### 7.4.8 High Patient Air Temp. (High Patient Air Temperature)

Property	Description
Alarm text	<b>High Patient Air Temp</b>
Priority	High
Alarm condition	A High Patient Air Temperature alarm will be given when the patient air temperature exceeds 43°C (109.4°F). If having a patient circuit with an active heated circuit, the patient air temperature is measured by the circuit's temperature sensor. Otherwise, the temperature is measured by the flow measurement sensor inside the Vivo 3.
Possible cause	<ul style="list-style-type: none"><li>• Blocked air inlets.</li><li>• Too high ambient temperature.</li></ul>
Ventilator action	The Vivo 3 will continue treatment. If a heated circuit or attachable humidifier is used, these will be turned off. The ventilator stops treatment and gives alarm for up to 2 minutes.
Reset criteria	The temperature goes below the limit again.

### 7.4.9 Flow Sensor Failure

Property	Description
Alarm text	<b>Flow Sensor Failure</b>
Priority	Medium
Alarm condition	No data or erroneous data from the flow sensor
Possible cause	
Reset criteria	Correct data from the sensor is received again. An information message remains until acknowledged by the user.
Ventilator action	The Vivo 3 will continue treatment but with the following limitations: <ul style="list-style-type: none"><li>• Monitoring of leakage is disabled.</li><li>• Volume measurements are disabled.</li><li>• The patient cannot trigger breaths (applies to assisted modes)</li></ul> An information message about the limitations is displayed on the screen.

#### 7.4.10 Internal Function Failure

Property	Description
Alarm text	<b>Int. Function Failure(01–18)</b>
Priority	High
Alarm condition	Failure of internal function that prevents treatment or normal operation of the ventilator. The error code that follows the alarm text indicates the kind of function failure.
Reset criteria	Correct function is restored.
Ventilator action	The ventilator will stop the treatment.
Action to take	Restart the Vivo 3. If the alarm persists or reoccurs: Take a note of the error code and contact your supplier of the Vivo 3 .

## 8 Cleaning and Maintenance

This chapter contains instructions for cleaning and maintenance actions that can be carried out by the care provider or by users with physical ability and working knowledge of the system.



### **WARNING!**

#### **Risk of Personal Injury**

- Repairs, upgrades and modifications must be carried out by technicians authorized by Breas Medical only and in accordance with instructions from Breas Medical
- The Vivo 3 must not be opened, repaired or modified by unauthorized personnel.

If subjected to unauthorized operations, Breas Medical is no longer responsible for the performance and safety of the device and all warranties will become invalid.

The patient-connected parts and the filters must be cleaned and replaced regularly to ensure correct function of the Vivo 3. All replaced parts must be disposed of in accordance with local environmental regulations regarding the disposal of used equipment and waste.

### 8.1 Cleaning the Vivo 3



### **WARNING!**

#### **Risk of Electric Shock**

Disconnect the power supply before cleaning the Vivo 3 according to the instructions in this manual.

Do not soak the Vivo 3 or immerse it into any fluids.

#### 8.1.1 Clean the Main Unit Externally

##### **Equipment**

- A lint free cloth.
- A mild soap solution or Ethanol 70%.

- 1 Turn off the Vivo 3 and disconnect the power supply.
- 2 Remove the patient circuit.
- 3 If any cable connected accessories (like the SpO<sub>2</sub> sensor or the accessory box) are used, disconnect them.
- 4 Clean the outside of the Vivo 3 using a lint free cloth moistened with a mild soap solution and / or ethanol 70%.
- 5 When the equipment is clean and dry, reconnect the patient circuit and any accessories that was disconnected during the cleaning.

### 8.1.2 Air Pathway Disinfection

The table below lists the parts that might get contaminated by exhaled gases or bodily fluids during normal use or single fault condition.

Condition	Parts
With bacteria filter	<ul style="list-style-type: none"><li>• Patient circuit</li><li>• Bacteria filter</li></ul>
Without bacteria filter	<ul style="list-style-type: none"><li>• Patient circuit</li><li>• Patient air outlet</li><li>• Air bypass unit/humidifier</li><li>• Blower unit</li><li>• Air inlet filters and filter holder</li></ul>

In case of contamination, the internal air pathways of the Vivo 3 may be disinfected up to 10 times by a maximum 60 minute long validated ozone gas process.

Low resistance bacteria filter, if used, should be replaced every 24 hours.

### 8.1.3 Clean the Patient Circuit



#### CAUTION!

The cleaning and replacement intervals should be established by the care provider, based on the care provider's infection control procedures and the instructions from the patient circuit's manufacturer.

Check the patient circuit regularly for damage. In case of damage, replace the circuit.



The patient circuit should be cleaned and replaced in accordance with the manufacturer's instructions, or by the care provider's instructions if additional cleaning or other replacement intervals are prescribed.


## 8.2 Clean and Replace Patient Air Inlet Filters

The Vivo 3 patient air inlet filters are located inside a magnetic filter holder at the back of the ventilator. The table below describes the filters and their minimum maintenance intervals.

Filter	Maintenance Intervals (minimum)*
Air inlet filter, grey (coarse)	<ul style="list-style-type: none"> <li>• Wash: every week.</li> <li>• Replace: every year or when assigning the Vivo 3 to a new patient.</li> </ul>
Air inlet filter, white (fine)	<ul style="list-style-type: none"> <li>• Replace: every 4th week or when assigning the Vivo 3 to a new patient.</li> </ul>


\* If the Vivo 3 is used in an environment with high grades of pollen or pollutions, shorter intervals might be required.

- 1 Turn off the Vivo 3 and place it on a clean dust free surface.
- 2 Pull out the filter holder and remove the filters.



- 3 If required by the interval or if visibly dirty, wash the grey coarse filter:
  1. Wash the filter using warm water and a mild soap.
  2. Rinse thoroughly.
  3. Dry the filter by first squeezing it in a towel and then letting it dry in the air. Do not wring the filter.
  4. Check that the filter is undamaged and completely dry before reinstalling it.

- 4 When reinstalling the air inlet filters in the filter holder: first install the grey coarse filter, then the white fine filter.



- 5 Reinstall the filter holder on the Vivo 3.

## 8.3 Change of Patients

If assigning the Vivo 3 to a new patient, first clean and prepare the Vivo 3 as described in this section.

If a humidifier is used, it is classified as single patient use.

### NOTE



If the Vivo 3 is used in a clinic, with frequent change of patients, a low resistance bacteria filter is located between the air outlet and the patient circuit to prevent patient cross contamination.

- 1 Clean the outside of the Vivo 3 according to 8.1.1 *Clean the Main Unit Externally*, page 109.
- 2 If a low resistance bacteria filter has been used between the patient air outlet and the patient circuit: Replace the filter.  
  
If a low resistance bacteria filter not has been used, the airways of the Vivo 3 may need disinfection to prevent patient cross contamination. See 8.1.2 *Air Pathway Disinfection*, page 110.
- 3 Replace the patient air inlet filters according to 8.2 *Clean and Replace Patient Air Inlet Filters*, page 110.
- 4 use a new patient circuit for the new patient.

## 8.4 Alarm Tests

The alarm tests should be performed when changing patient, if the ventilator's function needs to be checked for any other reason, or at least every 12 months.

### 8.4.1 Mandatory Test

#### 8.4.1.1 Low Pressure Alarm Test

This test ensures safe use of the Vivo 3.

- 1 Connect the Vivo 3 to a test lung.
- 2 Take notes of the current settings, if they are to be used also after the test.
- 3 Adjust the setting as follows:

Setting	Value
Ventilation mode	PSV
IPAP	15 cmH <sub>2</sub> O
EPAP	5 cmH <sub>2</sub> O



Rise time	9
Insp. trigger	9
Exp. trigger	3
Min insp. time	Off
Max insp. time	Off
Breath rate	12 bpm
Insp. time	2.0 S

All alarm settings shall be set to Off, if possible.

- 4 Start the treatment.
- 5 Set the **High Pressure Alarm** level to Off.
- 6 Set the **Low Pressure Alarm** level to 20 cmH<sub>2</sub>O.  
⇒ The Low Pressure Alarm shall be given.

## 8.4.2 Optional Tests

In this procedure, methods for additional alarm tests are described. These tests are optional and not needed to ensure safe use of the Vivo 3.

### 8.4.2.1 High Pressure Alarm

- 1 Connect the Vivo 3 patient circuit to a test lung and a CPAP device.
- 2 Set High Pressure alarm to 10 cmH<sub>2</sub>O.
- 3 Set the CPAP device treatment pressure to 12 cmH<sub>2</sub>O.
- 4 Adjust the Vivo 3 settings as follows:
  - Ventilation Mode: T
  - Target Volume: Off
  - IPAP: 10 cmH<sub>2</sub>O
  - EPAP: 5 cmH<sub>2</sub>O.
  - Breath Rate: 12 bpm
- 5 Start the treatment on both the Vivo 3 and the CPAP device.
- 6 Wait until the High Pressure alarm is activated.
- 7 Stop the treatment. The test is completed.

### 8.4.2.2 Low Pressure Alarm and the Disconnection Alarm

1. Start the treatment
2. Disconnect the patient circuit.
3. Wait 15 seconds.  
⇒The Low Pressure Alarm and/or the Disconnection Alarm will be given.

### 8.4.2.3 Low Breath Rate Alarm

- 1 Connect the Vivo 3 to a patient circuit.
- 2 Set the Low Breath alarm to 10 bpm.
- 3 Adjust the Vivo 3 settings as follows:
  - Ventilation Mode: T
  - Target Volume: Off
  - IPAP: 10 cmH<sub>2</sub>O
  - EPAP: 5 cmH<sub>2</sub>O.
  - Breath Rate: 12 bpm
- 4 Start the treatment on the Vivo 3 device.
- 5 Wait until the Low Breath Rate alarm is activated.
- 6 Stop the treatment. The test is completed.

### 8.4.2.4 High EPAP Alarm

- 1 Connect the Vivo 3 patient circuit to a test lung and a CPAP device.
- 2 Set the CPAP device treatment pressure to 10 cmH<sub>2</sub>O.
- 3 Adjust the Vivo 3 settings as follows:

Setting	Value
Ventilation Mode	PCV+A
Target Volume	Off
IPAP	15 cmH <sub>2</sub> O
EPAP	5 cmH <sub>2</sub> O
Breath Rate	12 bpm
Insp. Time	1.5 s

Rise Time	5
Insp. Trigger	Off

- 4 Start treatment on both the Vivo 3 and the CPAP device.
- 5 Wait three breaths before the High EPAP alarm.
- 6 Stop treatment. Test completed.

#### **8.4.2.5 Low EPAP Alarm**

1. Make sure Low EPAP alarm is set On and start the treatment.
2. Disconnect the patient circuit from the ventilator.
3. Wait for the duration of 3 breaths before the Low EPAP alarm shall be given.
4. Stop the treatment, the test is completed.

#### **8.4.2.6 High MV Alarm**

1. Connect the patient circuit and test lung.
2. Start the treatment.
3. Note the MVe value in the ventilator monitor.
4. Set the High MV Alarm limit below the monitored MVe.
5. Wait approximately 15 seconds before the High MV alarm shall be given.
6. Stop the treatment, the test is completed.

#### **8.4.2.7 Low MV Alarm**

1. Connect patient circuit and test lung.
2. Start treatment.
3. Note MVe value in ventilator monitor.
4. Set Low MV Alarm limit above monitored MVe.
5. Wait approximately 15 seconds before the Low MV alarm shall be given.
6. Stop treatment, test completed.

### **8.5 Repair**

Repairs of the Vivo 3 must be carried out by authorized personnel only, in accordance with instructions from Breas Medical.

Authorised service workshops can order the Vivo 3 service manual that contains all technical documentation required for the maintenance of the Vivo 3.

## 9 Technical Specifications

### 9.1 Ventilator Size and Weight

Property	Value
Dimensions (WxHxD)	166 x 185 x 200 mm
Weight	1,8 kg

### 9.2 Power Supply

#### Mains Power Supply

Property	Value
Mains Power Supply	100–240 V AC tolerance: +10%/-20%, 50 to 60 Hz, max 1.2 A. The approved AC/DC supply listed in 10 <i>Accessories</i> , page 130 must be used.
Protection against electric shock	Class II ME Equipment

#### External DC Power Supply

Property	Value
External DC Supply	12–24 V isolated DC The approved DC/DC supply listed in 10 <i>Accessories</i> , page 130 must be used.

#### Internal Battery

Property	Value
Battery type	4 cell Li-Ion battery
Nominal Voltage	14.4 V
Capacity	1.5 Ah
Operational time	2 hours
Expected life:	300 charging cycles



#### CAUTION!

Check the number of charging cycles regularly by selecting **Others** —> **Device information**. If 300 cycles have been exceeded, a warning appears here.




## 9.3 Environmental Conditions

Environmental Condition	Specification
Normal Operation Temperature	+5°C to +40°C
	<b>Precautions</b>
	<ul style="list-style-type: none"><li>• Be careful not to position the Vivo 3 in an extra warm place, such as in direct sunlight or above a radiator.</li><li>• Caution should be exercised if the room temperature is higher than 36°C (97°F). The air flow for breathing produced by the Vivo 3 can be as much as 4°C (7°F) higher than the ambient room temperature.</li><li>• Ventilator enclosure: 52°C NC, 53C SFC, touch time &lt; 1 min (NC=normal condition, SFC= single fault condition)</li><li>• Accessory box enclosure: 42°C NC, touch time &lt; 1 min</li><li>• Y-cable box: 42°C NC, touch time &lt; 1 min</li><li>• Humidifier enclosure: 46°C NC, 55C SFC, touch time &lt; 10 min</li><li>• Heated circuit: 42°C NC, 59C SFC, touch time &gt;= 10 min</li></ul>
Extended Operation Temperature	-20°C to +5°C
	The Vivo 3 is operational during the extended operation temperature for 4 hours, if: <ul style="list-style-type: none"><li>• The Vivo 3 is first started within the normal operation temperature span.</li><li>• The Vivo 3 is placed in it's protective cover.</li><li>• This condition happens maximum once a day.</li><li>• The ambient air is dry and still.</li></ul>




### WARNING!

The humidifier is heated only when the ventilator is operating on mains, so don't have the humidifier connected during mobile ventilator use.


Environmental Condition	Specification
<b>Transport and Storage Temperature</b>	<ul style="list-style-type: none"> <li>+5°C to +45°C (Maximum 90 days)</li> <li>-25°C to +70°C (Maximum 30 days)</li> </ul> <p>If stored in temperatures outside normal operation conditions, let the Vivo 3 acclimate before taking it to use.</p> <div style="border-left: 2px solid black; padding-left: 10px; margin-top: 10px;"> <p><b>CAUTION!</b></p> <p> The Vivo 3 must not be stored in a warm place, such as direct sunlight or close to a radiator.</p> <p> When the ventilator is brought from minimum/maximum allowed storage temperature, ensure that it is warmed up/cooled down for one hour before starting it.</p> </div>
<b>Humidity</b>	RH from 15% to 95%, non-condensing.
<b>Ambient Pressure Range</b>	70 to 106 kPa This corresponds to ~ 315 m below sea level to ~3000 m above sea level
<b>Ventilator ingress protection</b>	<p>IP 22</p> <p>Mechanical ingress protection: protected from touch by fingers and objects greater than 12.5 mm</p> <p>Liquid ingress protection: The device withstands dripping water (equivalent to 3 mm rainfall /minute) when not tilted more than 15 degrees from vertical. The protection has been tested for 10 minutes (2.5 minutes in every tilt direction).</p> <div style="border-left: 2px solid black; padding-left: 10px; margin-top: 10px;"> <p><b>CAUTION!</b></p> <p> There is a silicone lid to protect the USB, SD card, and communication ports. The IP22 classification is applicable only when this lid is in place. However, the accessory box can be connected with retained IP22 classification, but then only the lower part of the silicone lid can be opened. Ensure that the silicone lid on the back of the ventilator is closed when no accessories are connected.</p> </div>

## 9.4 Pneumatic Diagram




Number	Description
1	Vivo 3
2	Patient circuit
3	Leakage port / Patient interface connection
4	Patient

### With Bypass Unit



Number	Description
1	Air inlet (Ambient room air intake)
2	Blower
3	Sensors
4	Air bypass unit
5	Air outlet (Breathing air to the patient)

### With Humidifier



Number	Description
1	Air inlet (Ambient room air intake)
2	Blower
3	Sensors
4	Humidifier
5	Air outlet (Breathing air to the patient)
6	Cartridge heater

The humidification output is controlled by the ventilator by regulating the power to the cartridge heater in the humidifier.

## 9.5 Technical Data

### Noise Levels

Property	Value
Static sound pressure level	< 30 dB(A) @ 10 cmH <sub>2</sub> O, 4 mm leak. Measured at a distance of 1 meter, according to standard.
Static sound pressure level, stand by mode	< 22 dB(A) SPL. Measured at a distance of 1 meter, according to standard.
Maximum sound pressure level	40 dB(A) SPL according to applicable standards. Measured at a distance of 1 meter, according to standard.
Maximum sound power level	48 dB (A) according to applicable standards. Measured at a distance of 1 meter, according to standard.



## Maximum Flow Rate at Continuous Pressure

Pressure (hPa)	Flow Rate (l/min)
4	76
8	105
12	158
16	157
20	149

## Maximum Pressure

Property	Value
Maximum pressure limit in the event of fault condition	60 cmH <sub>2</sub> O

## CPAP Maximum Dynamic Pressure Variations

Device with standard tube / Device with worst case VBS (external humidifier with 22 mm tube and bacterial filter)

Pressure [cmH <sub>2</sub> O]	10 bpm	15 bpm	20 bpm
4	0.3/1.0	0.5/1.2	0.7/1.6
8	0.3/1.0	0.5/1.2	0.7/1.6
12	0.3/1.0	0.5/1.2	0.7/1.6
16	0.3/1.0	0.5/1.2	0.7/1.6
20	0.3/1.0	0.5/1.2	0.7/1.6

## Inlet Air Filter Specification

Specification for Air inlet filter, white (part no 007202):

Performance Characteristic	Value	Tested in Accordance with
NaCl Penetration	<7.35 %	BS EN 13274-7 NaCl at 16 cm/sec
Air flow resistance (pressure drop) flow	12.5 Pa typical 20 Pa max 0.3 Pa max	at 82.5 l/min (100 cm <sup>2</sup> ) at 0.2 cm/sec

## Ventilatory Breathing System Characteristics

The ventilator is verified to maintain specified accuracies over these ranges.

Property	Value
Resistance	0.4 to 1.7 cmH <sub>2</sub> O at 40 l/min
Compliance	Max 1.1 cmH <sub>2</sub> O

### 9.6 Data Accuracy



#### NOTE

All stated tolerances include measurement uncertainty (see separate table). The accuracies have been tested with all allowed configurations including worst case configuration (22 mm patient circuit, bacterial filter and external humidifier). Stated tolerances only disclose the maximum tolerance.

Parameter	Tolerance	Remark
Pressure at patient end	$\pm (0.5 \text{ cmH}_2\text{O} + 5 \%)$	Steady state
P <sub>peak</sub>	$\pm (0.5 \text{ cmH}_2\text{O} + 10 \%)$	Monitored in display
EPAP	$\pm (0.5 \text{ cmH}_2\text{O} + 10 \%)$	Monitored in display
Leakage	$\pm 10 \%$	Monitored in display (*) BTPS
MVe	$\pm (10\text{ml} + 10\%) \times \text{Total Breath rate}$	Monitored in display (*) BTPS
V <sub>te</sub>	$\pm (10\text{ml} + 10\%)$	Monitored in display (*) BTPS
% in TgV	$\pm 1 \%$	Monitored in display
Total Breath rate	$\pm 1 \text{ bpm}$	Monitored in display
Spontaneous Breath rate	$\pm 1 \text{ bpm}$	Monitored in display
% Spont breath rate	$\pm 1 \%$	Monitored in display
I:E	$\pm 0.1 \text{ unit}$	Monitored in display
Inspiration Time	$\pm 0.1 \text{ sec}$	Monitored in display
Rise Time	$\pm 10 \%$ or $\pm 0.1 \text{ sec}$ (Which-ever is greatest)	Monitored in display

(\*) BTPS — Body Temperature and Pressure Saturated.

## Measurement Uncertainty

Parameter	Uncertainty
Pressure	$\pm 0.75\%$ of reading or $\pm 0.1$ cmH <sub>2</sub> O
Leakage (flow)	$\pm 10\%$ or $\pm 0.1$ sl/min
Tidal volume V <sub>te</sub>	$\pm 1.75\%$ or $\pm 0.20$ sml
Minute volume M <sub>Ve</sub>	$\pm 2.5\%$ or $0.02$ sl/min
Breath rate	$\pm 2.5\%$ or $\pm 1$ bpm
Inspiration time, T <sub>i</sub>	$\pm 0.02$ s
I:E	$\pm 2.5\%$

### NOTE



sl/min = standard liter per minute (calculated using STP conditions of 21 °C and 1013 mbar).

## Filtering/Smoothing Techniques

Measurement	Description
Pressure	Low pass average time constant 16 ms
Inspiration trigger	Differential mass flow resolution 4 ms
Expiration trigger	Flow low pass filtering with level sensing
SpO <sub>2</sub>	No data post-processing done by the ventilator

## 9.7 Emission and Immunity Declaration

According to IEC 60601-1-2:2014.

The performance of all functions of the Vivo 3 is considered as essential performance for the purpose of immunity testing.

## Vivo 3 Essential Performance

The Vivo 3 will deliver ventilation at the patient-connection port within its published accuracy specifications and within the alarm limits set by the operator, or generate any of the following alarms:

- High pressure
- Low pressure
- Low inspired minute volume
- high inspired minute volume
- Low breath rate
- High leakage
- Low last power source
- Power failure
- Disconnection
- High EPAP
- Low EPAP

The Vivo 3 will provide SpO<sub>2</sub> and pulse rate values within its published accuracy specifications and generate an alarm upon a low SpO<sub>2</sub> condition. The Vivo 3 will provide indication when the SpO<sub>2</sub> value or pulse rate is potentially incorrect, and generate an alarm condition to indicate when the SpO<sub>2</sub> value update period has exceeded 30 seconds.

Under the immunity test conditions, the following allowances are acceptable:

- Temporary degradation of performance that does not adversely affect basic safety or essential performance.
- Any temporary degradation of SpO<sub>2</sub> performance following transient immunity test exposure shall recover from any disruption within 30 seconds.

Additionally, the following shall not be allowed:

- permanent damage or unrecoverable loss of function
- changes in programmable parameters or settings
- reset to default settings
- change of operating mode
- initiation of unintended operation
- negative false alarm condition

### 9.7.1 Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Vivo 3 is intended for use in the electromagnetic environment specified below. The customer or the user of the Vivo 3 should assure that it is used in such an environment.

Immunity test	Compliance level	Electromagnetic environment — Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	The relative humidity should be at least 5%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±2 kV for input/output lines	Mains power quality should be that of a typical commercial, hospital and residential environment.
Surge IEC 61000-4-5	± 0.5, ± 1 kV, ± 2 kV line to line	Mains power quality should be that of a typical commercial, hospital and residential environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial, hospital and residential environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% $U_T$ , 0.5 cycle (multiple phase analysis); 5% $U_T$ , 1 cycle; 70% $U_T$ , 25/30 cycles (50/60 Hz); 0% $U_T$ , 5 sec (50/60 Hz);	The Vivo 3 runs on internal battery during voltage dips, short interruptions and voltage variations on power supply input lines.

## NOTE



$U_T$  is the mains voltage prior to application of the test level.

## WARNING!



Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Vivo 3, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Immunity test	IEC 60601 test level	Recommended separation distance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> (6 V <sub>rms</sub> inside ISM/ASR bands)	$d=0.35*\sqrt{P}$ m at 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	$d= 0.6*\sqrt{P}$ m at 80 MHz to 800 MHz $d= 1.2*\sqrt{P}$ m at 800 MHz to 2.5 GHz <b>Equation description:</b> P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with this symbol.

- At 80 MHz and 800 MHz, the higher frequency range applies.
  - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Vivo 3 is used exceeds the applicable RF compliance level above, the Vivo 3 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Vivo 3.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

### 9.7.2 Guidance and Manufacturer's Declaration – Electromagnetic Emission

The Vivo 3 are intended for use in the electromagnetic environment specified below. The customer or the user of the Vivo 3 should assure that it is used in such an environment.

Immunity test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Group 1	The Vivo 3 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Vivo 3 are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emission IEC 61000-3-3	Complies	

### 9.7.3 Recommended separation distances between portable and mobile RF communications equipment and the Vivo 3

The Vivo 3 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Vivo 3 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Vivo 3 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to the frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 0.35 \cdot \sqrt{P}$ m	80 MHz to 800 MHz $d = 0.6 \cdot \sqrt{P}$ m	800 MHz to 2.5 GHz $d = 1.2 \cdot \sqrt{P}$ m
0.01	0.035	0.06	0.12
0.1	0.11	0.19	0.36
1	0.35	0.60	1.2
10	1.1	1.9	3.6
100	3.5	6.0	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

#### Notes

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

#### 9.7.4 Recommended separation distances between external power conductors and the Vivo 3

Rated maximum current in conductor (A)	Separation distance (m)
	50-60 Hz: $d = I/2\pi H = I/188$
1	0.005
10	0.05
30	0.16

For conductors rated at a maximum current not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation  $d = I/2\pi H$ , where  $I$  is the maximum current rating of the conductor in amperes (A) according to the transmitter manufacturer;  $H$  is the Vivo 3 immunity compliance level to electromagnetic fields in the 50-60 Hz frequency span (30 A/m).

## 9.8 Delivery Settings

Alarm	Default Setting
Low Pressure	10 cmH <sub>2</sub> O
High Pressure	25 cmH <sub>2</sub> O
High Pressure	Off
High Minute Volume	Off
Low Breath Rate	Off
High Breath Rate	Off
Apnea	Off
Low Pulse Rate	Off
High Pulse Rate	Off
Low SpO <sub>2</sub>	90%
High SpO <sub>2</sub>	Off
Disconnection (of mask)	On
Rebreathing (blocked mask leakage)	On
Low EPAP	Off
High EPAP	Off
Alarm Sound Level	5



<b>Setting</b>	<b>Default Value</b>
Time format	24h
Date format	dd/mm/yyyy
Pressure unit	cmH <sub>2</sub> O
Backlight setting	Delayed (2min/30s depending on power source and mode)
Alarm clock enabled	Off
Alarm clock time	00:00
Alarm clock sound level	3
Ventilation mode	PSV
IPAP	10 cmH <sub>2</sub> O
EPAP	5 cmH <sub>2</sub> O
CPAP	10 cmH <sub>2</sub> O
Breath Rate	12 bpm
Inspiration Time	1.5 sec
Min Insp. Time	Off
Max Insp. Time	Off
Inspiratory Trigger	Auto
Expiratory Trigger	Auto
Rise Time	Auto
Max IPAP	Set IPAP
Min IPAP	Set IPAP
Target Volume	100 ml
Ramp Up	Off
Ramp Down	Off
Ramp pressure	Set EPAP (PCV, PSV, S, S/T, T) or set CPAP (CPAP)
Humidifier On/Off	Off
Humidifier Level	5
Heated Circuit On/Off	Off
Heated Circuit Level	1

## 10 Accessories

The accessories described in this section, together with the *medical electric equipment* defined in chapter 3.1 *Main Components*, page 32, constitute the Vivo 3 *medical electric system*.



### CAUTION!

#### Responsibility for System

Accessory equipment connected to the analogue and digital interfaces must be certified according to the respective IEC standards (e.g. IEC 60950 or IEC 62368 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all configurations must comply with the valid version of the system standard IEC 60601-1. Anybody who connects additional equipment to the signal input part or signal output part is configuring a medical system, and is therefore responsible for ensuring that the system complies with the requirements of the valid version of the system standard IEC 60601-1. If in doubt, consult the technical service department or your local representative.


A PC, when connected to the ventilator, shall comply with IEC 62368-1, IEC 60950-1 or IEC 60601-1.

### 10.1 Accessories/Consumables

#### Circuit: Single limb 22 mm, disposable

**Function:** Delivers air to the patient, applied part


**Part No:** 005060



#### Circuit: Single limb 15 mm, disposable

**Function:** Delivers air to the patient, applied part


**Part No:** 006712



### Heated Circuit

**Function:** Delivers air to the patient, applied part. Prevents rain-out.

**Part No:** 006990



### Leakage Port


**Function:** Providing a leakage for clearing exhaled gases.

**Part No:** 007243 (10 pieces)

### Attachable humidifier

**Function:** Humidifies the patient air. For non-invasive use only.


**Part No:** 006977



### Air Bypass Unit

**Function:** Directs the patient air flow, if the attachable humidifier is not used.


**Part No:** 006983



### Filter Holder

**Function:** Holder for air inlet filters


**Part No:** 007598



### **Air Inlet Filter, Grey**

**Function:** Coarse air inlet filter, user replaceable part. Long life (washable).


**Part No:** 007203 (5 pieces)



### **Air Inlet Filter, White**

**Function:** Fine air inlet filter, user replaceable part. Disposable.

**Part No:** 007202 (5 pieces)




### **Low resistance bacterial filter (303 Respirgard-II Filter)**

**Function:** Filter air at ventilator outlet

#### **Characteristics**

- Resistance: 1.8 cmH<sub>2</sub>O @ 60 l/m
- Deadspace: 30 ml
- BFE (Bacterial Filtration Efficiency): 99.9%
- VFE (Viral Filtration Efficiency): 99.8 %



**Part No:** 004185

## **10.2 Accessories**

### **RRC CAR70M DC/DC Power Supply**

**Function:** Power supply adapter for the Vivo 3 shall be used for connecting the Vivo 3 to an external DC source.

**Part No:** 006995

### **Y Cable**

**Function:** Power supply cable, for connecting to both AC and DC power supply.

**Part No:** 007006

### **Protective Cover**

**Function:** Shock protection

**Part No:** 007014

## Lightweight Mobility Bag

**Function:** For mobile use of the Vivo 3 in hospital, institutions and home care environments.

**Part No:** 007380

## Trolley

**Function:** Mobile ventilator use

**Part No:** 007384

## Mounting Bracket

**Function:** Bracket to mount the Vivo 3 to the trolley


**Part No:** 006998

## Accessory Box

**Function:** For connecting measurement and communication accessories:

- Nurse call cable or Remote alarm
- SpO<sub>2</sub> sensor (Might also be connected directly to the Vivo 3, if no other measurement or communication accessories are used.)


**Part No:** 007000



## Remote alarm with cable

**Function:** Monitor Vivo 3 alarms remotely

**Part No:** 10 m: 006348, 25 m: 006349



### **Nurse call cable**

**Function:** Connect the ventilator to a hospital nurse call system


**Part No:**

NO: 006365

NC: 006364

10 k $\Omega$ , NO: 006363


10 k $\Omega$ , NC: 006362



### **SpO<sub>2</sub> module**

**Function:** Connection interface

**Part No:** 006369




### **SpO<sub>2</sub> sensor**

**Function:** Finger Clip SpO<sub>2</sub> sensor

**Part No:**

Adult: 006589


Paediatric: 006590



### **SpO<sub>2</sub> sensor**

**Function:** Multisite SpO<sub>2</sub> sensor


**Part No:** 006591



### **SpO<sub>2</sub> Adapter Cable**

**Function:** Connection cable. For use of SpO<sub>2</sub> without accessory box.

**Part No:** 007079



### **PC Software**

**Function:** Support software for follow-up on patient treatment.

### **USB Cable**


**Function:** USB cable for transferring data between a PC and the Vivo 3.

**Part No:** 005757

### **Memory card**

**Function:** Storage and transfer of settings, patient data and usage data


**Part No:** 006705



### **Memory card reader/writer**

**Function:** Read/write memory card

**Part No:** 002185



### **Oxygen Adapter, Low Pressure**

**Function:** Connection of oxygen supply.

**Part No:** 005032





# Appendices

# A Patient Settings Record

This page can be copied and used for noting the patient's settings.

Patient: \_\_\_\_\_

Clinic: \_\_\_\_\_

Vivo SN: \_\_\_\_\_

Set by: \_\_\_\_\_

Date: \_\_\_\_\_

## Treatment Settings

Ventilation  
Mode: \_\_\_\_\_

IPAP: \_\_\_\_\_

EPAP: \_\_\_\_\_

CPAP: \_\_\_\_\_

Breath Rate: \_\_\_\_\_

Insp. Time: \_\_\_\_\_

Min Insp.  
Time: \_\_\_\_\_

Max Insp.  
Time: \_\_\_\_\_

Insp. Trigger: \_\_\_\_\_

Exp. Trigger: \_\_\_\_\_

Rise Time: \_\_\_\_\_

Max Pressure: \_\_\_\_\_

Min  
Pressure: \_\_\_\_\_

Target  
Volume: \_\_\_\_\_

Ramp Up: \_\_\_\_\_

Ramp Down: \_\_\_\_\_

Humidifier: \_\_\_\_\_

Humidifier  
Level: \_\_\_\_\_

Heated  
Circuit: \_\_\_\_\_

Heated Cir-  
cuit Level: \_\_\_\_\_

## Notes

---

---

---

## **B FAA Compliance**

To whom it may concern:

The US Department of Transportation (DOT) Final Rule, “Nondiscrimination on the Basis of Disability in Air Travel” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009 provides important requirements for the accommodation of passengers with respiratory assistive devices (Ventilators, Respirators and CPAP machines).

In line with these requirements, respiratory assistive devices may be used onboard an aircraft, without further testing by the carrier, provided they have been tested for Electromagnetic Compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M.

Breas Medical has successfully completed testing for the ventilator System. The ventilator System complies with RTCA/DO-160, Section 21, Category M and can be considered FAA compliant.

Some airlines may require advance notification before travel, and devices may need to be operated by battery. Breas Medical recommends that customers check with their airline.



# Index

## A

AC Voltage	
supply .....	116
Accessories .....	130
using .....	73
AHI .....	45
Air bypass unit .....	32
Air outlet.....	35
Air pathway	
disinfection .....	110
Alarm.....	44
Ambient Pressure Compensation Lost ...	106
Apnea .....	99
Critically Low Battery. Connect to	
Mains .....	106
Disconnection.....	100
Flow sensor failure.....	107
High breath rate .....	97
High EPAP .....	101
High minute volume .....	98
High Patient Air Temp (High Patient Air	
Temperature).....	107
High pressure.....	95
High pressure limitation.....	104
High pulse rate.....	103
High SpO <sub>2</sub> .....	102
Internal Function Failure.....	108
Low Battery .....	106
Low breath rate.....	97
Low EPAP .....	101
Low minute volume .....	98
Low pressure .....	96
Low pulse rate.....	103
Low SpO <sub>2</sub> .....	102
Power Failure Alarm .....	104
Rebreathing.....	99
resetting .....	93
SpO <sub>2</sub> artifact .....	105
SpO <sub>2</sub> sensor failure / disconnected .....	104
Alarm condition	
identifying.....	91
Alarm message.....	92
Alarm signal	
LED .....	92
Alarm sound	
pausing.....	93
Alarm tests.....	112
Alarm time.....	44
Alarm volume .....	44

Alarms .....	90
physiological.....	95
Ambient Pressure Compensation Lost	
alarm.....	106
Apnea alarm.....	99
Average usage hours.....	45
<b>B</b>	
Battery .....	36
Battery percentage.....	43
Battery time left .....	43
Battery,(internal)	
specifications.....	116
Board revision.....	43
Bootloader version.....	43
Breaths	
minimum number .....	55
Buttons .....	34
<b>C</b>	
Carry bag .....	32
Caution	
Icon .....	14
Cleaning.....	27, 109
Clinical mode and Home mode	
switch.....	69
Clock.....	44
Command	
execute .....	40
Communciation port.....	35
Continuous Positive Airway Pressure	
(CPAP).....	54
Continuous Positive Airway Pressure,	
CPAP .....	53
Contraindications .....	13
Critically Low Battery. Connect to Mains	
alarm.....	106
<b>D</b>	
Data	
transferring .....	87
Data accuracy.....	122
Date.....	44
Date format .....	44
Days with usage .....	45
DC Voltage	
supply .....	116
Device information.....	43
Device memory.....	43
Device settings .....	42
Disconnection alarm.....	100
Disinfection	

air pathway.....	110	Inspiratory time.....	51, 55
Display.....	34	Inspiratory trigger.....	56
Display light.....	42	Interface board revision.....	43
Disposal.....	21	Internal Function Failure alarm.....	108
<b>E</b>		Invasive use.....	25
Electromagnetic compatibility.....	20	<b>L</b>	
Electromagnetic emission.....	126	Languages.....	42
Electromagnetic immunity.....	124	Light intensity.....	42
Electrostatic discharge.....	20	Log data.....	88
Environmental conditions.....	117	Low Battery alarm.....	106
Erase memory data.....	43	Low breath rate alarm.....	97
Expiration Positive Airway Pressure (EPAP).....	54	Low EPAP alarm.....	101
Expiratory trigger.....	57	Low pressure alarm.....	96
<b>F</b>		Low pulse rate.....	103
Filter		Low SpO <sub>2</sub> alarm.....	102
usage.....	24	<b>M</b>	
Filters		Mains power cord.....	32
clean, replace.....	110	Mains power supply.....	32
Firmware version.....	43	Maintenance.....	27, 109
Front panel.....	33	Memory card	
<b>H</b>		transferring data.....	88
Heated circuit		Memory card port.....	35
setting.....	62	Memory data.....	44
using.....	79	Menus.....	37
Heating.....	26	Min daily use.....	45
High breath rate alarm.....	97	Mobile use	
High EPAP alarm.....	101	safety.....	30
High Patient Air Temperature Alarm (High Patient Air Temp) alarm.....	107	Monitor treatment.....	71
High pressure alarm.....	95	<b>N</b>	
High pulse rate alarm.....	103	Note	
High SpO <sub>2</sub> alarm.....	102	Icon.....	14
Home mode and Clinical mode		Nurse call connection	
switch.....	69	using.....	83
Humidification.....	26	<b>O</b>	
Humidifier		Operating time.....	45
setting.....	61	Oxygen.....	82
Humidifier, attachable		usage.....	28
using.....	74	Oxygen concentration.....	30
Hypoventilation.....	97	Oxygen port.....	35
<b>I</b>		Oxygen saturation	
Icon		monitoring.....	73
Caution.....	14	<b>P</b>	
Note.....	14	Pathway	
Reference.....	14	internal disinfection.....	110
Warning.....	14	Patient air inlet filter.....	32
Inspiratory pressure.....	53	Patient air outlet.....	35
		Patient circuit.....	32
		connect.....	65

Patients		Time	44
changing	112	Time format	44
Placing	63	Tolerances	
Power Failure Alarm alarm	104	settings	122
Power management	36	Total days	45
Power port	35	Total usage hours	45
Power supply		Transit operation	
connect	64	safety	30
Pressure		Treatment	
max at Tgv	58	start	69
min at Tgv	59	Treatment data	40
Pressure Controlled Ventilation (Assisted),		Treatment modes	51
PCV+A	51	<b>U</b>	
Pressure Support Ventilation, PSV	51	Undesirable side effects	13
Pressure unit	42	Usage time	45
Preuse check	66	Internal battery	116
Product name	43	USB cable	
Protective cover	85	transferring data	89
<b>R</b>		USB port	35
Ramp down	60	Using instructions	69
Ramp up	59	<b>V</b>	
Rebreathing alarm	99	Ventilation mode	53
Reference		Voltage	
Icon	14	supply	116
Remote alarm unit		Volume	
using	84	target	58
Repair	115	<b>W</b>	
Reset compliance data	45	Warning	
RF communications equipment	127	Icon	14
Rise time	51	Wheel chair	
<b>S</b>		safety	30
S/T	52	<b>Y</b>	
Save memory data on card	43	Y-cable	86
SD card port	35		
Settings	51, 53		
adjust	67		
change	39		
delivery	128		
Setup pages	41		
Sound level, alarm	94		
SpO <sub>2</sub> sensor			
using	80		
Start date	45		
Stop the treatment	70		
Switch off	71		
Switch on	69		
Symbols	46		
<b>T</b>			
Technical specifications	116		