

BREAS

**User manual
Nippy 4+**

Table of Contents

| | | |
|-------|--|----|
| 1 | Introduction..... | 9 |
| 1.1 | What is the Nippy 4+? | 9 |
| 1.2 | Intended Use | 11 |
| 1.3 | Contraindications | 12 |
| 1.4 | Intended Environment..... | 12 |
| 1.5 | Intended Users..... | 12 |
| 1.5.1 | Respiratory Health Care Specialists..... | 12 |
| 1.5.2 | Lay Operators | 13 |
| 1.5.3 | Service Personnel..... | 13 |
| 1.6 | About this Manual | 14 |
| 1.6.1 | Audience | 14 |
| 1.6.2 | Icons in this Manual | 14 |
| 1.7 | Manufacturer Information..... | 15 |
| 2 | Safety Information | 17 |
| 2.1 | General User Precautions..... | 17 |
| 2.1.1 | Requirements for Life Supporting Treatment | 20 |
| 2.2 | Electrical Safety | 21 |
| 2.3 | Environmental Conditions..... | 22 |
| 2.4 | Usage of Patient Circuit | 23 |
| 2.5 | Usage of Filters | 26 |
| 2.6 | Humidification | 26 |
| 2.7 | Cleaning and Maintenance | 28 |
| 2.8 | Usage of Oxygen | 29 |
| 3 | Product Description | 31 |
| 3.1 | Main Components..... | 31 |
| 3.2 | Front View | 33 |
| 3.3 | Side Views..... | 34 |
| 3.3.1 | Detaching and Reattaching the Side Panels | 35 |
| 3.4 | Equipment Designation..... | 38 |
| 3.4.1 | Additional Symbols | 39 |
| 4 | Preparing the Nippy 4+ for Use | 41 |
| 4.1 | Checking the Nippy 4+ before First Use | 41 |
| 4.2 | Placing the Nippy 4+ | 41 |
| 4.3 | Connecting the Nippy 4+ to Mains..... | 43 |
| 4.4 | Connecting the Patient Circuit | 44 |
| 4.5 | Inspecting the Nippy 4+ before Use | 45 |
| 4.6 | Adjusting the Nippy 4+ Patient Settings..... | 46 |
| 4.7 | Performing the Pre-use Test | 47 |
| 5 | How to Use the Nippy 4+ | 48 |
| 5.1 | Switch the Nippy 4+ On and Off..... | 48 |

| | | |
|--------|--|----|
| 5.1.1 | Switch On and Enter Operating Mode | 48 |
| 5.1.2 | Stop Treatment | 49 |
| 5.1.3 | Turn off / Enter Sleep Mode | 49 |
| 5.2 | Navigating the Display | 50 |
| 5.3 | Symbols Used on the Display | 51 |
| 5.4 | Display Overview | 52 |
| 5.4.1 | Home Mode..... | 52 |
| 5.5 | Profiles | 52 |
| 5.5.1 | Selecting a Profile | 52 |
| 5.6 | Treatment Settings..... | 53 |
| 5.6.1 | Changing a Setting | 53 |
| 5.6.2 | Settings | 53 |
| 5.7 | The Alarms List | 63 |
| 5.8 | The Monitored Values Pane..... | 63 |
| 5.9 | The Menu | 65 |
| 5.9.1 | Opening the Menu | 66 |
| 5.9.2 | Humidification Settings | 66 |
| 5.9.3 | User Preferences | 66 |
| 5.9.4 | Compliance Data (Optional) | 67 |
| 5.9.5 | Pre-Use Test | 67 |
| 5.9.6 | Alarm/Event Logs..... | 68 |
| 5.10 | Transferring Data between the Ventilator and a PC | 69 |
| 5.10.1 | Transferring Data with a Memory Card..... | 69 |
| 5.10.2 | Transferring Data with a Data Cable..... | 70 |
| 5.11 | Using Batteries..... | 70 |
| 5.11.1 | Power Source Priority | 71 |
| 5.11.2 | Charging the Batteries | 71 |
| 5.11.3 | Battery Icons | 72 |
| 5.11.4 | Internal Battery..... | 73 |
| 5.11.5 | Click-in Battery..... | 73 |
| 5.11.6 | Battery Operating Time (Internal and Click-in) | 74 |
| 5.11.7 | Storing the Internal Battery and the Click-in Battery..... | 75 |
| 5.11.8 | External DC..... | 75 |
| 5.12 | Using Accessories | 76 |
| 5.12.1 | Connecting and Disconnecting the Cables..... | 76 |
| 5.12.2 | Using the ventilator with a Nurse Call System..... | 77 |
| 5.12.3 | Using the ventilator with the FiO ₂ Sensor | 77 |
| 5.12.4 | Using the ventilator with the Remote Alarm..... | 79 |
| 5.12.5 | Using the ventilator with the EtCO ₂ Sensor | 79 |
| 5.12.6 | Using the ventilator with the PtcCO ₂ Cable | 84 |
| 5.12.7 | Using the Ventilator with the SpO ₂ module..... | 85 |

| | | |
|---------|--|-----|
| 5.12.8 | Using the Ventilator with the Effort Belts..... | 85 |
| 5.12.9 | Using the ventilator with the Remote Start/Stop..... | 87 |
| 5.12.10 | Using the ventilator with the Protective Cover..... | 88 |
| 5.12.11 | Using the Nippy 4+ with the Trolley..... | 89 |
| 5.12.12 | Using the Click-in Humidifier..... | 90 |
| 5.12.13 | Installing the Water Chamber..... | 91 |
| 5.12.14 | Detaching the Water Chamber..... | 93 |
| 5.12.15 | Using the Patient Circuit with Heated Wire..... | 96 |
| 6 | Alarms..... | 98 |
| 6.1 | Alarm Function..... | 98 |
| 6.1.1 | Alarm Indication..... | 98 |
| 6.1.2 | Audible Signal Pause and Reactivation..... | 99 |
| 6.1.3 | Alarm Reset..... | 100 |
| 6.2 | Operator's Position..... | 100 |
| 6.3 | Physiological Alarms..... | 100 |
| 6.3.1 | High Flow Alarm..... | 101 |
| 6.3.2 | Low Flow Alarm..... | 102 |
| 6.3.3 | High Pressure Alarm..... | 103 |
| 6.3.4 | Low Pressure Alarm..... | 104 |
| 6.3.5 | High EPAP Alarm..... | 105 |
| 6.3.6 | Low EPAP Alarm..... | 105 |
| 6.3.7 | High V_{t_i} (High Inspired Tidal Volume Alarm)..... | 106 |
| 6.3.8 | Low V_{t_i} Alarm (Low Inspired Tidal Volume)..... | 106 |
| 6.3.9 | High MV_i Alarm (High Inspired Minute Volume Alarm)..... | 107 |
| 6.3.10 | Low MV_i (Low Inspired Minute Volume Alarm)..... | 107 |
| 6.3.11 | High V_{t_e} Alarm (High Expired Tidal Volume)..... | 108 |
| 6.3.12 | Low V_{t_e} Alarm (Low Expired Tidal Volume)..... | 108 |
| 6.3.13 | High MV_e (High Expired Minute Volume Alarm)..... | 109 |
| 6.3.14 | Low MV_e Alarm (Low Expired Minute Volume)..... | 109 |
| 6.3.15 | High Breath Rate Alarm..... | 110 |
| 6.3.16 | Low Breath Rate Alarm..... | 111 |
| 6.3.17 | Apnoea Alarm..... | 112 |
| 6.3.18 | Disconnection Alarm..... | 112 |
| 6.3.19 | Rebreathing Alarm..... | 114 |
| 6.3.20 | Obstruction Alarm..... | 115 |
| 6.3.21 | High FiO_2 Alarm..... | 115 |
| 6.3.22 | Low FiO_2 Alarm..... | 116 |
| 6.3.23 | High SpO_2 Alarm..... | 116 |
| 6.3.24 | Low SpO_2 Alarm..... | 117 |
| 6.3.25 | High $EtCO_2$ Alarm..... | 117 |
| 6.3.26 | Low $EtCO_2$ Alarm..... | 118 |

| | | |
|--------|---|-----|
| 6.3.27 | High InspCO ₂ Alarm (High Inspired CO ₂) | 118 |
| 6.3.28 | High Pulse Rate Alarm | 119 |
| 6.3.29 | Low Pulse Rate Alarm | 119 |
| 6.3.30 | PtcCO ₂ Alarm | 120 |
| 6.4 | Technical Alarms | 121 |
| 6.4.1 | Power Fail Alarm | 121 |
| 6.4.2 | High Patient Air Temp. (High Patient Air Temperature) | 121 |
| 6.4.3 | Low Patient Air Temp. (Low Patient Air Temperature Alarm) | 122 |
| 6.4.4 | Low Last Power Source Alarm | 122 |
| 6.4.5 | Crit. Low Last Power Source Alarm | 122 |
| 6.4.6 | Lost Mains Power Alarm | 123 |
| 6.4.7 | Exhalation Valve Control Error Alarm | 123 |
| 6.4.8 | SpO ₂ Disconnected (SpO ₂ Sensor Failure/Disconnection Alarm) | 124 |
| 6.4.9 | SpO ₂ Signal Lost Alarm | 124 |
| 6.4.10 | Poor SpO ₂ Signal | 124 |
| 6.4.11 | CO ₂ Disconnected (CO ₂ Sensor Failure/Disconnection Alarm) | 125 |
| 6.4.12 | CO ₂ Accuracy Error Alarm | 125 |
| 6.4.13 | Check CO ₂ Adapter Alarm | 125 |
| 6.4.14 | CO ₂ Sensor Error Alarm | 126 |
| 6.4.15 | FiO ₂ Disconnected (FiO ₂ Sensor Failure/Disconnection Alarm) | 126 |
| 6.4.16 | Ambient Pressure Compensation Lost Alarm | 126 |
| 6.4.17 | Temperature Comp. Lost (Ambient Temperature Compensation Lost Alarm) | 127 |
| 6.4.18 | Humidity Comp. Lost (Humidity Compensation Lost Alarm) | 127 |
| 6.4.19 | LED Failure Alarm | 127 |
| 6.4.20 | Low Alarm Battery Alarm | 128 |
| 6.4.21 | Alarm Battery Error Alarm | 128 |
| 6.4.22 | Internal/Click-In Battery Hot Alarm | 129 |
| 6.4.23 | Heated Circuit Temp. Alarm | 129 |
| 6.4.24 | High Humidifier Temp. Alarm | 130 |
| 6.4.25 | Humidifier Fault Alarm | 130 |
| 6.4.26 | Heated Circuit Fault Alarm | 131 |
| 6.4.27 | Internal Function Failure | 131 |
| 6.4.28 | Air Temp. Sensor Fail Alarm | 131 |
| 6.4.29 | Internal Error Alarm | 132 |
| 6.4.30 | Database Integrity Fail Alarm | 132 |
| 6.4.31 | Cooling Fan Error Alarm | 132 |
| 6.4.32 | Clock Failure Alarm | 133 |

| | | | |
|-----|--------|--|-----|
| | 6.4.33 | Internal Temp High Alarm | 133 |
| | 6.4.34 | Humidifier/Bypass Loose Alarm | 133 |
| 6.5 | | Alarm Test | 134 |
| | 6.5.1 | Alarm Signal Test | 134 |
| | 6.5.2 | Mandatory Alarm Tests | 134 |
| | 6.5.3 | Optional Alarm Tests | 137 |
| 7 | | Cleaning and Maintenance | 139 |
| | 7.1 | Cleaning the Nippy 4+ | 139 |
| | | 7.1.1 Main Unit | 140 |
| | | 7.1.2 Air Pathway Disinfection | 140 |
| | | 7.1.3 Patient Circuit..... | 140 |
| | 7.2 | Cleaning and Replacing the Air Filters | 141 |
| | | 7.2.1 Coarse Filter (Grey, Washable) | 141 |
| | | 7.2.2 Fine Filter (White, Disposable) | 142 |
| | 7.3 | Change of Patients..... | 142 |
| | 7.4 | Regular Maintenance | 142 |
| | 7.5 | Service and Repair..... | 143 |
| | 7.6 | Storage..... | 143 |
| | 7.7 | Disposal | 144 |
| 8 | | Technical Specifications | 145 |
| | 8.1 | System Description | 145 |
| | | 8.1.1 Pneumatic Diagram for the ventilator | 146 |
| | 8.2 | Data | 146 |
| | | 8.2.1 Worst Case Accuracy | 146 |
| | | 8.2.2 Settings Specifications..... | 146 |
| | | 8.2.3 Monitored Values Specifications..... | 150 |
| | | 8.2.4 Power Supply | 152 |
| | | 8.2.5 Environmental Conditions..... | 152 |
| | | 8.2.6 Other | 153 |
| | 8.3 | Emission and Immunity Declaration | 154 |
| | | 8.3.1 Nippy 4+ Essential Performance | 154 |
| | | 8.3.2 Guidance and Manufacturer's Declaration – Electromagnetic Immunity..... | 155 |
| | | 8.3.3 Guidance and Manufacturer's Declaration – Electromagnetic Emission..... | 158 |
| | | 8.3.4 Recommended separation distances between portable and mobile RF communications equipment and the ventilator | 158 |
| | | 8.3.5 Recommended separation distances between external power conductors and the ventilator..... | 159 |
| | | 8.3.6 Recommended separation distances between external power conductors and the ventilator..... | 160 |
| | 8.4 | Delivery Settings | 160 |
| 9 | | Accessories and Parts | 164 |

| | | |
|-----|-----------------------------|-----|
| 9.1 | Breas Accessories List..... | 164 |
| 10 | Patient Settings | 176 |
| 11 | FAA Compliance..... | 178 |

1 Introduction

WARNING!



Risk of Personal Injury

The Nippy 4+ must only be used:

- For the intended treatment in accordance with this manual and with the instructions given by the responsible clinical personnel.
- In accordance with the operating conditions specified in this manual.
- In original and unmodified shape and only with accessories specified or approved by Breas Medical.

Every other use may lead to risk of personal injury!

CAUTION!



Read this manual thoroughly so that you completely understand how the Nippy 4+ is operated and maintained before taking it into use, to ensure correct usage, maximum performance and serviceability. Non-professional caregivers (e.g. family members) should consult the medical equipment provider's respiratory therapist if they have any questions about the function, proper use, operation, service or maintenance of the Nippy 4+.



Breas Medical reserves the right to make changes to this product without any prior notification.

1.1 What is the Nippy 4+?

The Nippy 4+ is a pressure and volume ventilator capable of delivering continuous or intermittent ventilatory support for patients who require invasive or non-invasive mechanical ventilation. The Nippy 4+ is capable of running 24 hours/day.

The Nippy 4+ can be operated in the following combinations of ventilation and breath modes:

- Pressure Support (PSV)
- Pressure Support with TgV (PSV+TgV)
TgV= Target Volume
- Pressure Control (PCV)
- Pressure Control with TgV (PCV+TgV)
TgV= Target Volume
- Mouthpiece - Pressure (PCV-MPV)
- SIMV-Pressure (SIMV-P)

SIMV= Synchronized Intermittent Mandatory Ventilation

- Volume Control (VCV)
- Mouthpiece - Volume (VCV-MPV)
- SIMV-Volume (SIMV-V)

SIMV= Synchronized Intermittent Mandatory Ventilation

- CPAP

Compatible Patient Circuits

The Nippy 4+ can be used with a leakage circuit, an MPV circuit or a circuit with active exhalation valve. See 9 *Accessories and Parts*, page 164 for detailed information about approved patient circuits.

For leakage circuits: The patient circuit shall comply to ISO 17510. The leakage should be at least 12 l/min at 4 cmH₂O, to prevent rebreathing of exhaled air. The recommended leakage is 20 to 50 l/min at 10 cmH₂O pressure.

Compatible Patient Interfaces

For invasive use, the patient interface may be a tracheostomy tube (cuffed or uncuffed).

For non-invasive use it may be a mask, mouthpiece or pillow interface. See the patient interface's instructions for use when selecting the interface to use.

Data Log

The Nippy 4+ has an internal memory with a data log that holds the following data:

- Running hours
- Technical alarms
- Settings
- Asset data
- Treatment hours
- Treatment settings
- Device serial number
- Physiological alarms
- Detailed log, containing at least 24 h data of clinical data (monitored values)
- Breath log, containing at least 30 day data of (monitored values)
- Usage log (containing at least 1 year data of non-clinical events, alarms and settings)

The internal memory data is maintained also during power failure. The data can be transferred to a computer, printed out, and analysed via Breas software products.



For more information about Breas software products, please contact your Breas representative.

Expected Service Life

The expected service life of the Nippy 4+ is 8 years.

1.2 Intended Use

The Nippy 4+ ventilator (with or without the SpO₂ and CO₂ sensor) is intended to provide continuous or intermittent ventilatory support for the care of individuals who require mechanical ventilation. Specifically, the ventilator is applicable for paediatric through adult patients weighing more than 10 kg (22 lbs.)

The Nippy 4+ with the SpO₂ is intended to measure functional oxygen saturation of arterial hemoglobin (%SpO₂) and pulse rate.

The Nippy 4+ with the CO₂ sensor is intended to measure CO₂ in the inspiratory and expiratory gas.

The device is intended to be used in home, institution, hospitals and portable applications such as wheelchairs and gurneys. It may be used for both invasive and non-invasive ventilation. Nippy 4+ is not intended to be used as an emergency transport or critical care ventilator.

1.3 Contraindications

- The use of the Nippy 4+ is contraindicated for patients who need to be ventilated with oxygen concentrations (FiO₂) higher than achievable when combining inlet from a low pressure oxygen source at 15 l/min with actual ventilator settings.
- Generally, after surgery, the surgeon should be consulted to avoid organ damage and help determine ventilator parameters that do not adversely affect hemodynamics or have a negative impact on the patient's health status.
- In case of facial surgery, make sure to choose a suitable patient interface, in order to avoid discomfort and injury.
- The Nippy 4+ is not intended to be used as an emergency transport ventilator or critical care ventilator.

Undesirable Side Effects

If the patient experiences chest discomfort, pain, severe headache or shortness of breath while using the Nippy 4+, a physician or responsible clinician should be contacted immediately.

The following side effects may occur during the course of therapy with the Nippy 4+, patients are advised to report any new or changing adverse effects to their physician:

- Nasal, mouth or throat dryness
- Nosebleeds
- Abdominal bloating
- Ear or sinus discomfort
- Eye irritation
- Skin rashes

1.4 Intended Environment

The Nippy 4+ is intended to be used in clinical settings (e.g. hospitals, sub-acute care institutions), public spaces and home environments as well as during portable applications such as wheelchairs, personal family vehicles, ground ambulances and civil aircraft, excluding helicopters.

It is not intended for use during emergency transportation.

1.5 Intended Users

This section describes the intended users of the Nippy 4+, their qualifications and their related documents.

1.5.1 Respiratory Health Care Specialists

Health care professionals such as physicians and respiratory therapists, assigned to form the clinical authority when it comes to operating mechanical ventilators.

They have a good understanding of the human respiratory system and a general understanding of mechanical ventilators.

They are allowed to change the clinical settings of a ventilator and prescribe new settings. They may also operate software applications for follow-up on patient's ventilator treatment.

Training

The respiratory health care specialists shall be trained to a good knowledge of the Nippy 4+, its capabilities and the settings that can be made. This training consists of reading the Clinician's manual in full and it shall be conducted before operating the Nippy 4+.

Related Documents

The Clinician's manual is intended for the respiratory health care specialists. It shall be available for training of new personnel and as reference when operating the Nippy 4+.

1.5.2 Lay Operators

Day-to-day caregivers, patients, relatives and other non-professional users that operate the Nippy 4+ within the prescribed settings.

They are allowed to operate the Nippy 4+ with the Home mode activated. The lay operator may also perform basic maintenance that doesn't require special equipment or a service environment.

In Home mode, the device is locked in order to limit settings accessibility and hide features/controls.

The User Manual contains the information intended for patients and lay operators.

Training

The lay operator shall be trained to basic knowledge of the Nippy 4+ and in the specific operations they are assigned to perform. The training shall be based on the user manual and the responsible clinical personnel shall assess the level of training required for each lay operator.

Related Documents

The User manual is intended for lay operator. It shall be available for the training and as reference when operating the Nippy 4+.

1.5.3 Service Personnel

Certified service personnel with responsibility to maintain the equipment in proper working order. They have a technical education and/or relevant experience of technical work on electrical equipment. If local or national regulations requests additional authorization or competence, these shall be complied to.

Certified service personnel may perform any repairs, upgrades or service operations that they have been certified to perform, as long as they have the required equipment and the operation is performed in an appropriate environment. They may also operate software applications for follow-up on ventilators usage and for troubleshooting.

Training and Certification

Service personnel shall be trained on the Nippy 4+ and certified by Breas for being allowed to perform any service, repairs or other operations on the Nippy 4+. The training consists of reading the services manual in complete. After completed training the certification test may be performed.

Related Documents

- The Service Manual.
- The Clinician's Manual.
- Service bulletins, available for certified service personnel on the Breas extranet.

1.6 About this Manual



CAUTION!

Always read this manual before setting up and using the Nippy 4+ or performing maintenance on the machine, to ensure correct usage, maximum performance and serviceability.

1.6.1 Audience









This manual is intended for patients and other lay users operating the Nippy 4+.



- Care providers, clinical personnel, physicians and others who require a working knowledge of the Nippy 4+ will find additional information on settings and functions in the Clinician's Manual.
- Service personnel may order the Service Manual that contains detailed technical information for maintenance, service, repair and disposal procedure.

1.6.2 Icons in this Manual

In this manual, icons are used to highlight specific information. The meaning of each icon is explained in the table below.

| Icon | Explanation |
|---|---|
|  | Warning! Risk of death or personal injury. |
|  | Warning! Risk of Cross-contamination. |
|  | Warning! Risk of electric shock. |
|  | Warning! Hot surface, risk of burns. |
|  | Warning! Flammable material, risk of fire. |
|  | Caution! Risk of equipment damage, loss of data, extra work, or unexpected results. |
|  | Note Information that may be valuable but is not of critical importance, tips. |
|  | Reference Reference to other manuals with additional information on a specific topic. |

1.7 Manufacturer Information

Legal Manufacturer



Postal Address

Breas Medical AB
Företagsvägen 1
SE-435 33 Mölnlycke Sweden

Web Address

www.breas.com

Email address

breas@breas.com

Phone

+46 31 868800, Order: +46 31 868820, Technical support: +46 31 868860

Fax

+46 31 868810

UK Sales and Support Information

Phone UK Head Office

01789293460

Fax UK Head Office

01789262470

General Enquiries and Ordering

orders@nippyventilator.com

Tech Support

techsupport@nippyventilator.com

2 Safety Information

2.1 General User Precautions

WARNING!



Risk of Personal Injury

The Nippy 4+ must only be used:

- For the intended treatment in accordance with this manual and with the instructions given by the responsible clinical personnel.
- In accordance with the operating conditions specified in this manual.
- In original and unmodified shape and only with accessories specified or approved by Breas Medical.

Every other use may lead to risk of personal injury!



Risk of Insufficient Ventilation

Usage outside the specified operating conditions may cause reduced performance. The Nippy 4+ must only be used in accordance with the operating conditions specified in this manual



Risk of Reduced Safety and Performance

Accessories that have not been tested with the Nippy 4+ might affect safety features and performance negatively. Only use the Nippy 4+ with accessories approved by Breas Medical.

Incompatible parts can result in degraded performance and change of pressure gradient.

If unapproved accessories are used, Breas Medical has no responsibility for the safe and effective use of the Nippy 4+.

The responsible organization is responsible for ensuring the compatibility of the ventilator and all parts used to connect to the patient before use.



Changes to the patient circuit, like adding or removing accessories or changing type or length of breathing tube, may affect both circuit compliance and alarm triggering conditions.

It's recommended to perform a pre-use test and re-test the alarm function after making changes to the patient circuit.



When a patient is treated, there must be a supervising person present during the treatment in order to take care of alarms and conditions that the patient cannot solve on their own.



Always have immediate access to an alternative means of ventilation, which is ready for use, to avoid patient death or serious injury. Failure to have an alternate means of ventilation can result in serious injury or patient death if ventilator fails.



To prevent death or serious injury, monitor the patient and the ventilator regularly in order to determine the need to provide emergency ventilation when an alarm sounds or the ventilator malfunctions.



If using the ventilator in a transport case, only use the protective cover listed as approved accessory to prevent adverse ventilator performance, which can consequently result in patient death.



Measured values for volume and expired CO₂ may differ from the actual patient values due to unintentional leaks.



Do not obstruct the gas intake port.



To prevent disconnection of the patient circuit during use, only tubes in compliance with ISO 5367 or ISO 80601-2-74 should be used.



The ventilator may not work properly if any part has been dropped, damaged or submerged in water.

WARNING!



Risk of Burns

Covering breathing tubes with a blanket or heating them with an overhead heater can affect the quality of the therapy or injure the patient.

WARNING!



Risk of Faulty Treatment

If the patient is admitted to a hospital or is prescribed any other form of medical treatment, always inform the medical staff that the patient is on mechanical ventilation treatment.



Risk of Faulty Treatment

Do not use the Nippy 4+ in the event of:

- Suspected damage to the device, including the occurrence of Internal Functional Failure alarms.
- Unexpected patient symptoms during treatment.
- Unexplainable or sudden changes of pressure, performance or sound during operation.
- Delivered air being abnormally hot or emitting an odor.

Contact your responsible care provider for an inspection.



Risk of Faulty Treatment

The responsible organization should periodically reassess the settings of the therapy for effectiveness.



The ventilator therapy settings must always be based on medical advice and must be carried out by authorised clinical personnel only. When changing treatment settings or changing to another device a clinical assessment must be performed to determine if blood gas measurement needs to be performed.



Before starting treatment, always perform the procedure 4.5 *Inspecting the Nippy 4+ before Use*, page 45.



Risk of Unnoticed Critical Conditions

- The alarm sound level should be set to a clearly audible level. Setting the alarm sound level below that of the ambient sound level can impede recognition of alarm conditions.

CAUTION!



Clinical personnel must read the Clinician's manual thoroughly and understand the ventilator operation before setting up and using the ventilator.



- Handle the ventilator with care.
- Do not use the ventilator while in the carry bag.
- Do not use the ventilator with nitric oxide, helium or helium mixtures.



Contact Injuries: Skin irritation may occur due to prolonged exposure to either a mask (if used) or the SpO₂ module.



Ensure that the cooling air intake vents are not blocked. If the vents are blocked, especially in hot use environments, the surface temperature of the patient circuit may rise above 41°C (106°F). In a 40°C (104°F) environment and with the cooling air intake vents blocked, surface temperatures as high as 50°C (122°F) can occur. Before an unsafe temperature is reached, the “High Patient Air Temp” alarm will occur. If this alarm occurs, assure that the ventilator cooling air intake path is free of obstruction and that the patient circuit surface is not heating the patient’s skin.

WARNING!



Risk of Electric Shock

Modifying or using the ventilator with accessories that are not specified or approved by Breas may cause cardiac arrhythmia.

The Nippy 4+ must only be used in original and unmodified shape and only with accessories specified or approved by Breas Medical.

Inadequate use of device or accessories may cause loss of treatment or decreased performance.

NOTE



Any serious incident that has occurred in relation to this device should be reported to the competent authority and the manufacturer.

2.1.1 Requirements for Life Supporting Treatment

Life supporting-supporting treatment requires that:

- An emergency equipment (e.g. resuscitation bag) is available.
- One of the following means of surveillance is used:
 - Using the EtCO₂ sensor accessory, or an external EtCO₂ monitor (capnometer). This surveillance method can be used for ventilation with active exhalation valve circuits as well as leakage port circuits.

The CO₂ sensor must be connected between the patient and the exhalation valve or leakage port to be able to measure exhaled gases. If using an external CO₂ monitor, it shall fulfil the ISO 80601-2-55 standard (Medical electrical equipment – Particular requirements for the basic safety and essential performance of respiratory gas monitors).

- Supervising the ventilator’s monitoring of exhaled volume. This surveillance method can be used for ventilation with leakage circuits only.

2.2 Electrical Safety



WARNING!

Risk of Electric Shock

High voltage contact may cause cardiac arrhythmia.

- Do not operate the Nippy 4+ if it has a damaged power cord, power supply or casing.
- To avoid electrical shock, only clean the Nippy 4+ according to instructions in this manual. Do not soak or immerse the Nippy 4+ into any fluids.
- Use the approved power supply units only.

Use of unapproved power supply units may compromise the electrical isolation and lead to risk of electric shock.

- Do not use more than one multiple portable socket-outlet or extension cord.
If a multiple portable socket-outlet is used, it must not be placed on the floor.
- The operator must not touch accessible contacts of connectors and the patient simultaneously.
- Nurse Call must only be connected to a safety extra low voltage system with an isolation from mains voltage which complies with the requirements of IEC 60601-1.



WARNING!

Risk of Faulty Treatment

Electromagnetic Interference may cause electrical equipment to malfunction.

- The aspects of electromagnetic compatibility must be considered.
 - The Nippy 4+ should not be used adjacent to or stacked with other equipment; if adjacent or stacked use is necessary, the Nippy 4+ should be observed to verify normal operation in that configuration.
 - Mobile or transportable radio transmitters may interfere with the Nippy 4+.
 - Further guidance for safe installation of the ventilator can be found in the chapter about emission and immunity declaration.
- If a portable AC power supply is used, make sure that the voltage variations are within the operating limits of the Nippy 4+.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Nippy 4+, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING!



Avoid touching the contacts within the ventilator click-in battery compartment. Under certain circumstances touch current limits per IEC 60601-1 may be exceeded.

2.3 Environmental Conditions

WARNING!



Risk of Intoxication

Do not use the Nippy 4+ in a toxic environment.

WARNING!



Risk of Fire

Do not use the Nippy 4+ in environments where explosive gases or flammable anesthetic agents present.

WARNING!



The delivered patient air can be as much as 4°C (7°F) higher than ambient temperature. Caution should be exercised if the room temperature is greater than 36°C (97°F).



Risk of Faulty Treatment

If a room humidifier is used, place it at least 2 meters away from the Nippy 4+.



Risk of faulty Treatment

The performance of the Nippy 4+ may deteriorate at altitudes or ambient temperatures outside the operation conditions specified in the chapter *Technical Specifications*.

- Do not use the ventilator while positioned in a warm place, such as direct sunlight or close to a radiator as this might lead to temperature outside the specifications.
- Do not use the ventilator in an hyperbaric chamber, as this would cause an ambient pressure outside the specifications.
- Do not use the ventilator immediately after storage or transport outside the recommended operating conditions.



Risk of Faulty Treatment



Do not use or store the Nippy 4+ in a magnetic resonance (MR) environment. Use of the Nippy 4+ in an MR environment may result in malfunction of the Nippy 4+ and pose unacceptable risk to the patient, medical staff or other persons.



Unsteady indicated values for delivered volumes or pressures and the occurrence of alarm conditions without apparent cause may be an indication of loss of performance due to electromagnetic disturbances. Follow the instructions above and the guidance provided in 8.3 *Emission and Immunity Declaration*, page 154 to mitigate the effects of electromagnetic disturbances.



CAUTION!

The ventilator, any accessories and all replaced parts, must be disposed of in accordance with the local environmental regulations regarding the disposal of used equipment and waste.

2.4 Usage of Patient Circuit



WARNING!

The ventilator supports leakage circuits, circuits with an active exhalation valve and circuits with mouthpiece interface. Failure to use a mask or accessory that minimizes rebreathing of carbon dioxide or permits spontaneous breathing can cause asphyxiation.



For the ventilator to deliver treatment according to settings, it is important that the selection of the patient circuit type is correctly set.



Risk of Insufficient Ventilation

Insufficient ventilation may cause transient hypoxia.

Before use:

- Make sure that the patient circuit and joined parts are undamaged and correctly connected, in order to avoid unwanted leakage.
- For leakage circuits: Make sure that the leakage port of the circuit or mask is not blocked or obstructed. This port prevents rebreathing by flushing the exhaled air.
- For active exhalation valve circuits: Check the function of the exhalation valve and that is not blocked or obstructed.
- The Nippy 4+ should be turned on and the function of the leakage port should be checked before use: The pressurized air from the Nippy 4+ causes a continuous flow of air through the leakage port, enabling flushing of exhaled air.



When the patient circuit is replaced or modified, check the alarm settings as changes to the patient circuit may affect the alarm triggering. Also consider performing a pre-use test for optimizing the therapy.



Risk of Suffocation

If the patient needs assistance to remove the patient interface, the patient shall not be left alone. This is to avoid the risk of re-breathing of CO₂ in case of accidental ventilator failure.

Do not breathe through the connected patient circuit unless the ventilator is turned on and operating properly.

WARNING!



Risk of Electric Shock

Do not use antistatic or electrically conductive hoses or tubing with the ventilator breathing system. This could result in electrical shock.

WARNING!



Patient connected parts and all filters must be replaced regularly to ensure correct function of the ventilator. All replaced parts must be disposed of according to local environmental regulations regarding the disposal of used equipment and parts.



By conducting a pre-use test (see 4.7 *Performing the Pre-use Test*, page 47) the compatibility of the complete patient circuit configuration with the ventilator can be verified. If a pre-use test is successfully performed the circuit configuration meets the required characteristics.



Risk of Suffocation

Periodically check for moisture in the patient circuit. When present, remove the moisture. Before attempting to dry the circuit, disconnect it from the Nippy 4+ to ensure no water flows back into the Nippy 4+. The frequency at which these checks must be performed will depend on the patient's condition and the device used. The responsible caregiver should assess this on an individual basis in accordance with the patient's needs.



Risk of Insufficient Ventilation

Insufficient ventilation may cause transient hypoxia. The use of equipment such as endotracheal tubes, oral/nasal tubes, adaptors etc. with small inner diameters or high resistance filters, humidifiers etc. increase the resistance in the patient circuit which may interfere with the operation of the patient disconnect function. It may also interfere with the device trigger function. This impact can be reduced by conducting a pre-use test.



Make sure that the exhalation valve or the leakage port never is blocked or obstructed.



Risk of Constriction

Entanglement with cables or tubing constricting airways may cause asphyxiation. Do not leave long lengths of air tubing or cables around the top of the bed. They could twist around the patient's head or neck while sleeping.



The ventilator is equipped with a rebreathing alarm. The alarm is not a substitute for operator vigilance in ensuring that the leakage port or exhalation valve remains clear at all times. Periodically check the patient circuit during therapy.



In general, as pressure decreases, the potential of rebreathing increases. Lower pressures produce less flow through the leakage port which may not clear all CO₂ from the circuit to prevent rebreathing.



Risk of Excessive Carbon Dioxide

Insufficient carbon dioxide removal may cause arterial acidemia. For reducing the risk of rebreathing CO₂, make sure that the leakage port is located as near the patient interface as possible. This is even more important for treatments with low pressure, as this reduces the flow through the leakage port.



WARNING!

Risk of Cross-Contamination

Patient circuits might get contaminated by exhaled gases. To avoid cross-contamination, always use a properly cleaned or a new patient circuit when the Nippy 4+ is to be used by a new patient.

NOTE



For masks and accessories, always follow the manufacturer's instructions.

2.5 Usage of Filters

WARNING!



Always use the ventilator with patient air inlet filters installed. Only use the ventilator with accessories recommended by Breas Medical.



Risk of Insufficient Ventilation

Insufficient ventilation may cause transient hypoxia.

Replace or clean the inlet filters as specified in the *Maintenance* chapter.

Using old or clogged filters may cause the Nippy 4+ to operate at higher temperatures than intended.

When operating the Nippy 4+, make sure that the air inlet and filters are not obstructed or occluded.



Risk of Insufficient Ventilation

Insufficient ventilation may cause transient hypoxia.

Do not use high resistance bacteria filter at the air outlet of the Nippy 4+. High resistance bacteria filters placed between the air outlet and the patient interface may interfere with the operation of the patient disconnect function. It may also interfere with the device trigger function.

WARNING!



Risk of Cross-Contamination

Deep tissue or mucosal contact with infectious agents may cause infections.

If the Nippy 4+ is used by several patients, a low resistance bacteria filter shall be used between the air outlet and the patient circuit, for preventing patient cross-contamination. Reuse of bacteria filter, patient circuit or mask may expose patients to contagious agents.

2.6 Humidification

WARNING!



When adding or removing an HME (Heat and Moisture Exchanger, artificial nose) or HCH (Hygroscopic Condenser Humidifier), always reassess the settings, including alarm settings, and perform a pre-use test.

**Risk of Suffocation**

When the attachable humidifier is installed, the Nippy 4+ must be located below the patient and on a flat surface. This is to prevent personal injury from accidental spillage or from excess water or condensation flowing down the patient tube and into the patient's mask. Extra cautions should be taken for patients who are unable to guard their airways or cannot pull off the mask.



When using the humidifier or a nebuliser any patient air filter will need more frequent replacement to prevent increased resistance or blockage.



The ventilator accuracy can be adversely affected by the gas added by the use of a pneumatic nebuliser.

WARNING!**Risk of electric shock**

If using the protective cover or the carry bag, first remove the attachable water chamber. Water spillage may cause electric shocks.

WARNING!

The use of a heated wire patient circuit decreases condensation in the patient circuit.



In case of invasive application, the use of an appropriate external heated humidifier or HME is recommended.



If the condensation in the patient circuit is excessive, the use of a heated humidifier may require the installation of a water trap in the circuit. The water trap prevents any condensed water in the patient circuit from running into the patient airways and causing personal injury.



Any external humidifier connected to the ventilator must comply with ISO 8185 or 80601-2-74.



Any HME connected to the ventilator must comply with ISO 9360.



Do not add any attachments or accessories to the humidifier that are not listed in the instruction for use of the humidifier or the humidifier might not function correctly affecting the quality of the therapy or injuring the patient.



The use of an HME or an external humidifier may require readjustment of the ventilator low-pressure alarm.



Certain HMEs and HCHs are sufficient to provide humidification when the ventilator is used invasively. Check specific suppliers' recommended use.

NOTE



The ventilator has been tested and validated with the Fisher & Paykel MR850 heated humidifier.

2.7 Cleaning and Maintenance

WARNING!



Risk of Electric Shock

Cleaning with excessive water or opening the device's casing without certified training may cause electric shocks.

The Nippy 4+ should be regularly cleaned and maintained in accordance with this operating manual.

WARNING!



Risk of Faulty Treatment

Service and Maintenance of the Nippy 4+ shall not be performed when the Nippy 4+ is in use.

WARNING!



Risk of Electric Shock

High voltage contact may cause cardiac arrhythmia.

Repairs and modifications must be carried out by authorized technicians only and in accordance with instructions from Breas Medical

- The Nippy 4+ must not be opened, repaired or modified by unauthorized personnel. If subjected to unauthorized operations, Breas Medical is no longer responsible for the device's performance and safety and all warranties will become invalid.
- The Nippy 4+ must not be modified or interconnected to unapproved equipment.

CAUTION!



Do not attempt to autoclave or sterilize the Nippy 4+.

2.8 Usage of Oxygen

When using the Nippy 4+ with oxygen, always follow the oxygen provider's instructions and use only medical grade oxygen complying with local regulations.

WARNING!



As this medical device uses an alternative small-bore connector design different from those specified in the ISO80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need to be taken by the user to mitigate these reasonable foreseeable risks.

WARNING!



Risk of fire

The presence of oxygen can speed up combustion of inflammable materials.



Risk of Fire

When oxygen is used with the Nippy 4+, the oxygen flow must be turned off when the Nippy 4+ is not operating. Oxygen delivered into the patient tubing may accumulate within the machine enclosure. Oxygen accumulated in the machine enclosure increases the risk of fire.

WARNING!



Do not use a humidifier between the oxygen source and the ventilator, in order to humidify the oxygen flow.
If humidification is required, use the attachable humidifier or an external humidifier after the patient air outlet.

WARNING!



Risk of Fire

Ventilate the room adequately. Do not smoke in a room where oxygen is being used.



Risk of Fire

Naked light bulbs and other sources of ignition must be kept a minimum of 2 meters (6 feet) away from the oxygen cylinder or any part of the patient circuit.



Risk of Fire

Do not use aerosols or solvents close to the oxygen supply, even when the oxygen supply is shut off.

WARNING!



Supplemental oxygen with a flow up to 15 l/min can be added by an oxygen source with rotameter such as oxygen cylinder, central oxygen supply system or an oxygen concentrator.



Risk of faulty Treatment

At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure delivered, the patient's breathing pattern, mask selection and leak rate.

To monitor the oxygen concentration, use the FiO₂ sensor accessory.



Supplemental oxygen flow and pressure must not exceed 15 l/min and 100 kPa.

CAUTION!



Supplemental oxygen is added before the volume measurement sensor and thereby included in the measurements. However, the oxygen concentration still has influence on the volume measurement of the delivered air.

This measurement is based on a normal oxygen concentration of 21%. If the oxygen concentration is higher, the actual inspired volume will deviate from the monitored volume as follows:

- 40% oxygen concentration: -2.5% deviation
- 60% oxygen concentration: -5% deviation
- 80% oxygen concentration: -7.5% deviation

3 Product Description

3.1 Main Components

This section describes the components of the Nippy 4+ medical electric equipment.




NOTE

- There might be local variations of the main components configuration.
- The standard Nippy 4+ and its packaging do not contain any natural rubber latex.

Carry bag

Function: Storage for transportation

Part No: 006014




Manual

Function: Product and usage information

Part No:

User's manual: 007303


Clinician's manual: 007302



Patient Circuit

Function: Delivers air to the patient (applied part)

Delivered patient circuit depends on the sales configuration. See 9 *Accessories and Parts*, page 164 for approved patient circuits.




Patient air inlet filter, fine, white, disposable

Function: Fine inlet air filtration.

Material: AS 100

NaCl Penetration: (0.65 μm NaCl @ 95 l/min) = <7.35%


Part No: 007103 (5pcs)



Power Supply

Function: Deliver power to the ventilator

Part No: 006396




Power cord

Function: Deliver power to the AC power supply


Part No:

UK: 003521




Nippy 4+ Main Unit

Main Unit




3.2 Front View



| No | Item | Function |
|-----|--------------------------|--|
| 1 | Alarms button | Access to alarm settings |
| 2 | Sensor | Ambient light sensor |
| 3-4 | Alarm (red & yellow) LED | Alarm indication: Red = High priority Yellow = Medium priority |
| 5 | Audio pause LED | Paused alarm sound indication |
| 6-8 | Plus, Set, Minus buttons | Function according to display Plus = Increase, go up Set = Enter / Navigation Minus = Decrease, go down |
| 9 | Mains LED | Power source indication: Mains |
| 10 | External DC LED | Power source indication: External DC |
| 11 | Click-in battery LED | Power source indication: Click-in battery |
| 12 | Internal battery LED | Power source indication: Internal battery |
| 13 | Menu button | Menu/Navigation |
| 14 | Settings, Mode buttons | Select settings, modes and profiles |


3.3 Side Views

Click-in compartment side



| No | Item | Function |
|----|----------------------|---|
| 1 | Audio pause | Pause the alarm sound |
| 2 | Start/Stop | Start/Stop ventilation treatment |
| 3 | Patient air outlet | Connection for patient circuit |
| 4 | Air bypass unit | Click-in airway/silencer for use without the click-in humidifier. (If the click-in humidifier is used, it replaces the air bypass unit) |
| 5 | Click-in compartment | Compartment for either of the accessories click-in humidifier or click-in battery. |
| 6 | Side panel | Cover |
| 7 | Carrying handle | Handle for lifting the ventilator |
| 8 | Cooling air outlet | Outlet internal cooling |
| 9 | Cooling air inlet | Inlet internal cooling |

Filter Side




| No | Item | Function |
|----|-------------------------------|--|
| 1 | Side panel | Cover |
| 2 | Memory card slot (SD card) | Memory download |
| 3 | Alarm beeper | Alarm Sounds Output |
| 4 | Patient air inlet | Air bypass unit in, replaceable filters |
| 5 | Internal battery | Compartment for the internal battery |
| 6 | FiO ₂ sensor hatch | Compartment for the optional FiO ₂ sensor |

3.3.1 Detaching and Reattaching the Side Panels


Detaching the Filter Side Panel

- 1 Lift the handle to access the release button (A).
- 2 Looking from behind, to dismount the filter side panel press the button above the panel (B). The panel is released.
- 3 Remove the panel. (C)




Reattaching the Filter Side Panel

- 1 Lift the handle to access the release button (A).
- 2 To mount the filter side panel, insert the tabs (B) on the lower side of the panel into the holes (C).
- 3 Press the side panel into the casing until it clicks in place at the button (D).




Detaching the Click-in Compartment Side Panel

- 1 Lift the handle to access the release button (A).
- 2 Press the button above the panel (B).
⇒The panel is released.
- 3 Remove the panel (C).



















Reattaching the Click-in Compartment Side Panel

- 1 Lift the handle to access the release button (A).
- 2 To mount the click-in compartment side panel, insert the tabs (B) on the lower side of the panel into the holes (C).
- 3 Press the side panel into the casing until it clicks in place at the button (D).















3.4 Equipment Designation









| No | Item/Symbol | Description | Colour |
|----|---|--|---|
| 1 |  | Electrical connector, for power to the heated patient circuit. | |
| 2 |  | Exhalation valve port | |
| 3 |  | Connection for low pressure/bleed-in oxygen source | |
| 4 |  | CO ₂ interface port |  |
| 5 |  | SpO ₂ interface port |  |
| 6 |  | Remote start/stop, Audio pause, and effort belt interface port |  |
| 7 |  | Remote alarm and Nurse call interface port |  |
| 8 |  | Network connection port | |
| 9 |  | USB data connection port | |
| 10 |  | Mains/External DC inlet |  |

3.4.1 Additional Symbols

This section describes symbols and markings that might appear on the parts or packaging of the Nippy 4+.

| Symbol | Description |
|---|--|
|  | Internal battery |
|  | Product number |
|  | Read user instructions. |
|  | Attention: Read the intended use in the manual. Read the Safety chapter in the manual. |
|  | This product must not be exposed to open fire. |
|  | This product should be recycled. |
|  | Read 7.7 <i>Disposal</i> , page 144 for information about recycling and disposal. |
| IP22 | Degree of protection provided by enclosure: IP22. See 8.2.5 <i>Environmental Conditions</i> , page 152 for detailed information. |
|  | Manufacturer |
|  | Serial number |
|  | This product is a Medical Device. |
|  | Date of Manufacture |
|  | IEC protection Class II: Double insulated equipment. |

| Symbol | Description |
|---|--|
|  | Indication of applied parts (IEC 60601-1 Type BF, Isolated Applied Part) |
| Rx Only | (Symbol only applicable in U.S.) Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. |
|  | Meets all requirements for CE marking according to relevant European health, safety and environmental protection legislation |
|  | Meets all requirements for CE marking according to relevant European health, safety and environmental protection legislation |
|  | Do not obstruct air inlets or outlet. |
|  | Single patient use |
|  | Hot Surface. Do not touch. (Heating plate in click-in compartment.) |

4 Preparing the Nippy 4+ for Use




WARNING!

Read 2 *Safety Information*, page 17 before setting up the Nippy 4+.


4.1 Checking the Nippy 4+ before First Use

When using the Nippy 4+ for the first time, follow the instructions below:

- 1 Check that all main components and ordered accessories have been delivered (refer to the packing note or the invoice, if available).



- 2 Ensure that the equipment is in good condition.
- 3 If stored more than 1 month, connect the Nippy 4+ to the power supply to recharge the internal battery.
- 4 Check that the grey and white air filters are installed.



4.2 Placing the Nippy 4+




WARNING!

Read 2.3 *Environmental Conditions*, page 22 carefully to make sure all conditions are met and considered.


- 1 Place the Nippy 4+ on a solid, flat, and clean surface.

The Nippy 4+ should be placed lower than the patient in order to prevent the device from falling on the patient, as well as preventing condensed water from reaching the patient.



Overnight, the Nippy 4+ should be placed close enough to the patient's bedside to allow movements during the sleep without pulling the Nippy 4+ of its surface.



- 2 Make sure that nothing can block the patient air inlet.



- 3 Make sure that nothing can block the cooling air inlet or outlet.




CAUTION!



- Do not place the Nippy 4+ on a soft surface that will prevent the air flow underneath the device.

Never cover the device.



- Always position the Nippy 4+ so the power supply lays on a surface without strain to the power cord. The power supply shall be easy to disconnect, if required to isolate the Nippy 4+ from the mains.


4.3 Connecting the Nippy 4+ to Mains

WARNING!




Read 2.2 *Electrical Safety*, page 21 carefully to make sure all conditions are considered and met.


- 1 Plug the power supply into the power inlet of the Nippy 4+.



- 2 Make sure a clicking sound is heard to ensure the power supply is completely inserted.



- 3 Connect the power supply's power cord to the mains supply.



To isolate the Nippy 4+ from the mains supply, disconnect the power supply.

4.4 Connecting the Patient Circuit



WARNING!


Read 2.4 *Usage of Patient Circuit*, page 23 carefully to make sure all conditions are considered and met.

The Nippy 4+ can be used with the following circuits:


- Single limb circuit with external active exhalation valve
- Single limb circuit with external leakage port
- Single limb circuit connected to patient interface with integrated leakage port
- Circuit with mouthpiece interface

Connect the Patient Circuit


- 1 Inspect the circuit for damages.
- 2 Connect the patient circuit to the patient air outlet on the ventilator.



- 3 If having a heated patient circuit, connect the heated wire electrical plug to the socket on the ventilator.



- 4 If having an active exhalation valve patient circuit, connect the pilot pressure tube at the back of the ventilator.



- 5 Connect the other end of the patient circuit to the patient interface.
If having a leakage circuit and a patient interface without integrated leakage port: Make sure to use a leakage port between the circuit and the patient interface.
- 6 If using a leakage circuit or an active exhalation circuit, make the correct setting for the type of circuit. See .
For MPV circuits, the type of circuit is set when activating the MPV mode.

4.5 Inspecting the Nippy 4+ before Use

Inspection of Device

- Check that there is no visible damage.
- Check that the surface is clean.

Inspection of Cables

- Check that all cables are recommended by Breas.
- Check that the cables are undamaged.
- Check that the cables are properly connected.

Inspection of Placement

- The Nippy 4+ shall be placed on solid flat surface below the patient level (see 4.2 *Placing the Nippy 4+*, page 41).
- Make sure that nothing can block the air inlet at the side.

Inspection before Use

1. Connect a patient circuit to the Nippy 4+.
2. Press the “Start/Stop” button to start Nippy 4+.
3. Ensure that the treatment settings and alarm settings are adjusted as prescribed.
4. If needed, perform a pre-use test. (The Nippy 4+ default settings are to prompt for a pre-use test when powering up.)
5. Press and hold the Start/Stop button until the progress bar is filled.
6. Check that the alarm LEDs flash and a is heard when the Nippy 4+ is started. This is the alarm signal test, if there is no signal, do not use the Nippy 4+ and contact your service provider.
7. Disconnect the power supply for more than 5 seconds. Check that the device switches to the internal battery (or click-in battery if connected) and that a *Lost mains power alarm* is shown on the screen together with an audible warning. If this is not the case contact your service provider.
8. Reconnect the power supply. Check that the device switches to the mains supply and that an information message is shown on the screen together with an audible signal.

4.6 Adjusting the Nippy 4+ Patient Settings



WARNING!

The configuration of the Nippy 4+ therapy settings must always be prescribed by a licensed physician and carried out by an authorised health care professional.

For detailed information about the treatment parameters of the Nippy 4+, see 5.6 *Treatment Settings*, page 53.

Follow the instructions below when setting up the Nippy 4+:

- Adjust the settings to find the best possible breathing comfort for each patient.
- If you have changed the ventilation mode, press **Select** and review the settings before pressing **Confirm**.
- Always document the patient settings.
- The ventilator always starts in the mode and with the settings that were active when it was switched off.

4.7 Performing the Pre-use Test

The pre-use test is used for detecting the type and characteristics of the patient circuit that is connected to the Nippy 4+. The resistance and compliance of the patient circuit are measured and calculated. This will be used to compensate for pressure drop in the patient circuit and the compliance of the patient circuit.

The patient shall not be connected during the pre-use test.

NOTE



If the pre-use test has not been performed, the Nippy 4+ will operate with default patient circuit compensation.

Starting the Pre-use Test Manually

- 1 On the **Menu**, select **Pre-use Test** and then **START PRE-USE TEST**.

Activating the Pre-use Test Prompt

- 1 On the **Menu**, select **Pre-use Test** and press the **Set** button.
- 2 Using the + and - buttons, set **Start Pre-use Test** to **On**.
- 3 Press the **Set** button to confirm the setting.

Pre-use Test Sequence

When performing a pre-use test, the instructions on the display will guide you through this sequence:

| Step | Action |
|------|--|
| 1 | Start of pre-use test. |
| 2 | Connect the patient circuit. |
| 3 | Make sure that nothing is blocking the patient end of the circuit. |
| 4 | Wait while the Nippy 4+ is checking the patient circuit resistance. If the resistance is not within the limits, the test will end without performing the following steps. The result will be displayed for review. |
| 5 | Block the end of the patient circuit with an air tight object. |
| 6 | Wait while the Nippy 4+ is checking the patient circuit compliance and leakage. |
| 7 | Test finished. Review the test result. |

5 How to Use the Nippy 4+

WARNING!





Read 2 *Safety Information*, page 17 before using the ventilator. When the ventilator is handed over to the patient, the physician in charge or hospital staff must instruct the patient in how the unit works.

5.1 Switch the Nippy 4+ On and Off


5.1.1 Switch On and Enter Operating Mode

- 1 Make sure the power supply is properly connected.
- 2 To start treatment and enter operating mode first press and hold the Start/Stop button on the Nippy 4+.
- 3 Release the Start/Stop button when the progress bar is filled.
- 4 Select Yes/No if asked to “Perform Pre-use Test”.




5.1.2 Stop Treatment


1 To stop treatment and enter standby mode, first press and hold the Start/Stop button on the front panel.



2 Release the Start/Stop button when the progress bar is filled.



3 Press "Yes" to stop the treatment.




5.1.3 Turn off / Enter Sleep Mode


If the Nippy 4+ is running on batteries, this procedure turns it off.

If the Nippy 4+ is connected to mains, this procedure puts it in sleep mode (all functions are off, except battery charging)

- 1 When the Nippy 4+ is in standby mode (no treatment is running), press the Start/Stop button on the front panel.
- 2 When asked to confirm the action, press "Yes".










5.2 Navigating the Display



| No. | Item |
|-----|--|
| 1 | Mode / Profile Use this button to select between user profiles (if profiles are configured by your clinician). |
| 2 | Treatment settings buttons Use these buttons to change the settings in their respective frame on the display. See 5.6 <i>Treatment Settings</i> , page 53 for more information. |
| 3 | Accessory/Function icons Indicates connected or activated accessories or functions. |
| 4 | Alarms settings button Use this button to change the alarm settings. |
| 5 | + button Use this button to increase a value when editing a setting, or to move up in the menu or on the alarm settings list. |
| 6 | Set button Use this button to select a setting to edit and to confirm a change. |
| 7 | - button Use this button to decrease a value when editing a setting, or to move down in the menu or on the alarm settings list. |
| 8 | Menu/More button Use this button to open the menu for accessing information and settings for the device and for comfort functions. |
| 9 | Alarms list Displays alarm settings. See 5.7 <i>The Alarms List</i> , page 63 for more information. To see an alarms history list, open to the Alarm/Event log from the Menu. |

| No. | Item |
|-----|--|
| 10 | Monitored Values pane. Displays read outs of monitored values and bar graphs of the current Pressure, Flow and Volume during treatment. Values from connected accessories are displayed with yellow text. |
| 11 | Home mode lock Indicates whether the settings are locked to Home mode. |
| 12 | Battery Icon Displays battery charge status by colour and percentage of fully charge. If having the click-in battery, it will have an icon of its own. |
| 13 | Settings list Displays treatment mode and treatment settings. See 5.6 <i>Treatment Settings</i> , page 53 for more information. |


5.3 Symbols Used on the Display

| Symbol | Description |
|---|---|
|  | Internal battery For battery level information, see 5.11 <i>Using Batteries</i> , page 70. |
|  | Click-in battery (accessory) For battery level information, see 5.11 <i>Using Batteries</i> , page 70. |
|  | Keypad lock activated |
|  | Keypad lock deactivated |
|  | Click-in Humidifier (accessory) The number in the drop indicates humidity setting. If the click-in humidifier is connected but not activated, the symbol is stroked through. |
|  | Heated circuit (accessory) The number in the symbol indicates the set temperature for the heated circuit. If the heated circuit is connected but not activated, the symbol is stroked through. |
|  | Multiple pages Press the MORE button to display the next page. |


5.4 Display Overview

5.4.1 Home Mode

In Home Mode, the Nippy 4+ start screen has the following layout:



Icon explanation:



Page



Sub menu

● Menu



Humidification Settings



User Preferences



Patient Monitor



Compliance Data (Optional)



Alarm/Event Logs



Device Information

5.5 Profiles

Three different profiles can be used for storing complete parameter and alarm settings. This function is suitable as a quick selection for a patient using different settings, for example at night or during daytime.

5.5.1 Selecting a Profile

- 1 Press the **Mode/Profile** button.
- 2 Press the - button to select the profile to use and then press the **Select** button.

The Settings list and the Alarm lists are now displayed. Note that the frames around the settings are red, indicating that the change of profile needs to be confirmed.

- 3 Review the treatment settings and press the **Confirm** button.

The profile is now saved and applied, indicated by blue frames around the settings.


To revert to the original settings, press the **Cancel** button.

5.6 Treatment Settings

5.6.1 Changing a Setting

- 1 Make sure the *Home mode* lock is off.
- 2 Press the **Setting** button for the frame that contains the setting to change.

If the frame contains several settings, press repeatedly until the setting to change is selected.



- 3 Use the + and - buttons to adjust the value and then press the **Set** button.

5.6.2 Settings

5.6.2.1 IPAP

| Item | Description |
|--------------------|---|
| Definition | The IPAP setting is used to define the airway pressure during the inspiratory phase. Minimum/maximum working IPAP is limited/achieved by a software control of blower speed vs. measured pressure. |
| Setting min | 4 cmH ₂ O |
| Setting max | 50 cmH ₂ O |
| Setting resolution | 0.5 below 10 cmH ₂ O, 1.0 above 10 cmH ₂ O |
| Default value | 15 cmH ₂ O |

5.6.2.2 EPAP

| Item | Description |
|-------------|--|
| Definition | The EPAP setting is used to define the airway pressure at the end of the expiratory phase. |
| Setting min | 2 cmH ₂ O Off (for patient circuits with active exhalation valve) |
| Setting max | 20 cmH ₂ O For pressure ventilation modes: The max setting is also limited by IPAP -2 cmH ₂ O and Min Pressure -2 cmH ₂ O. |

| Item | Description |
|--------------------|---|
| Setting resolution | 0.5 below 10 cmH ₂ O,1.0 above 10 cmH ₂ O |
| Default value | 5 cmH ₂ O |

5.6.2.3 Breath Rate

| Item | Description |
|--------------------|---|
| Definition | The Breath Rate setting defines the minimum number of breaths the ventilator will deliver as long as no inspiratory trigger effort from the patient is detected. The cycles will be ventilator-initiated breaths. The combination of the Breath Rate and Inspiratory Time setting is limited by the I:E ratio range 1:9.9 to 2:1. |
| Setting max | 50 bpm |
| Setting resolution | 1 bpm |
| Default value | 12 bpm |

5.6.2.4 Backup Rate

| Item | Description |
|--------------------|---|
| Definition | The Backup Rate setting defines the minimum number of breaths the ventilator will deliver as long as no inspiratory trigger effort from the patient is detected. The cycles will be ventilator-initiated breaths. The combination of the Breath Rate and Inspiratory Time setting is limited by the I:E ratio range 1:9.9 to 2:1. |
| Setting min | 4 bpm 0 bpm (MPV modes) |
| Setting max | 50 bpm |
| Setting resolution | 1 bpm |
| Default value | 12 bpm |

5.6.2.5 SIMV Rate

| Item | Description |
|--------------------|--|
| Definition | The SIMV Rate setting is used in the SIMV ventilation modes, for defining the minimum frequency of mandatory, ventilator-controlled breaths. The mandatory breaths can be either triggered by an inspiratory effort from the patient, or ventilator-initiated. The SIMV Rate setting determines the SIMV cycle time. The combination of the SIMV Rate and Inspiratory Time setting is limited by the I:E ratio 2:1. |
| Setting min | 4 bpm |
| Setting max | 40 bpm |
| Setting resolution | 1 bpm |
| Default value | 12 bpm |

5.6.2.6 Insp. Time (Inspiratory Time)

| Item | Description |
|--------------------|--|
| Definition | The Inspiratory Time setting defines the length of each inspiration from start of inspiration to cycling off to expiration. The combination of the Inspiratory Time and Breath Rate settings is limited by the I:E ratio 2:1. |
| Setting min | 0.3 s |
| Setting max | 5 s |
| Setting resolution | 0.1 s |
| Default value | 1.5 s |

5.6.2.7 Backup Insp. Time (Backup Inspiratory Time)

| Item | Description |
|--------------------|---|
| Definition | The Backup Inspiratory Time setting defines the length of each inspiration delivered during ventilator-triggered backup ventilation, initiated by the set Backup Rate. The combination of the Backup Inspiratory Time and Backup Rate setting is limited by the I:E ratio 2:1. |
| Setting min | 0.3 s |
| Setting max | 5 s |
| Setting resolution | 0.1 s |
| Default value | 1.5 s |

5.6.2.8 Sigh Rate

This setting is available in the menu, after enabling the Sigh function

| Item | Description |
|--------------------|--|
| Definition | The Sigh rate sets the frequency of which breaths with an increased pressure or volume are delivered to the patient. If the High Pressure alarm or the High Tidal Volume alarm is given, the Sigh function will be disabled as long as the alarm condition remains. |
| Setting min | Off, every 10 breaths . |
| Setting max | Every 250 breaths. |
| Setting resolution | 10 breaths |
| Default value | Off |

NOTE



In pressure modes (during the sigh breath), the high pressure alarm will automatically be set 10 cmH₂O above set sigh pressure (max 55 cmH₂O).

5.6.2.9 Sigh %

This setting is available in the menu, after enabling the Sigh function

| Item | Description |
|--------------------|---|
| Definition | Sigh % sets the increased % of the set pressure is delivered to the patient. |
| Setting min | 125% of actual set pressure or volume . |
| Setting max | 200% of actual set pressure or volume . The maximum sigh is also limited by the max allowed set volume. |
| Setting resolution | 25% |
| Default value | 125% |

NOTE



In pressure modes (during the sigh breath), the high pressure alarm will automatically be set 10 cmH₂O above set sigh pressure (max 70 cmH₂O).

5.6.2.10 Rise Time

| Item | Description |
|--------------------|--|
| Definition | The Rise Time setting controls the speed of the pressure/volume increase from start of inspiration to the set pressure or volume. A low setting will give a faster increase and therefore a longer plateau at the set value. A high setting will give a slower increase and therefore a shorter plateau. |
| Setting min | 1 (PSV and PCV modes) 50% of the inspiration time (Min. 0.3 s) (VCV modes) |
| Setting max | 9 (PSV and PCV modes) 90% of the inspiration time (Min. 0.3 s) (VCV modes) Off (VCV modes) |
| Setting resolution | 1 step(PSV and PCV modes) 10% of the inspiration time (VCV modes) |

5.6.2.11 Insp. Trigger (Inspiratory Trigger)

| Item | Description |
|--------------------|--|
| Definition | The inspiratory trigger defines the patient's effort required to initiate a ventilator assisted breath. When the patient starts a breath, an increasing flow is created in the patient circuit. If the patient's effort reaches the set inspiratory trigger level an inspiration is initiated. If the patient cannot trigger a breath, the ventilator will deliver breaths according to the set Backup Rate or Breath Rate. The Assisted breath modes are turned off if the inspiratory trigger is set to Off. |
| Setting min | 1 |
| Setting max | 9 Off (Not PSV nor MPV modes) |
| Setting resolution | 1 (Setting 1 is the most sensitive and 9 is the least sensitive) |
| Default value | 3 |

5.6.2.12 Sup. Pressure (SIMV mode)

| Item | Description |
|--------------------|---|
| Definition | The Support Pressure setting is used in the SIMV ventilation modes, for defining the inspiratory pressure for the support breaths triggered by the patient. |
| Setting min | EPAP +2 cmH ₂ O |
| Setting max | Max IPAP |
| Setting resolution | 0.5 below 10 cmH ₂ O, 1.0 above 10 cmH ₂ O |
| Default value | 15 cmH ₂ O |

5.6.2.13 Exp. Trigger (Expiratory Trigger)

| Item | Description |
|--------------------|--|
| Definition | The Expiratory Trigger setting defines the moment when the ventilator will cycle from the inspiratory to the expiratory phase. |
| Setting min | 1 (10% decrease of peak flow) |
| Setting max | 9 (90% decrease of peak flow) |
| Setting resolution | 1 (Setting 1 is the most sensitive and 9 is the least sensitive) |

5.6.2.14 Max Insp. Time (Maximum Inspiratory Time)

| Item | Description |
|--------------------|--|
| Definition | The Maximum Inspiratory Time setting defines a maximum length for each inspiration. If the Maximum Inspiratory Time is set to Off, the length of the inspiration and/or minimum inspiratory time is dependent on the set Expiratory Trigger. |
| Modes | PSV, PSV(TgV) |
| Setting min | 0.3 s |
| Setting max | 5 s Off |
| Setting resolution | 0.1 s |
| Default value | Off |

5.6.2.15 Min Insp. Time (Minimum Inspiratory Time)

| Item | Description |
|--------------------|--|
| Definition | The Minimum Inspiratory Time setting defines a minimum length for each inspiration. If the Minimum Inspiratory Time is set to Off, the length of the inspiration is dependent on the set Expiratory Trigger. |
| Setting min | <ul style="list-style-type: none"> • 0.3 s • Off |
| Setting max | 3 s |
| Setting resolution | 0.1 s |
| Default value | Off |

5.6.2.16 Target Volume

NOTE



If Target Volume is used with a patient circuit with an active exhalation valve, leakage may be misinterpreted by the ventilator as an increase of tidal volume. This will lead to a decrease of the Inspiratory Pressure (the Inspiratory Pressure will not be lower than the set Min Pressure). This may result in hypoventilation as the true delivered tidal volume will decrease both as a result of the leakage and the decrease in Inspiratory Pressure. This does not occur if a patient circuit with leakage port is used.

| Item | Description |
|--------------------|--|
| Definition | The Target Volume setting defines the tidal volume that the ventilator will aim for while ventilating the patient in a pressure mode. To aim for the preset volume, the ventilator will adapt the Inspiratory Pressure between two adjustable pressure limits: Min Pressure and Max Pressure. When Target Volume is active, the mode field on the ventilator display will indicate "(TgV)". |
| Setting min | Off 100 ml |
| Setting max | 2000 ml |
| Setting resolution | 10 below 500 ml,50 above 500 ml |
| Default value | Off |

5.6.2.17 Tidal Volume

| ITEM | DESCRIPTION |
|--------------------|---|
| Definition | The Tidal Volume setting defines the volume that will be delivered by the Nippy 4+ for each breath when using volume control modes. |
| Setting min | 100 ml |
| Setting max | 2000 ml |
| Setting resolution | 10 ml below 500 ml 50 ml above 500 ml |

5.6.2.18 Flow Pattern

| ITEM | DESCRIPTION |
|---------------|---|
| Definition | Flow pattern sets the characteristics of the air flow in VCV modes. |
| Settings | Square (constant flow during the inspiratory phase) Decelerating (flow decreases linearly, may prevent air hunger) |
| Default value | Square |

5.6.2.19 Max Pressure

| Item | Description |
|--------------------|--|
| Definition | The Max Pressure setting is only used when Target Volume is activated. Max Pressure defines the upper pressure limit up to where the ventilator can increase the pressure to reach the set Target Volume. If Target Volume is not reached at Max Pressure, the ventilator will continue to ventilate at this Max Pressure setting. |
| Setting min | Current Min pressure setting |
| Setting max | 50 cmH ₂ O |
| Setting resolution | 0.5 below 10 cmH ₂ O, 1.0 above 10 cmH ₂ O |
| Tolerance | ±0.5 cmH ₂ O or ±5 %, whichever is the greatest. |
| Default value | 15 cmH ₂ O |

5.6.2.20 Min Pressure

| Item | Description |
|--------------------|---|
| Definition | The Min Pressure setting is only used when Target Volume is activated. Min Pressure defines the lower pressure limit down to where the ventilator can decrease the pressure to maintain the set Target Volume. If the actual volume is above Target Volume at Min Pressure, the ventilator will continue to ventilate at this Min Pressure setting. |
| Setting min | 4 cmH ₂ O |
| Setting max | Current Max pressure setting |
| Setting resolution | 0.5 below 10 cmH ₂ O, 1.0 above 10 cmH ₂ O |
| Tolerance | ±0.5 cmH ₂ O or ±5 %, whichever is the greatest. |
| Default value | 15 cmH ₂ O |

5.6.2.21 CPAP

| Item | Description |
|--------------------|---|
| Definition | The CPAP setting defines the pressure that will be applied to the airways in CPAP mode. |
| Setting min | 4 cmH ₂ O |
| Setting max | 20 cmH ₂ O |
| Setting resolution | 0.5 below 10 cmH ₂ O, 1.0 above 10 cmH ₂ O |
| Max bias error | ±0.8 cmH ₂ O or ±4 %,whichever is the greatest. |
| Default value | 10 cmH ₂ O |

5.6.2.22 Humidifier

| Item | Description |
|---------------|--|
| Definition | Humidifier allows the user to start or stop the heated humidification. The click-in water chamber needs to be connected before the setting can be turned On. |
| Setting min | Off |
| Setting max | On |
| Default value | Off |

5.6.2.23 Humidifier Setting

| Item | Description |
|--------------------|---|
| Definition | The Humidifier Setting defines the level of humidity of the air delivered to the patient. |
| Setting min | 1 |
| Setting max | 5 |
| Setting resolution | 1 |
| Default value | 3 |

5.6.2.24 Heated Circuit Temp

| Item | Description |
|--------------------|--|
| Definition | Heated Circuit Temp setting will define the temperature of the heated circuit. |
| Setting min | 16 °C (61 °F) |
| Setting max | 30 °C (86 °F) |
| Setting resolution | 0.5 ° |
| Default value | 27 °C (81 °F) |

5.6.2.25 Circuit Heating


| Item | Description |
|---------------|---|
| Definition | Circuit Heating allows the user to start or stop the heating of the circuit. The Heated Circuit needs to be connected before the setting can be turned On. |
| Setting min | Off |
| Setting max | On |
| Default value | Off |

5.7 The Alarms List

The Alarms list displays the alarm settings. Which alarms that are available depends on the current treatment mode and settings.

5.8 The Monitored Values Pane

This section describes the monitored treatment values that are displayed on the start screen



| Value | Description |
|------------------|--|
| RR | Respiratory rate |
| Vte / Vti | Tidal Volume (the volume delivered with each breath) <ul style="list-style-type: none"> Expiratory volume (Vte) is displayed for leakage circuits. Vte is a calculated value. Inspiratory volume (Vti) is displayed for active exhalation valve circuits and for mouthpiece circuits. |
| MVe / MVi | Minute Volume, calculated as Tidal Volume multiplied with the Total Breath Rate. <ul style="list-style-type: none"> Expiratory volume (MVe) is displayed for leakage circuits. Inspiratory volume (MVi) is displayed for active exhalation valve circuits and for mouthpiece circuits. |
| LEAK | The total leakage (intentional and unintentional) as calculated at expiratory pressure level. |
| HR | Heart rate An SpO ₂ sensor needs to be in place to measure and display this value. |
| SpO ₂ | Displays the patient's oxygen saturation. An SpO ₂ sensor needs to be connected to measure and display this value. |

| Value | Description |
|--------------------|---|
| FiO ₂ | Displays the fraction of inspired oxygen as measured at the air outlet of the Nippy 4. An FiO ₂ sensor needs to be connected to measure and display this value. |
| EtCO ₂ | Displays the end-tidal carbon dioxide, measured on the last portion of the exhaled volume that is passing through the EtCO ₂ sensor. An EtCO ₂ sensor needs to be connected to measure and display this value. |
| PtcCO ₂ | Displays transcutaneous CO ₂ pressure from an external PtcCO ₂ monitor. An PtcCO ₂ monitor needs to be connected to measure and display this value. |
| Pressure bar graph | Displays the pressure during treatment. <ul style="list-style-type: none"> To the left of the bar a mark discloses the highest pressure during last breath. To the right of the bar marks discloses the set values for IPAP and EPAP. Red lines in the bar indicates alarm levels. |
| Flow bar graph | Displays the flow during treatment. <ul style="list-style-type: none"> To the left of the bar a mark discloses the highest flow during last breath. |
| Volume bar graph | Displays the air volume delivered during treatment. <ul style="list-style-type: none"> To the left of the bar a mark discloses the total volume delivered during last breath. To the right of the bar marks discloses the set target volume value (if used). Red lines in the bar indicates alarm levels |


In clinical mode, curves, trends and additional values can be viewed from the Patient Monitor page, see .

5.9 The Menu


This section contains information about the items on the menu.

5.9.1 Opening the Menu

- 1 Press the **Menu** button.




- 2 In the menu, use the **+** and **-** buttons to select the item to open and then press the **Set** button to open it.



5.9.2 Humidification Settings

This menu item lets you make changes to the humidifier and the heated circuit settings.

- 1 Press the **Menu** button and then select **Humidification Settings**.
- 2 Use the **+** and **-** buttons to select the setting to change and click the **Edit** button.
- 3 Use the **+** and **-** buttons to change the value and then click the **Edit** button to leave the editing mode for the specific setting.
- 4 Click **Close** to save the settings when done.



5.9.3 User Preferences

This menu item lets you view or make changes to the user preferences:


- Time and Date
- Display Light
 - On (will keep the display lit up regardless of use)
 - Auto (will adjust the light intensity depending on the ambient light)

— Delayed (the display is dimmed after 30 seconds or more depending on the mode and battery setup. If any button is pressed or any alarm occurs, the display light will return to normal again).

- Light Intensity (setting range: 1-9, where 1 is the lowest and 9 is the highest light intensity setting).
- Alarm sound Level (setting range: 1-5, where 1 is the lowest and 5 is the highest. Make sure to use a setting where the alarm sound is clearly audible.)
- Pressure Unit
- CO₂ Unit
- Compliance in Home Mode (Show/Hide)

Change User Preferences

- 1 Press the **Menu** button and then select **User Preferences**.
- 2 Use the **+** and **-** buttons to select the setting to change and click the **Edit** button.
- 3 Use the **+** and **-** buttons to change the value and then click the **Edit** button to leave the editing mode for the specific setting.
- 4 Click **Close** to save the settings when done.



5.9.4 Compliance Data (Optional)

Compliance data is statistics of usage, such as usage hours and days used.

From this menu item you view the compliance data, if the Nippy 4+ has been configured to show it.

5.9.5 Pre-Use Test

This menu item is only available in clinical mode and when the Nippy 4+ is in standby mode.


This menu item contains:

- Start Pre-use Test On/Off (Selects whether a pre-use test shall be prompted every time the Nippy 4+ is powered on).

The default setting for pre-use test is On.

- START PRE-USE TEST (starts a pre-use test immediately).

- 1 Press the **Menu** button and then select **Pre-use Test**.
- 2 Use the **+** and **-** buttons to select the item and then click the **Select** button.
- 3 If configuring the pre-use test prompt, use the **+** and **-** buttons to change the value and then click the **Select** button to leave the editing mode.
- 4 Click **Close** when done.



5.9.6 Alarm/Event Logs

This menu item displays the alarm and event logs.



- Red: High priority alarms
- Yellow: Medium priority alarms
- Blue: Messages
- White: Currently selected row

Click **Close** when done.

5.10 Transferring Data between the Ventilator and a PC

WARNING!



Read the chapter 2.2 *Electrical Safety*, page 21 carefully to make sure all conditions are considered and met.

CAUTION!



Do not eject the memory card or disconnect the USB cable while the Nippy 4+ is transferring data. Doing so may result in loss of data and/or damaged equipment.

NOTE



In order to view and present patient data, Breas software must be used.



Instructions on how to manage data in the Breas software can be found in the software help.


5.10.1 Transferring Data with a Memory Card

NOTE




The *Home mode* lock must be off for copying and transferring data to the memory card.

- 1 Insert the memory card in the memory card slot on the side of the Nippy 4+. Make sure the memory card is properly inserted.



- 2 Press the Menu button and navigate to the Device Memory page. Menu > Advanced Settings > Device Memory.

- 3 Select **Save Memory Data on Card** and press the **Select** button. Confirm to save the data and wait while the data is being saved.




- 4 Remove the memory card from the ventilator and insert it to the computer. You need Breas software to read the data on the card.

5.10.2 Transferring Data with a Data Cable

With a USB cable, real-time data can also be received and sent between the ventilator and a PC.

- 1 Connect the USB cable to the ventilator. Make sure it is fitted correctly.



- 2 Connect the other end of the cable to a PC running Breas PC software.

WARNING!



The PC must be placed outside the patient area (i.e. more than 2 meters from the patient).

5.11 Using Batteries

Since all batteries, in general, degenerate over time, the recommendations below will ensure that the battery capacity of the Nippy 4+ is maximized during its lifetime.

The internal and click-in batteries in the ventilator are of the Lithium-ion type, which is a high performance battery. It has a long expected lifetime, low weight in relation to its capacity and low self discharge.



See the Nippy 4+ Service Manual on how to perform service on the batteries.

5.11.1 Power Source Priority

1. AC power (Mains)
2. External DC
3. Click-in battery
4. Internal battery

If the AC power source fails, the ventilator will switch to either the external DC (if installed), or the click-in battery (if attached) or the internal battery and show a message in the display window.

NOTE



How to test the Internal Battery:

The switch-over to internal battery can be tested by disconnecting the AC power cord, and confirming the behaviour described below is observed.

- Internal battery power source LED will be illuminated
- Medium priority “Lost Mains Power” alarm will be triggered
- Information message “Switched to Internal Battery” will be posted



How to test the Click-In Battery:

The switch-over to click-in battery can be tested by disconnecting the AC power cord while having a click-in battery connected, and confirming the behaviour described below is observed.

- Click-in battery power source LED will be illuminated
- Medium priority “Lost Mains Power” alarm will be triggered
- Information message “Switched to Click-In Battery” will be posted



How to test the External DC (“alternative Supply Mains”):

The switch-over to external DC can be tested by disconnecting the AC power cord while having an external DC source connected, and confirming the behaviour described below is observed.

- External DC power source LED will be illuminated
- Medium priority “Lost Mains Power” alarm will be triggered
- Information message “Switched to External DC” will be posted

5.11.2 Charging the Batteries



CAUTION!

Do not charge the ventilator while placed in the carry bag or other types of closed or non-ventilated spaces.



Charging of batteries is only started when the state of charge is below 95%.

The internal and click-in batteries are automatically charged when connecting the Nippy 4+ to the mains supply. To ensure that the batteries are fully charged, a maintaining charging cycle will be performed.

The batteries are not charged when connecting the Nippy 4+ to an external DC supply. While charging, the battery level will be animated. The batteries are only charged if the internal temperatures are between 0 to 45°C (32 to 113°F). High power consuming settings in combination with high ambient temperatures may make the battery temperature rise above 45°C (113°F).

Behaviour of the Ventilator while Internal Battery or Click-in Battery is Charging

The battery icon will be animated (filling from bottom to top).




Charging Times



| Battery | Charger | Time |
|------------------|--------------------------|------|
| Internal battery | Nippy 4+ | 2 h |
| Click-in battery | Nippy 4+ | 4 h |
| Click-in battery | Click-in battery charger | 3 h |

Times are based on charging empty batteries.

5.11.3 Battery Icons

When running on battery, the battery status is indicated by the following symbols:

| Symbol | Description |
|---|---|
|  | Internal battery Green symbol indicates over 50 % state of charge. |
|  | Click-in battery Green symbol indicates over 50 % state of charge. |
|  | Medium State of charge (20 % – 50 %) |

| Symbol | Description |
|---|----------------------------------|
|  | Low state of charge (below 20 %) |
|  | Malfunctioning battery |

5.11.4 Internal Battery

The internal battery is intended as a backup power source if the primary power source fails. It can also be used as a temporary power source. For example during transportation between one stationary power source to another.


5.11.5 Click-in Battery

The click-in battery is intended as a power source during transportation, or if the primary mains power source fails.


The click-in battery can be replaced during treatment, provided that the internal battery is charged.

Connect the Click-in Battery


- 1 Release the side cover by pressing the button under the handle.




- 2 Open and remove the side panel.



3 Insert the click-in battery.



4 Close the side panel. Make sure there is a clicking sound to secure the side panel.



When removing the battery, press down the latch at the bottom of the battery compartment and tilt the ventilator sideways. Make sure to close the side panel after removing the click-in battery.

5.11.6 Battery Operating Time (Internal and Click-in)

The operation time is dependent on the battery condition, its capacity, the ambient air temperature and the Nippy 4+ pressure setting. These data are based on new and fully charged batteries.

| Parameter | Example |
|---------------------------------|-----------------------|
| Environmental Conditions | |
| Ambient temperature | 20°C (68°F) |
| Ventilator Settings | |
| Mode | PCV |
| IPAP | 20 cmH ₂ O |
| EPAP | 4 cmH ₂ O |
| Breath Rate* | 20 bpm |
| Insp. Time* | 1.0 s |
| I:E | 1:2 |

| Parameter | Example |
|------------------|---------------------------|
| Insp. Trigger | Off |
| Rise Time | 1 |
| Target Volume | Off |
| Display Light* | Off |
| Light Intensity* | - |
| | |
| Monitored Value | |
| Tidal Volume | 800 ml |
| Resistance | 5 hPa (l/s) ⁻¹ |
| Compliance | 50 ml (hPa) ⁻¹ |

*: These ventilator settings affect the operating time significantly

| Battery | Operating Time |
|------------------|----------------|
| Internal Battery | 2.5 h |
| Click-in Battery | 6.5 h |

5.11.7 Storing the Internal Battery and the Click-in Battery

Storage longer than 1 month should be initiated with half-charged batteries in order to maintain maximum capacity. Optimal storage temperature is 5 to 30°C (41 to 86°F).

5.11.8 External DC



WARNING!

Do not connect the ventilator to a wheelchair unless the operating manual for the wheelchair permits this as this can affect the ventilator performance and consequently result in patient death.



CAUTION!


Only use a battery charger compliant to IEC 606011 if you are charging a battery that at the same time is connected to the ventilator.

The ventilator can be operated from:

- a 12 V external DC source using the 12/24 V converter.
- a 24 V external DC source using the external battery cable.

With an external DC source connected, the Nippy 4+ will automatically switch over to the external DC source if the mains power cord is removed or if the mains power supply fails. The external DC voltage level is shown under “Device Information” in the menu.

- 1 Connect the external DC cable to the Nippy 4+. Make sure that it is fitted correctly.




- 2 Connect the other end of the cable to the DC source.


5.12 Using Accessories

5.12.1 Connecting and Disconnecting the Cables


- 1 Insert the cable in the appropriate port.



- 2 Make sure to insert the connector with the marking pointing upwards.



- 3 Pull the connector sleeve, not the cable itself or cable restrainer to release the connector.




5.12.2 Using the ventilator with a Nurse Call System

The ventilator can be connected to a nurse call system using the nurse call cable. When connected, the ventilator alarms will also be forwarded to the nurse call system.

5.12.2.1 Connect the ventilator to a Nurse Call System

- 1 Connect the nurse call cable at the back of the ventilator.



- 2 Test the connection by triggering an alarm on the ventilator and verify that the nurse call system activates.

5.12.3 Using the ventilator with the FiO₂ Sensor

The FiO₂ sensor can be used to monitor and store FiO₂ measurements. The FiO₂ sensor measures the fraction of inspired oxygen in the air channel of the ventilator. The FiO₂ measurements will be stored in the data memory which can be downloaded to a PC and viewed in Breas software.

| Usage | Time |
|-------------------------|--|
| Operating temperature | 10 to 40°C (50 to 104°F) |
| Operating pressure | 700 to 1250 mbar |
| Expected operating life | <3 years (in ambient air) or 500,000 Vol.% h, whichever comes first. |
| Shelf life | < 6 months (recommended) |




CAUTION!


Note that the operating conditions for the FiO₂ sensor is different from the ventilator system conditions. If the sensor is used outside its operating conditions the FiO₂ measurements might deviate.

5.12.3.1 Installing the FiO₂ Sensor

- 1 Place the ventilator so the bottom is accessible.
- 2 Remove the hatch for FiO₂ sensor.
Use a torx TX10 screwdriver.



- 3 Insert the FiO₂ sensor with the electrical contact side in.



- 4 Reinstall the hatch.
- 5 Calibrate the FiO₂ sensor in the advanced settings of the main menu.



When installed, the ventilator automatically detects the sensor, also after powering off/on and after power failure.

5.12.3.2 Calibrating the FiO₂ sensor

The FiO₂ sensor should be calibrated when first used and then at least once a month.



FiO₂ calibration can be performed from the “FiO₂/CO₂ Calibration” page under the Advanced settings section of the main menu.

5.12.4 Using the ventilator with the Remote Alarm



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for Remote Alarm.

The Remote Alarm enables care providers and clinical personnel to monitor the ventilator alarms remotely. The Remote Alarm forwards alarms from the ventilator. When an alarm sounds, the care provider or clinical personnel must attend to the patient quickly.

When installing a remote alarm system, check that it operates as intended before starting the treatment.

5.12.5 Using the ventilator with the EtCO₂ Sensor

The EtCO₂ sensor can be connected to the patient breathing circuit and to a Nippy 4+ in order to monitor and store CO₂ measurements. The CO₂ measurements will be stored in the ventilator data memory which can be downloaded to a PC and viewed in the ventilator PC software.



More information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for the EtCO₂ sensor.

5.12.5.1 Safety Information



WARNING!

Read the instructions thoroughly so that you completely understand how the EtCO₂ sensor is operated before taking it into use, to ensure correct usage and maximum performance.

Breas Medical reserves the right to make changes to this product without any prior notification.



Do not use a damaged CO₂ sensor or adapter.
The CO₂ sensor is intended to be used by authorized and trained medical personnel only.



The CO₂ sensor is intended only as an adjunct in patient assessment. It must be used in conjunction with other assessments of clinical signs and symptoms.




Masks' dead-space, patient's volumes and unintentional leakage may influence the CO₂ measurements.



Used airway adapters shall be disposed of in accordance with local regulations for medical waste.



Measurements can be affected by mobile and RF communications equipment. It should be assured that the CO₂ sensor is used in the electromagnetic environment specified in 8.3 *Emission and Immunity Declaration*, page 154.



Do not place the airway adapter between the endotracheal tube and an elbow as this may allow patient secretions to block the adapter windows and result in incorrect operation.



Incorrect CO₂ zeroing will result in false gas readings.
Replace the airway adapter if rain-out/condensation occurs inside the airway adapter.



Only use airway adapters distributed by Breas Medical.
Do not apply tension to the CO₂ sensor cable.



To keep secretions and moisture from pooling on the windows, always position the CO₂ sensor in a vertical position with the green LED pointing upwards.

WARNING!



Disposable airway adapters shall not be reused. Reuse of the single use adapter can cause cross infection.

CAUTION!



If an intentional leakage port is used, make sure that the CO₂ sensor is placed between the patient interface and the leakage port.
If a patient interface with integrated leakage is used, the monitored CO₂ values may be influenced.



The CO₂ sensor should be placed as close to the patient interface as possible. However, a HME (if used) shall be placed between the patient interface and the CO₂ sensor. This will protect the airway adapter from secretions and effects of water vapour and eliminates the need of changing the airway adapter.

NOTE




The CO₂ monitoring automatically compensates for changes in ambient barometric pressure. The CO₂ monitor shall fulfil the ISO 80601-2-55 standard (Medical electrical equipment – Particular requirements for the basic safety and essential performance of respiratory gas monitors).


5.12.5.2 How to Connect the EtCO₂ Sensor

- 1 Connect the CO₂ sensor cable to the CO₂ connection port on the ventilator (according to the instruction 5.12.1 *Connecting and Disconnecting the Cables*, page 76).
A green LED indicates that the CO₂ sensor is ready to use.


- 2 Snap the CO₂ sensor probe on top of the airway adapter. It will click into place when properly sealed.



- 3 Perform a CO₂ zeroing procedure.
- 4 Connect the airway adapter to the patient circuit.



- 5 Make sure to position the CO₂ sensor with the LED pointing upwards.



When installed, the ventilator automatically detects the sensor, also after powering off/on or after power failure.

WARNING!



The CO₂ sensor is not intended to be in contact with the patient body.

CO₂ Zeroing

CO₂ zeroing is recommended when changing the airway adapter. Besides from that, zeroing only needs to be performed when an offset in monitored CO₂ values is observed, or when a CO₂ sensor accuracy unspecified message is displayed.

| LED Status | Description |
|----------------------|---------------------|
| Steady green light | System OK |
| Flashing green light | Zeroing in progress |
| Steady red light | Sensor error |
| Flashing red light | Check adapter |

Maintenance

No periodical maintenance is required for the CO₂ sensor.

To verify the CO₂ sensor readings, a gas span check shall be performed every year, preferably when the ventilator is sent for service.



See the ventilator service manual for how to perform the gas span check.

WARNING!



Do not under any circumstances attempt to service or repair the CO₂ sensor yourself. If you do so, the manufacturer will no longer be responsible for the performance and safety of the CO₂ sensor.

Cleaning

WARNING!



- Always be careful when cleaning to ensure that you do not damage any equipment.
- Fluid must not be allowed to enter into the CO₂ sensor.
- Always clean the T-piece with plug when to be used by a new patient. All parts that come into contact with the respiration gas must be cleaned.



- Remove the airway adapter before cleaning.
- Do not sterilise the CO₂ sensor.
- Do not autoclave the CO₂ sensor.

Clean the outside of the CO₂ sensor using a lint-free cloth moistened, but not wet, with ethanol or isopropyl alcohol (< 70%).

Disposal

The CO₂ sensor must be disposed of in accordance with the local environmental regulations regarding the disposal of used equipment and waste.

5.12.6 Using the ventilator with the PtcCO₂ Cable

An external monitor for Transcutaneous CO₂ Pressure (PtcCO₂) may be connected to the ventilator by an accessory PtcCO₂ cable. For information about available PtcCO₂ cables, see 9 *Accessories and Parts*, page 164.

NOTE



Both the PtcCO₂ cable and the EtCO₂ sensor connects at the ventilator's yellow CO₂ port. Only one CO₂ measuring device can be connected at a time.

When connected, the ventilator will:

- Display the monitored values and include them in trend views.
- Store monitored values in the internal memory. The PtcCO₂ values will also be included in the data that can be downloaded and analysed with Breas PC software.
- Repeat CO₂ alarms from the external PtcCO₂ monitor.
- Automatically detect the sensor, also after powering off/on or after power failure

5.12.7 Using the Ventilator with the SpO₂ module



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for SpO₂ module.

The SpO₂ module enables connection to an SpO₂ sensor for measuring of functional oxygen saturation of arterial haemoglobin (SpO₂) and pulse rate. The SpO₂ module can be connected to the Nippy 4+ in order to monitor and store SpO₂ measurements.

The SpO₂ measurements will be stored in the data memory which can be downloaded to a PC and viewed in the Breas PC software.

When installed, the Nippy 4+ automatically detects the sensor, also after powering off/on and after power failure.

5.12.8 Using the Ventilator with the Effort Belts

NOTE



The effort belt communication box and the Remote Start/Stop use the same port on the ventilator. Only one of the accessories can be connected at a time.

The Effort Belt Communication Box


Up to two effort belts may be connected to the ventilator by the accessory Effort belt communication box.

When the effort belt communication box is connected, the ventilator will:

- Perform a start up test when an effort belt is connected.
- Check the effort belt unit for internal failures.
- Include real time wave forms from the effort belts on the Effort monitoring page.
- Store the effort belt measurements in the internal memory. The effort belt measurements will also be included in the data that can be downloaded and analysed with Breas PC software.
- Automatically detect the belt, also after powering off/on or after power failure

Connecting the Effort Belt

1. Connect the black pins of the wire set to the effort belt.
2. Connect the key-shaped connector of the wire set to the communication box.
3. Connect the Effort belt communication box to the ventilator.



Effort Belt Connection Status

The connection status of the effort belts are indicated by LEDs on the effort belt communication box and by the effort belt symbol on the display

| Status | Indication |
|--|--|
| Connection OK | <ul style="list-style-type: none">• Communication box Green LED for each effort belt.• Nippy 4+ Display White belts in the effort belt symbol. |
| Communication box connected but not belt | <ul style="list-style-type: none">• Communication box Red LED for unconnected belt. Green LED for connected belt.• Nippy 4+ Display White belt in the effort belt symbol for connected belt. Red belt in the effort belt symbol for unconnected belt. |
| Communication box not connected | <ul style="list-style-type: none">• Communication box No light from LEDs.• Nippy 4+ Display No effort belt symbol. |

5.12.9 Using the ventilator with the Remote Start/Stop

NOTE




The effort belt communication box and the Remote Start/Stop use the same port on the ventilator. Only one of the accessories can be connected at a time.

5.12.9.1 Connecting the Remote Start/Stop

- 1 Connect the Remote Start/Stop cable to the ventilator.



Information about safety, warnings, product description, installation, usage, cleaning, maintenance and technical specifications can be found in the user instruction for Remote Start/Stop.




5.12.10 Using the ventilator with the Protective Cover

The protective cover is intended for additional protection of the ventilator during transportation, and in hospital, institutional and home care environments. It can be used while the ventilator is operating, for example mounted on a wheelchair, in a personal vehicle, or carried by hand.

The protective cover protects the ventilator from environmental impact such as shock, water spill, sunlight, dust and dirt, under normal handling.

The protective cover has the following functions:



1. Transparent window, for accessing front panel and buttons
2. Port for patient circuit
3. Cooling air inlet
4. Port for cables and O₂ inlet
5. Mounting straps
6. Patient air inlet

CAUTION!



Do not cover the air inlets or outlets.


5.12.11 Using the Nippy 4+ with the Trolley

Intended Use

The intended use of the trolley system is to allow patient mobility while receiving ventilator treatment. The trolley shall only be used in indoor, hospital environment. The trolley system consists of a trolley base and a mounting bracket.

This section describes how to use the Nippy 4+ and a trolley with mounting bracket.

Mount and dismount the Nippy 4+ as shown in the picture:



The bottom plate is mounted to the trolley using two screws.

Be careful when handling the trolley with the ventilator mounted, in order to avoid any risk of the trolley falling. The trolley can be tipped 10° and return to vertical position, when loaded in accordance with the weight specifications below.

WARNING!



The maximum total weight of the trolley and added accessories is 37 kg (82 lbs). (Trolley base weight = 12 kg (26 lbs), maximum externally added load = 25 kg (55 lbs).)

- The maximum load of the trolley basket is 0.9 kg (2 lbs).
- The maximum load of the IV-pole is 3 kg (6.5 lbs).
- The maximum load of the trolley rail is 9 kg (20 lbs).
- The maximum load of the E-cylinder holder is 7.9 kg (17.5 lbs).

No maintenance is required.

5.12.12 Using the Click-in Humidifier

WARNING!



Read the chapter “Humidification” on page 24 before using the Nippy 4+ with the humidifier.

CAUTION!



The click-in humidifier and the circuit heating operates on the AC power source only. If the AC power source fails and the internal or the external battery activates, the click-in humidifier and the circuit heating will be turned off automatically.

NOTE



Maximum Humidifier Start-Up Period in Normal Use: 30 minutes

The click-in humidifier is intended to humidify the patient air. It is intended for non-invasive use only. The click-in water chamber is for single patient use only. Reusing a water chamber for a new patient might cause a risk of cross-contamination. The Nippy 4+ shall not be moved with a filled water chamber installed.

5.12.12.1 Prerequisites

- The water chamber must be installed in order to access the humidifier setting on the ventilator menu, both in clinical and home mode. If the water chamber is disconnected and reconnected after usage, the ventilator will remember the humidity setting used.
- The click-in humidifier only operates during treatment. When the ventilator is in standby mode, the humidification is paused.

5.12.13 Installing the Water Chamber

CAUTION!




Do not switch on the humidifier without a filled water chamber in order to avoid burn or damage to the humidifier's electronics.

NOTE



If the ventilator is equipped with click-in battery, remove it before installing the water chamber.

Follow the instructions in the illustration below to install the water chamber to the ventilator




CAUTION!



Always make sure the water chamber is in correct position before use.
Store the airway bypass unit in a clean and dust free environment.

5.12.14 Detaching the Water Chamber



CAUTION!



Always insert the air bypass unit after disconnecting the water chamber.

WARNING!



Always stop treatment before detaching or attaching the water chamber. Make sure the Nippy 4+ with the attached water chamber is placed lower than the patient and on a flat surface. This is to prevent personal injury from accidental spillage or from excess water or condensation flowing down the patient tube and into the patient's interface.



Never add or pour out water from the water chamber when it is attached to the ventilator. If there is water outside of the water chamber after filling, dry it using a lint-free cloth before reconnecting it to the ventilator.

WARNING!



To avoid burn injury, be careful not to touch the heater plate or the heated water in the water chamber when the humidifier is switched on or has not yet cooled down. Wait 10 minutes for the heater plate and water to cool.

5.12.14.1 Adding Water to the Water Chamber

CAUTION!



Use only distilled or sterilised water or boiled, chilled tap water in the humidifier water chamber. This is to reduce mineral deposits and maximize the life of the water chamber.



Do not fill the water chamber with hot water.



Do not overfill the water chamber. Fill only the water chamber to the maximum level indicated on the water chamber.



Always ensure the lid with seal is properly mounted after filling and reassembling the water chamber. Also check that the water chamber is correctly docked in place and locked to the ventilator.



Avoid to remove the seal from the lid at normal, daily usage.



Make sure all parts are dry before the ventilator is connected to the mains and put into operation.

NOTE



Duration of Operation between Humidifier refills

Default setting (3): 16 hours and 40 minutes


Max setting (5): 8 hours and 40 minutes

- 1 Detach the water chamber, see 5.12.14 *Detaching the Water Chamber*, page 93.

- 2 Fill water to the chamber, by filling through one of the airway connections.


Make sure not to fill above the Max indication. A water chamber filled to the maximum level contains approximately 350 ml

You can also remove the lid and fill water through the top of the chamber.



5.12.14.2 Opening the Water Chamber

The water chamber lid shall be opened when manually emptying or cleaning the water chamber



CAUTION!



Always make sure the lid of the water chamber is totally sealed.

5.12.14.3 Activating the Humidification

- The water chamber shall be filled with water and attached.
 - The ventilator shall be connected to the mains power supply
1. In the **Main menu**, select **Humidification Settings**.
 2. Select **Humidifier Setting** and set the level of humidification. 1 is the lowest level 5 is the highest level.
 3. Select **Humidifier** and set it to **On**.
 4. The humidifier is now activated and will start to operate when the treatment starts.

5.12.14.4 Cleaning the Water Chamber

- 1 Open the water chamber as described in 5.12.14.2 *Opening the Water Chamber*, page 95.
- 2 Clean the parts of the water chamber either by hand using a mild detergent or in a dishwasher.
- 3 If there are mineral deposits inside the water chamber, dissolve them using warm water and citric acid for 30 minutes.

For disinfecting the water chamber, use any of the agents listed below. The water chamber will withstand at least 20 disinfections without degradation.

| Disinfection Agent | Duration |
|--------------------------|------------|
| Gigasept® FF 5% solution | 15 minutes |
| Steranios 2% solution | 10 minutes |

5.12.15 Using the Patient Circuit with Heated Wire

The ventilator may be used with the accessory *Patient Circuit, Heated Wire with Cable Connector*.

When the heated circuit is used, the time for the patient air temperature to reach the set temperature from a starting temperature of $(23 \pm 2)^\circ\text{C}$ may be up to 3 minutes.

Prerequisites

The wire heating only operates during treatment. When the ventilator is in standby mode, the wire heating is paused.



Read the User Instruction for the Patient Circuit, Heated Wire with Cable Connector before using the patient circuit.

5.12.15.1 Connecting the Patient Circuit

Connect the circuit as described in 4.4 *Connecting the Patient Circuit*, page 44.

When the circuit is connected, continue with activating the circuit heating.

5.12.15.2 Activating the Circuit Heating

The ventilator shall be connected to the mains power supply

- 1** In the **Main menu**, select **Humidification Settings**.
- 2** Select **Heated Circuit Temp** and set the temperature according to the respiratory therapist's prescription.
- 3** Select **Circuit Heating** and set it to **On**.

The circuit heating is now activated and will start to operate when the treatment starts.

6 Alarms

WARNING!



The adjustable alarm settings should be re-evaluated whenever a change in settings is made on the ventilator.

CAUTION!



Never leave a patient unattended during an alarm condition.



Setting alarm limits to extreme values could put the patient at risk. Permitted distributed alarm systems are Nippy 4+ remote alarm with cable and Nippy 4+ nurse call cables provided by Breas Medical only.

NOTE



The alarm settings are maintained during an extended power failure.

This chapter describes the alarm functions used for the ventilator.

6.1 Alarm Function

The alarm function of the ventilator consists of the alarm LEDs on the front panel, an audible alarm, and messages on the display (see the front panel section for an overview of the position of the LEDs).

6.1.1 Alarm Indication


As soon as an alarm condition is detected, the ventilator main unit and the remote alarm unit (if connected) will alarm without delay.

When an alarm condition arises, the alarm is indicated in three ways:

Colour LED on the panel

Indicates the priority of the active alarm condition.

- High priority: red colour, flashing twice per second.
- Medium priority: yellow colour, flashing every 2 seconds.




Alarm text in display

Displays the name of the active alarm condition and a guiding text.

The priority of the alarm is indicated by the background colour

- Red = High priority
- Yellow = Medium priority



Audible signals

- **High priority:** 3 signals followed by 2 more. The signal sequence is repeated with a 0.5 second pause and thereafter a 3 second pause.



- **Function failure:** Same signal as the high priority alarm or a constant signal, depending on the kind of function failure.
- **Medium priority:** 3 signals, with a lower frequency than the high priority alarm. The signal sequence repeats after a 6 second pause.



- **Information:** 1 signal with a low frequency. The signal is repeated after a 5 second pause and stopped after 5 sequences.



Alarm signal sound pressure: Adjustable from 45 to 85 dB(A) measured at 1 m. Accuracy: ± 5 dB(A).



The power failure alarm sounds in the case of power failure.


If the external DC falls below the warning limit and it is the last power source, the Low External DC warning is displayed.

If a battery that is the last power source falls below the warning limit, the Low Last Power Source alarm is set.

6.1.2 Audible Signal Pause and Reactivation

The audible signal can be paused for 60 seconds by pressing the Audio Pause button. The audible signal can be reactivated by pressing the Audio Pause button again.

If a new alarm condition occurs during the audio pause period, the audible signal will be reactivated.



6.1.3 Alarm Reset

An alarm will automatically be reset once the cause of the alarm has been corrected.

In the alarm descriptions, read the *Possible cause* information and perform corrective actions, if applicable.



WARNING!

If an alarm condition cannot be corrected, discontinue use and refer the ventilator for service.

6.2 Operator's Position

To receive the audible part of an alarm, the operator's position should be within audible range from the ventilator, depending on the set audible alarm level.

To receive the visual part of an alarm and its priority, the operator's position should be within a distance of 4 metres (13 feet) from the ventilator, and within an angle of 30° to the normal of the ventilator display.

6.3 Physiological Alarms

The physiological alarms of the ventilator are related to the treatment parameters of the ventilator.

6.3.1 High Flow Alarm

| Property | Description |
|--------------------|---|
| Alarm text | High Flow Alarm |
| Priority | High |
| Alarm condition | A high flow alarm will be given when the total flow exceeds the set High Flow alarm limit for 3 consecutive breaths during inspiration. The alarm is reset after a full breath with a flow below the alarm limit. |
| Possible cause | <ul style="list-style-type: none">• Unintentional leaks from the patient interface or breathing circuit• Mismatch between IPAP/CPAP and alarm setting• Coughing during inspiration• Changes in airway resistance and / or compliance |
| Setting range | <ul style="list-style-type: none">• 10 l/min to 200 l/min• Off |
| Setting resolution | 5 l/min |
| Setting display | The alarm setting is also displayed by a red line in the Flow bar graph. |
| Ventilator action | TheNippy 4+ will continue treatment with the current settings, and try to compensate for unintentional leakages. |

6.3.2 Low Flow Alarm

| Property | Description |
|--------------------|--|
| Alarm text | Low Flow Alarm |
| Priority | High |
| Alarm condition | A low flow alarm will be given when the total flow remains below the Low Flow alarm setting for more than 10 seconds or during 3 consecutive breaths. The alarm is reset when the flow exceeds the alarm limit again. |
| Possible cause | <ul style="list-style-type: none">• An Obstruction in the breathing circuit• Mismatch between IPAP/CPAP and alarm setting• An obstructed CO₂ leak valve, reducing the total flow• Changes in airway resistance and or compliance |
| Setting range | <ul style="list-style-type: none">• 3 l/min to 180 l/min• Off |
| Setting resolution | Below 10 l/min: 1 l/min Above 10 l/min: 5 l/min |
| Setting display | The alarm setting is displayed by a red line in the Flow bar graph. |
| Ventilator action | The Nippy 4+ will continue treatment with the current settings. |

6.3.3 High Pressure Alarm

| Property | Description |
|--------------------|--|
| Alarm text | High Pressure |
| Priority | High |
| Alarm condition | A High Pressure alarm will be given when the patient pressure reaches the set High Pressure alarm limit for three consecutive breaths. It will also be given if pressure exceeds 76 cmH ₂ O. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between pressure setting and alarm setting.• Coughing during inspiration.• Changes in airway resistance and or compliance. |
| Reset criteria | A full breath is performed with maximum pressure below the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. The actual breath is however terminated if the High Pressure alarm limit is reached. |
| Setting range | <ul style="list-style-type: none">• 5 cmH₂O to 70 cmH₂O <p>Note that the High pressure alarm cannot be set lower than the value set for the Low pressure alarm.</p> |
| Setting resolution | Below 10 cmH ₂ O: 0.5 cmH ₂ O Above 10 cmH ₂ O: 1.0 cmH ₂ O |
| Setting Display | The High Pressure alarm setting is displayed by a red line in the pressure bar graph. |

6.3.4 Low Pressure Alarm

| Property | Description |
|--------------------|--|
| Alarm text | Low Pressure |
| Priority | High |
| Alarm condition | A Low Pressure alarm will be given when the Nippy 4+ pressure fails to reach the low pressure alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Disconnection of patient circuit.• Mismatch between pressure setting and alarm setting.• Leakage from the mask or other components of the patient circuit. |
| Reset criteria | The pressure rises above the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• 1 cmH₂O to 50 cmH₂O <p>Note that the Low pressure alarm cannot be set higher than the value set for IPAP or the High pressure alarm.</p> |
| Setting resolution | Below 10 cmH ₂ O: 0.5 cmH ₂ O Above 10 cmH ₂ O: 1.0 cmH ₂ O |
| Setting display | The Low Pressure alarm setting is displayed by a red line in the pressure bar graph. |

6.3.5 High EPAP Alarm

| Property | Description |
|-------------------|---|
| Alarm text | High EPAP |
| Priority | Medium |
| Alarm condition | A High EPAP alarm will be given when the measured EPAP is 30% above the set value for more than 15 seconds |
| Possible cause | <ul style="list-style-type: none">• Blocked leakage port.• Too short expiratory time.• Changes in airway resistance and or compliance.• Malfunction of the exhalation valve.• Blocked exhalation valve. |
| Reset criteria | EPAP has gone below the alarm limit (lower than 30% above the set value). |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• On• Off |

6.3.6 Low EPAP Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Low EPAP |
| Priority | Medium |
| Alarm condition | A Low EPAP alarm will be given when the measured EPAP is 30% below the set value for more than 60 seconds |
| Possible cause | <ul style="list-style-type: none">• Excessive leakage.• Malfunction of the exhalation valve. |
| Reset criteria | EPAP has gone above the alarm limit (higher than 30% below the set value). |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• On• Off |

6.3.7 High V_{t_i} (High Inspired Tidal Volume Alarm)

| Property | Description |
|--------------------|--|
| Alarm text | High V_{t_i} |
| Priority | Medium |
| Alarm Condition | A High Inspired Tidal Volume alarm will be given when the monitored Inspired Tidal Volume exceeds the set limit for the High Inspired Tidal Volume alarm for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Inspired Tidal Volume and alarm setting.• Pressure settings causing the Inspired Tidal Volume to exceed the set alarm level.• Leakage from the mask or other components of the patient circuit.• Mismatch between selected and used patient circuit. |
| Reset criteria | When inspired tidal volume is below set alarm limit |
| Setting range | <ul style="list-style-type: none">• 100 ml to 2500 ml• Off |
| Setting resolution | 10 below 600 ml, 100 above 600 ml |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.8 Low V_{t_i} Alarm (Low Inspired Tidal Volume)

| Property | Description |
|--------------------|---|
| Alarm text | Low V_{t_i} |
| Priority | High |
| Alarm Condition | A Low Inspired Tidal Volume alarm will be given when the monitored Inspired Tidal Volume fails to reach the set limit for the Low Inspired Tidal Volume alarm for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Inspired Tidal Volume and Alarm setting.• Changes in airway resistance and or compliance. |
| Setting range | <ul style="list-style-type: none">• 50 ml to 2000 ml• Off |
| Setting resolution | 10 below 600 ml, 100 above 600 ml |
| Reset criteria | A full breath above set alarm limit |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.9 High MV_i Alarm (High Inspired Minute Volume Alarm)

| Property | Description |
|--------------------|--|
| Alarm text | High MV_i |
| Priority | Medium |
| Alarm condition | A High Inspired Minute Volume alarm will be given when the monitored inspired minute volume exceeds the set limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Breath Rate, Inspired Tidal Volume settings and the alarm setting.• Increased Breath Rate.• Leakage around the mask or within one of the components of the circuit. |
| Reset criteria | When inspired minute volume is below set alarm limits |
| Setting range | <ul style="list-style-type: none">• 1.0 to 40 l/min• Off |
| Setting resolution | 0.5 l/min |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.10 Low MV_i (Low Inspired Minute Volume Alarm)

| Property | Description |
|--------------------|---|
| Alarm text | Low MV_i |
| Priority | High |
| Alarm condition | A Low Inspired Minute Volume alarm will be given when the monitored minute volume does not reach the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Breath Rate and Inspired Tidal Volume settings and the alarm setting.• Changes in airway resistance and or compliance.• Decreased Breath Rate. |
| Setting range | <ul style="list-style-type: none">• 1.0 l/min to 30 l/min• Off |
| Setting resolution | 0.1 l up to 1.0 l, 0.5 l above 1.0 l. |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.11 High V_{t_e} Alarm (High Expired Tidal Volume)

| Property | Description |
|--------------------|--|
| Alarm text | High Vte |
| Priority | Medium |
| Alarm condition | A High Expired Tidal Volume alarm will be given when the monitored Expired Tidal Volume exceeds the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Expired Tidal Volume and alarm setting.• Mismatch between selected and used patient circuit. |
| Setting range | <ul style="list-style-type: none">• 100 ml to 2500 ml• Off |
| Setting resolution | 10 below 600 ml, 100 above 600 ml |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.12 Low V_{t_e} Alarm (Low Expired Tidal Volume)

| Property | Description |
|--------------------|--|
| Alarm text | Low Vte |
| Priority | High |
| Alarm Condition | A Low Expired Tidal Volume alarm will be given when the monitored Expired Tidal Volume fails to reach the set limit for the Low Expired Tidal Volume alarm for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Expired Tidal Volume and Alarm setting.• Changes in airway resistance and or compliance.• Leakage around the mask or within one of the components of the circuit. |
| Reset criteria | Full breath above set alarm limit |
| Setting range | <ul style="list-style-type: none">• 50 ml to 2000 ml• Off |
| Setting resolution | 10 below 600 ml, 100 above 600 ml |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.13 High MV_e (High Expired Minute Volume Alarm)

| Item | Description |
|--------------------|---|
| Alarm text | High MVe |
| Priority | Medium |
| Alarm condition | A High Expired Minute Volume alarm will be given when the monitored expired minute volume exceeds the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Breath Rate, Tidal Volume settings and the alarm setting.• Increased Breath Rate. |
| Setting range | <ul style="list-style-type: none">• 1.0 to 40 l/min• Off |
| Setting resolution | 0.5 l/min |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.14 Low MV_e Alarm (Low Expired Minute Volume)

| Property | Description |
|--------------------|--|
| Alarm text | Low MVe |
| Priority | High |
| Alarm condition | A Low Expired Minute Volume alarm will be given when the monitored minute volume is below the alarm limit for more than 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between Breath Rate and Tidal Volume settings and the alarm setting.• Changes in airway resistance and or compliance.• Decreased Breath Rate.• Leakage around the mask or within one of the components of the circuit. |
| Setting range | <ul style="list-style-type: none">• 1.0 l/min to 30 l/min• Off |
| Setting resolution | 0,5 l |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.15 High Breath Rate Alarm

| Property | Description |
|--------------------|--|
| Alarm text | High Breath Rate |
| Priority | Medium |
| Alarm condition | A High Breath Rate alarm will be given when the alarm limit has been exceeded for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between the Breath Rate setting and the alarm setting.• Increased Breath Rate.• Too sensitive setting of the inspiratory trigger setting. |
| Reset criteria | The breath rate goes below the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• 10 bpm to 70 bpm• Off |
| Setting resolution | 1 bpm. |

6.3.16 Low Breath Rate Alarm

| Property | Description |
|--------------------|---|
| Alarm text | Low Breath Rate |
| Priority | High |
| Alarm condition | A Low Breath Rate alarm will be given when the delivered total breath rate is below the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Mismatch between the Breath Rate setting and the alarm setting.• The patient cannot trigger breaths because the inspiratory trigger setting is too high.• Decrease in the patient's spontaneous breathing.• Circuit disconnection. |
| Reset criteria | The breath rate goes above the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• 4 bpm to 40 bpm (non-MPV modes)• 1 bpm to 40 bpm (MPV modes)• Off |
| Setting resolution | 1 bpm. |

6.3.17 Apnoea Alarm

| Property | Description |
|--------------------|---|
| Alarm text | Apnoea |
| Priority | High |
| Alarm condition | An Apnoea alarm will be given when no patient-triggered breath is detected for the set period of time. The Apnoea alarm is only available if the Inspiratory trigger is activated. |
| Possible cause | <ul style="list-style-type: none">• Patient stopped breathing.• Patient decreases spontaneous breathing.• Circuit disconnection.• Inspiratory Trigger is set too high. |
| Reset criteria | Inspiratory effort detected by the Nippy 4+. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• 5 to 60 s.(non-MPV modes)• 15 to 900 s. (MPV modes)• Off |
| Setting resolution | 5 s below 15 s. 15 s above 15 s. MPV modes: 15 s below 60 s. 60 s above 60 s. |

6.3.18 Disconnection Alarm



CAUTION!

No single alarm can reliably detect all disconnections due to the number of possible combinations of therapy settings, circuit configurations and patient interfaces. To verify that patient disconnection can be detected, including if the patient interface becomes accidentally detached from the patient, it is advised to test the functionality of the Disconnection Alarm upfront with the complete set-up as used during treatment, including items such as filters, circuit, connectors, interface (mask, cannula etc.)

| Property | Description |
|-------------------|--|
| Alarm text | Disconnection |
| Priority | High |
| Alarm condition | A Disconnection alarm will be given when the measured flow exceeds the expected leakage flow at the set Pressure for 15 seconds. |
| Possible cause | <ul style="list-style-type: none"> • Too high leakage in the patient circuit. • The patient has removed the mask. • Circuit disconnection. • Pilot pressure tube disconnection |
| Reset criteria | The leakage is back within limits. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings |
| Setting range | <ul style="list-style-type: none"> • On • Off |

6.3.19 Rebreathing Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Rebreathing (with leakage circuit) Rebreathing (with active exhalation valve circuit) |
| Priority | High (with leakage circuit) Medium (with active exhalation valve circuit) |
| Alarm condition | <p>Leakage Circuit A Rebreathing alarm will be given if the intentional leakage is too low for more than 15 seconds.</p> <p>Exhalation valve circuit A Rebreathing alarm will be given if the exhalation valve is obstructed for more than 10 consecutive breaths.</p> <p>MPV circuit A Rebreathing alarm will be given if air returns into the ventilator for more than 10 consecutive breaths.</p> |
| Possible cause | <ul style="list-style-type: none"> • Obstructed or occluded patient circuit. • Incorrect patient circuit. • Patient exhales through mouthpiece. • Obstructed or removed CO₂ port from leakage circuit. • Disconnected pilot pressure line from rear of device or from the active exhalation valve. |
| Reset criteria | The leakage is back within limits. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none"> • On • Off |

6.3.20 Obstruction Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Obstruction |
| Priority | High |
| Alarm condition | An Obstruction alarm will be given if the inspiratory breathing tube becomes blocked and remains blocked for 2 consecutive breaths. |
| Ventilator action | With each breath cycle, upon detection of an obstruction the ventilator will reduce the airway pressure to the set EPAP. Treatment will resume with the start of the next breath cycle. |
| Reset Criteria | When the monitored compliance and resistance become normal after a breath. |
| Setting Range | <ul style="list-style-type: none">• High• Low• Off |

6.3.21 High FiO₂ Alarm

| Property | Description |
|--------------------|--|
| Alarm text | High FiO₂ |
| Priority | Medium |
| Alarm condition | A High FiO ₂ alarm will be given when the measured FiO ₂ exceeds the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Increased oxygen inflow.• Decreased minute ventilation. |
| Reset criteria | FiO ₂ goes below the alarm limit |
| Setting range | <ul style="list-style-type: none">• 21% to 100%• Off |
| Setting resolution | 1% |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.22 Low FiO₂ Alarm

| Property | Description |
|--------------------|--|
| Alarm text | Low FiO₂ |
| Priority | High |
| Alarm condition | A Low FiO ₂ alarm will be given when the measured FiO ₂ is below the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Decreased oxygen inlet.• Disconnection of oxygen inlet.• Increased minute ventilation.• High leakage. |
| Setting range | <ul style="list-style-type: none">• 21% to 100%• Off |
| Setting resolution | 1% |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.3.23 High SpO₂ Alarm

| Property | Description |
|--------------------|--|
| Alarm text | High SpO₂ |
| Priority | Medium |
| Alarm condition | A High SpO ₂ alarm will be given when the measured SpO ₂ exceeds the alarm limit for 30 seconds. |
| Possible cause | Too high flow of bleed-in oxygen. |
| Reset criteria | The SpO ₂ value goes back below the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | <ul style="list-style-type: none">• 90 % to 100 %• Off |
| Setting resolution | 1 % |

This alarm requires a connected SpO₂ sensor.

6.3.24 Low SpO₂ Alarm

| Property | Description |
|--------------------|---|
| Alarm text | Low SpO₂ |
| Priority | High |
| Definition | A Low SpO ₂ alarm will be given when the measured SpO ₂ is below the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Too low flow of bleed-in oxygen.• Oxygen inlet is disconnected.• Delivered tidal volumes are too small. |
| Setting range | 85% to 100% |
| Setting resolution | 1% |
| Ventilator action | The ventilator will continue treatment with the same settings. |

This alarm requires a connected SpO₂ sensor.

6.3.25 High EtCO₂ Alarm

| Property | Description |
|--------------------|---|
| Alarm text | High EtCO₂ |
| Priority | High |
| Alarm condition | A High EtCO ₂ alarm will be given when the measured EtCO ₂ exceeds the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Alarm limit is set too low.• Breath Rate too low.• Delivered Tidal Volume too low.• Excessive dead space between patient and leakage port.• Leak port occluded. |
| Setting range | 1 to 99 mmHg Off |
| Setting resolution | 1 mmHg |
| Ventilator action | The ventilator will continue treatment with the same settings. |

This alarm requires a connected EtCO₂ sensor

6.3.26 Low EtCO₂ Alarm

| Property | Description |
|--------------------|--|
| Alarm text | Low EtCO₂ |
| Priority | Medium |
| Alarm condition | A Low EtCO ₂ alarm will be given when the measured EtCO ₂ is below the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Alarm limit is set too high.• Ventilator disconnection.• Excessive leakage in the Patient circuit/Interface.• Partial obstruction of the airways.• Breath Rate too high.• Delivered Tidal Volume too high.• Self triggering of the ventilator. |
| Setting range | 1 to 99 mmHg Off |
| Setting resolution | 1 mmHg |
| Ventilator action | The ventilator will continue treatment with the same settings. |

This alarm requires a connected EtCO₂ sensor

6.3.27 High InspCO₂ Alarm (High Inspired CO₂)

| Property | Description |
|--------------------|--|
| Alarm text | High InspCO₂ |
| Priority | High |
| Alarm condition | A High Inspired CO ₂ alarm will be given when the measured inspired CO ₂ exceeds the alarm limit for 30 seconds. |
| Possible cause | <ul style="list-style-type: none">• Alarm limit is set too low.• Excessive dead space between patient and exhalation valve/leakage port.• Leakage port/valve occluded. |
| Setting range | 1 to 99 mmHg Off |
| Setting resolution | 1 mmHg |
| Ventilator action | The ventilator will continue treatment with the same settings. |

This alarm requires a connected EtCO₂ sensor

6.3.28 High Pulse Rate Alarm

| Property | Description |
|--------------------|---|
| Alarm text | High Pulse Rate |
| Priority | Medium |
| Alarm condition | A High Pulse Rate alarm will be given when the measured pulse rate exceeds the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Insufficient ventilatory support.• Too low flow of bleed-in oxygen.• The EPAP value is set too high.• Bad positioning of the finger probe. |
| Reset criteria | The pulse rate goes back below the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | 30 to 230 bpm (beats per minute) Off |
| Setting resolution | 5 bpm (beats per minute) |

This alarm requires a connected SpO₂ sensor.

6.3.29 Low Pulse Rate Alarm

| Property | Description |
|--------------------|---|
| Alarm text | Low Pulse Rate |
| Priority | High |
| Alarm condition | A low pulse rate alarm will be given when the measured pulse rate goes below the alarm limit for 15 seconds. |
| Possible cause | <ul style="list-style-type: none">• Bad positioning of the finger probe.• Too low flow of bleed-in oxygen.• Insufficient ventilatory support. |
| Reset criteria | The pulse rate goes back above the alarm limit. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |
| Setting range | 30 to 230 bpm (beats per minute) Off |
| Setting resolution | 5 bpm (beats per minute) |

This alarm requires a connected SpO₂ sensor.

6.3.30 PtcCO₂ Alarm

| Property | Description |
|-------------------|--|
| Alarm text | PtcCO₂ Outside Limits |
| Priority | High |
| Alarm Condition | A PtcCO ₂ alarm will be given when PtcCO ₂ is outside alarm limits. Check the PtcCO ₂ monitor. |
| Possible cause | <ul style="list-style-type: none">• External PtcCO₂ monitor is outside its alarm limits.• Breath Rate needs adjustment• Delivered Tidal Volume needs adjustment.• Excessive dead space between patient and exhalation valve/leakage port.• Leak port/valve occluded.• Ventilator disconnection.• Excessive leakage in the Patient circuit/Interface.• Partial obstruction of the airways.• Self triggering of the ventilator. |
| Ventilator action | The ventilator will continue treatment with the same settings. |

This alarm requires a connected PtcCO₂ sensor.

6.4 Technical Alarms

6.4.1 Power Fail Alarm

| Property | Description |
|-------------------|---|
| Alarm text | The alarm is given audibly with a tone and the display is blinking with the alarm message Power Fail |
| Priority | High |
| Alarm condition | The Power Fail alarm is given if the last power source fails to provide enough power for running the ventilator. |
| Possible cause | The last available power source cannot deliver power to the ventilator. Battery discharged or battery failure. |
| Reset criteria | External power supply connected to ventilator. |
| Ventilator action | The Nippy 4+ stops the treatment, and gives the Power Fail alarm for at least 2 minutes. If power is restored within the alarm time, the ventilator will automatically resume treatment with current settings. When powered up again, the power failure will be logged. |

6.4.2 High Patient Air Temp. (High Patient Air Temperature)

| Property | Description |
|-------------------|--|
| Alarm text | High Patient Air Temp |
| Priority | High |
| Alarm condition | A High Patient Air Temperature alarm will be given when the patient air temperature exceeds 43°C (109.4°F). |
| Possible cause | <ul style="list-style-type: none">Blocked air inlets.Blocked cooling air outlets.Too high ambient temperature. |
| Ventilator action | The ventilator stops treatment and gives alarm for up to 2 minutes. If power is restored within the alarm time, the ventilator will automatically resume treatment with current settings. |
| Reset criteria | The temperature goes below the limit again. |

6.4.3 Low Patient Air Temp. (Low Patient Air Temperature Alarm)

| Property | Description |
|-------------------|---|
| Alarm text | Low Patient Air Temp |
| Alarm condition | A Low Patient Air Temperature alarm will be given when the patient air temperature is below the preset limit -30°C (-22°F). |
| Priority | High |
| Possible cause | Too low ambient temperature |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.4.4 Low Last Power Source Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Low Last Power Source |
| Priority | Medium |
| Alarm condition | This alarm will be given when the last battery source (internal battery) has 15 minutes of operating time left with current settings. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. |

6.4.5 Crit. Low Last Power Source Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Crit. Low Last Power Source |
| Alarm condition | A Crit. Low Last Power Source alarm will be given when the last battery source (internal battery or click-in battery) has 5 minutes of operating time left with current settings. |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | Connection of “higher” power source. |

6.4.6 Lost Mains Power Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Lost Mains Power |
| Alarm condition | A Mains Power Lost alarm will be given when the ventilator switched from Mains to another power source due to Mains power is lost. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. An information message will be shown on the screen. |
| Reset | Confirmation by user or mains reconnected. |

6.4.7 Exhalation Valve Control Error Alarm

| Property | Description |
|-----------------|---|
| Alarm text | Exhalation Valve Control Error |
| Alarm condition | An Exhalation Valve Control Error alarm will be given when the ventilator fails to control the internal /external exhalation valve. |
| Priority | High |
| Possible cause | <ul style="list-style-type: none">• Exhalation valve occluded• Exhalation valve control tube disconnected• Insert not properly mounted• Internal function failure of the exhalation valve controls |
| Reset | The pilot pressure gets a normal value. |

6.4.8 SpO₂ Disconnected (SpO₂ Sensor Failure/Disconnection Alarm)

| Property | Description |
|-------------------|---|
| Alarm text | SPO2 Disconnected |
| Alarm condition | An SpO ₂ Sensor Failure/Disconnection alarm will be given when an error signal or no signal from the SpO ₂ sensor has been detected for 2 seconds. Check the SpO ₂ sensor. |
| Priority | High |
| Possible cause | The SpO ₂ electronics cable has been disconnected and subsequently no communication (possibly due to disconnection) for 2 seconds. Failure in the SpO ₂ sensor. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | Confirmation by user or reconnected/changed. |

6.4.9 SpO₂ Signal Lost Alarm

| Property | Description |
|-------------------|---|
| Alarm text | SPO2 Signal Lost |
| Alarm condition | SpO ₂ signal lost. |
| Priority | High |
| Possible cause | Signal lost reported by SpO ₂ electronics (due to patient removing the probe from finger, or sensor detached from SpO ₂ electronics). |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | User presses OK or electronics cable is disconnected by the user, or the sensor is reconnected to the finger. |

6.4.10 Poor SpO₂ Signal

| Property | Description |
|-------------------|---|
| Alarm text | Poor SPO2 Signal |
| Alarm condition | A Poor SpO ₂ signal alarm will be given when the SpO ₂ signal is not correct. Check the SpO ₂ sensor. |
| Priority | High |
| Possible cause | Artifact or low perfusion reported by SpO ₂ electronics |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | OK message from SpO ₂ electronics or SpO ₂ electronics disconnected by user or SpO ₂ Signal Lost alarm is triggered. |

6.4.11 CO₂ Disconnected (CO₂ Sensor Failure/Disconnection Alarm)

| Property | Description |
|-------------------|--|
| Alarm text | CO2 Sensor Disconnected |
| Alarm condition | A CO ₂ Sensor Failure/Disconnection alarm will be given when communication between the ventilator and the CO ₂ sensor has been lost for 2 seconds. Check the CO ₂ sensor. |
| Priority | High |
| Possible cause | <ul style="list-style-type: none">• CO₂ Sensor disconnected.• Failure in the CO₂ sensor. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | Confirmation by user or reconnected/changed. |

6.4.12 CO₂ Accuracy Error Alarm

| Property | Description |
|-------------------|--|
| Alarm text | CO2 Accuracy Error |
| Alarm condition | A CO ₂ Accuracy Error alarm will be given when an accuracy error in the CO ₂ measurement has occurred. |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | OK message from sensor or sensor disconnected by user. |

6.4.13 Check CO₂ Adapter Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Check CO2 Adapter |
| Alarm condition | A Check CO ₂ Adapter alarm will be given when the airway adapter is not attached correctly to the CO ₂ sensor. Check/replace the airway adapter. |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | OK message from sensor or sensor disconnected by user. |

6.4.14 CO₂ Sensor Error Alarm

| Property | Description |
|-------------------|---|
| Alarm text | CO2 Sensor Error |
| Alarm condition | A CO ₂ Sensor Error alarm will be given when an error in the CO ₂ sensor has occurred. Replace the CO ₂ sensor. CO ₂ monitoring cannot be performed in this condition. |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | OK message from sensor or sensor disconnected by user. |

6.4.15 FiO₂ Disconnected (FiO₂ Sensor Failure/Disconnection Alarm)

| Property | Description |
|-------------------|--|
| Alarm text | FiO2 Disconnected |
| Alarm condition | An FiO ₂ Sensor Failure/Disconnection alarm will be given when no signal from the FiO ₂ sensor has been detected for 2 seconds. Check the FiO ₂ sensor. |
| Priority | High |
| Possible cause | <ul style="list-style-type: none">• FiO₂ Sensor disconnected.• Communication with the FiO₂ sensor failed. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | Confirmation by user or reconnected/changed. |

6.4.16 Ambient Pressure Compensation Lost Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Pressure Comp Lost |
| Priority | Medium |
| Alarm condition | An Ambient Pressure Compensation Lost alarm will be given when the automatic ambient pressure compensation functionality is out of order. |
| Ventilator action | The Nippy 4+ will continue treatment according to the current settings. Normal atmospheric pressure at sea level will be used as approximation for the temporary ambient pressure compensation. If used at other altitude, delivered and measured pressures may deviate. |
| Reset | Reset of ventilator. |

6.4.17 Temperature Comp. Lost (Ambient Temperature Compensation Lost Alarm)

| Property | Description |
|-------------------|---|
| Alarm text | Temperature Comp. Lost |
| Alarm condition | An Ambient Temperature Compensation Lost alarm will be given when the automatic ambient temperature compensation is out of order. There is no communication with the air temperature sensor or the value is out of range (less than -30°C (-22°F) or more than 70°C (158°F). |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. The accuracy of the volume measurement may be impaired. |
| Reset | Ambient temperature inside valid range. |

6.4.18 Humidity Comp. Lost (Humidity Compensation Lost Alarm)

| Property | Description |
|-------------------|---|
| Alarm text | Humidity Comp. Lost |
| Alarm condition | An Humidity Compensation Lost alarm will be given when the automatic humidity compensation is out of order. 50% relative humidity is used for temporary compensation. If the ventilator is used at other humidities, delivered and measured pressure and flow may deviate. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. The accuracy of the volume measurement may be impaired. |
| Reset | Air humidity sensor values (RH and temperature) inside valid range. |

6.4.19 LED Failure Alarm

| Property | Description |
|-------------------|--|
| Alarm text | LED Failure |
| Alarm condition | A LED Failure alarm will be given when one or more LED indicators on the front panel are broken. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | Power-on reset of ventilator (or repair). |


6.4.20 Low Alarm Battery Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Low Alarm Battery |
| Alarm condition | An alarm for <i>Low Alarm Battery</i> will be given if the alarm battery is not charged enough to have power for a <i>Power Fail</i> alarm for at least 2 minutes. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings and start charging the alarm batteries. |
| Reset | When alarm energy storage level is sufficient to give an alarm for at least 2 minutes. |

6.4.21 Alarm Battery Error Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Alarm Battery Error |
| Alarm condition | Unable to communicate with super capacitor and read super capacitor status. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | When triggering condition is removed. |

6.4.22 Internal/Click-In Battery Hot Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Internal Battery — Internal Battery Hot Click-In Battery — Click-In Battery Hot |
| Alarm condition | An alarm for Internal/Click-In Battery Overheat in Discharge will be given when the internal or click-in battery reaches 55° C (131°F). |
| | <div style="border-left: 2px solid black; padding-left: 10px;"> <p>NOTE</p> <p> The battery electronics by manufacture stops discharge at 60°C (140°F).</p> </div> |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. Battery discharging will be disabled (by the battery electronics) once the temperature gets to 60°C (140°F). (If the battery is last power source, the ventilator will stop running). |

6.4.23 Heated Circuit Temp. Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Heated Circuit Temp. |
| Alarm condition | A Heated Circuit temp alarm will be given when the measured temperature of the heated wire is outside the tolerance. |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Priority | Medium |
| Reset | Heated wire measured temp tolerance is inside limits. |

6.4.24 High Humidifier Temp. Alarm

| Property | Description |
|-------------------|--|
| Alarm text | High Humidifier Temp. |
| Alarm condition | A Humidifier High Temperature alarm will be given if the humidifier heater plate temperature exceeds 76°C (169°F) for more than 2 seconds. |
| Priority | Medium |
| Ventilator action | The ventilator will turn off the click-in humidifier and then continue treatment with the same settings. A message with option to turn on the humidifier again will be displayed. |
| Reset | The alarm is dismissed when the humidifier temperature drops below 76°C (169°F), set humidifier temperature). |

6.4.25 Humidifier Fault Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Humidifier Fault |
| Alarm condition | <ul style="list-style-type: none">• All humidifier enabling conditions have been satisfied for 10 minutes, and• No humidifier setting changes have been made for 10 minutes, and• Heater plate temperature < 50°C (122 °F)• Humidifier set temperature > Ambient temperature ,and• The heater plate temperature is more than 5°C (41 °F) below the set temperature, or the heater plate temperature < -20°C (68 °F) or greater than 400°C (752 °F) |
| Priority | Medium |
| Ventilator action | The ventilator will turn off the humidifier and continue treatment with the same settings. The humidifier must be restarted manually when the cause of the alarm is resolved. |

6.4.26 Heated Circuit Fault Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Heated Circuit Fault |
| Alarm condition | A Heated Circuit Fault alarm will be given if a fault in the heated circuit electronics or temperature sensor is detected. |
| Priority | Medium |
| Ventilator action | The ventilator will turn off the heated circuit and continue treatment with the same settings. The heated circuit must be restarted manually when the cause of the alarm is resolved. |
| Reset | The alarm is dismissed when the heated circuit setting is changed to OFF, or the treatment is stopped. The power to the heated circuit is re-enabled when all enabling conditions are satisfied. |

6.4.27 Internal Function Failure

| Property | Description |
|-------------------|--|
| Alarm text | Int. Function Failure |
| Priority | High |
| Alarm condition | Failure of internal function that prevents treatment or normal operation of the ventilator. The error code that follows the alarm text indicates the kind of function failure. All Internal Function Failure alarm error codes are defined and explained in the ventilator Service Manual. |
| Reset criteria | Restart the ventilator. |
| Ventilator action | The ventilator will stop the treatment and shut down. |
| Action to take | Restart the Nippy 4+. If the alarm persists or reoccurs: Take a note of the error code and contact your supplier of the Nippy 4+ . |

6.4.28 Air Temp. Sensor Fail Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Air Temp Sensor Fail |
| Alarm condition | The alarm is given in case of swivel boot temperature sensor communication failure or sensor reporting temperatures out of range (below -30°C (-22°F) or above 60°C (140°F)). |
| Priority | Medium |
| Ventilator action | The ventilator will continue treatment with the same settings. |

6.4.29 Internal Error Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Internal Error |
| Priority | High |
| Alarm Condition | An internal Error alarm will be given when the ventilator has an internal error, followed by an error code for the specific failure. All Internal Error alarm error codes are defined and explained in the ventilator Service Manual. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset action | Power off and restart the ventilator. |

6.4.30 Database Integrity Fail Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Database Integrity Failed |
| Priority | High |
| Alarm Condition | This alarm is given when the database integrity check fails. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset action | Rebuild the database and restart the ventilator. |

6.4.31 Cooling Fan Error Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Cooling Fan Error |
| Alarm Condition | The Cooling Fan Error alarm shall be given when the cooling fan runs too slow. |
| Priority | High |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | When cooling fan speed is above 275 rpm. |

6.4.32 Clock Failure Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Clock Failure |
| Priority | High |
| Alarm condition | The alarm shall be given when the real time clock value is invalid. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset action | Restart the ventilator. |

6.4.33 Internal Temp High Alarm

| Property | Description |
|-------------------|--|
| Alarm text | Internal Temp High |
| Priority | High |
| Alarm condition | The Internal High Temp alarm shall be given when the ventilator internal temperature is high. The internal temp high alarm is triggered when PTU/Sensor board temperature is higher than 65°C (149°F), or main board temperature is higher than 65°C (149°F), or motor temperature is higher than 85°C (185°F). |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset | When the triggering conditions are resolved. |

6.4.34 Humidifier/Bypass Loose Alarm

| Property | Description |
|-------------------|---|
| Alarm text | Humidifier/Bypass Loose |
| Priority | Medium |
| Alarm condition | The Humidifier/Bypass Loose alarm shall be given when the air bypass/humidifier latch is stuck in the down position for 5 secs. |
| Ventilator action | The ventilator will continue treatment with the same settings. |
| Reset action | Reinsert the air bypass unit/humidifier and make sure the latch closes. |

6.5 Alarm Test

6.5.1 Alarm Signal Test

When starting treatment, an automatic alarm signal test is performed. Check that the test is performed successfully, this is indicated by:

- A short beep indicating functional audio signaling.
- The alarm LED first lights yellow, then red, indicating functional visual signaling.
- The audio pause LED lights yellow.
- In about a second, both LEDs are turned off.

If the test fails, do not use the Nippy 4+. Contact your supplier of the Nippy 4+ for a technical check.

6.5.2 Mandatory Alarm Tests

This alarm test should be performed when changing patient, if the ventilator's function needs to be checked for any other reason, or at least every 12 months.

The alarm test should be included in the regular inspections during maintenance.

To perform the alarm test, follow the instructions below:

Alarm Test Preparation

- 1 Connect the ventilator patient circuit to a test lung.
- 2 Connect the ventilator to Mains power supply.
- 3 Start the ventilator.
- 4 Adjust the settings as follows:

| Setting | Value |
|------------------|------------------------------------|
| Ventilation Mode | Pressure Support Ventilation (PSV) |
| Patient Mode | Adult |
| IPAP | 15 cmH ₂ O |
| EPAP | 5 cmH ₂ O |
| Rise Time | 9 |
| Insp. Trigger | 9 |
| Exp. Trigger | 3 |
| Min Insp. Time | Off |

| | |
|-------------------|--------|
| Max Insp. Time | Off |
| Backup Rate | 12 bpm |
| Backup Insp. Time | 2.0 s |
| Target Volume | Off |

5 All alarm settings shall be set to Off if possible.

6 Start the treatment.

6.5.2.1 High and Low Flow Alarm Tests

- 1 Set the high flow alarm to 20 l/min.
⇒ The high flow alarm shall be given.
- 2 Set the high flow alarm to Off
- 3 Set the low flow alarm to 150L/min.
⇒ The low flow alarm shall be given.

6.5.2.2 High and Low Pressure Alarm Tests

- 1 Set the high pressure alarm to 10 cmH₂O.
⇒ The high pressure alarm shall be given.
- 2 Set the high pressure alarm to 55 cmH₂O.
- 3 Set the low pressure alarm to 20 cmH₂O.
⇒ The low pressure alarm shall be given.
- 4 Set the low pressure alarm to 1.0 cmH₂O.

6.5.2.3 Expiratory Tidal Volume Alarm (V_{t_e}) Tests

This alarm test applies if having a patient circuit with intentional leakage.

- 1 Set up the ventilator as described in *Alarm Test Preparation*, page 134.
- 2 Set the high V_{t_e} alarm to 150 ml.
⇒ The high V_{t_e} alarm shall be given.
- 3 Set the high V_{t_e} alarm to Off.
- 4 Set the low V_{t_e} alarm to 400 ml.
The low V_{t_e} alarm shall be given.

6.5.2.4 Inspiratory Tidal Volume Alarm (V_t) Tests

This alarm test applies if having a patient circuit with exhalation valve or a patient circuit with mouthpiece.

- 1 Set up the ventilator as described in *Alarm Test Preparation*, page 134.
- 2 Set the high V_t alarm to 150 ml.
⇒ The high V_t alarm shall be given.
- 3 Set the high V_t alarm to Off.
- 4 Set the low V_t alarm to 400 ml.
The low V_t alarm shall be given.

6.5.2.5 EtCO₂ Related Alarm Test

This alarm test applies if the EtCO₂ accessory is used.

- 1 Connect the EtCO₂ sensor with an attached airway adapter to the Nippy 4+.
- 2 Disconnect the airway adapter from the CO₂ sensor.
⇒ The check CO₂ adapter alarm shall be given.
- 3 Connect the airway adapter to the CO₂ sensor again.

6.5.2.6 SpO₂ Related Alarm Tests

These tests apply if the SpO₂ accessory is used,

- 1 Connect SpO₂ sensor to device and to your finger.
- 2 Set the low SpO₂ alarm to 85%.
- 3 Set the high SpO₂ alarm to be 90%.
- 4 Start treatment and wait 30 s.
⇒ High SpO₂ alarm should be given.
- 5 Stop treatment.
- 6 Set the high SpO₂ alarm to off.
- 7 Set the low SpO₂ alarm to be 100%.
- 8 Start treatment and wait 30 s.
⇒ Low SpO₂ alarm should be given.
- 9 Stop Treatment.
- 10 Set the low SpO₂ alarm to 85%.
- 11 Set the low pulse rate alarm to off.
- 12 Set the high pulse rate alarm to 30 bpm.
- 13 Start treatment and wait 30 s.

⇒ High pulse rate alarm should be given.

- 14 Stop treatment.
- 15 Set the high pulse rate alarm to off.
- 16 Set the low pulse rate alarm to be 230 bpm.
- 17 Start treatment and wait 30 s.
⇒ Low pulse rate alarm should be given.
- 18 Stop Treatment.
- 19 Set the low pulse rate alarm to off.

6.5.2.7 Power Related Alarm Tests

- 1 If having the Click-in battery installed, disconnect it.
- 2 Check that the internal battery is fully charged and disconnect the Mains power supply while the treatment is running.
⇒ The Lost mains power alarm shall be given.
- 3 Start a timer and record the time for the following alarms to appear.
 - the Low last power source alarm shall be given after a run time of at least 1 hour and 40 minutes. The ventilator shall continue to run for at least 15 minutes more.
 - When the Critical low last power source alarm is given, the ventilator shall run for at least five more minutes.
 - When the Power Fail alarm is given, the treatment stops. The alarm shall continue to sound for about two more minutes.

The ventilator shall have been able to run at least 2 hours on the internal battery before the Power fail alarm is given

6.5.3 Optional Alarm Tests

In this chapter, methods for additional alarm tests are described. These tests are optional and not needed to ensure safe use of the ventilator.

6.5.3.1 High EPAP Alarm

- 1 Connect the ventilator patient circuit to a test lung and a CPAP device.
- 2 Set the CPAP device treatment pressure to 10 cmH₂O.
- 3 Adjust the ventilator settings as follows:

| Setting | Value |
|------------------|-----------------------|
| Ventilation Mode | Pressure Control |
| IPAP | 15 cmH ₂ O |

| | |
|---------------|----------------------|
| EPAP | 5 cmH ₂ O |
| Breath Rate | 12 bpm |
| Insp. Time | 1.5 s |
| Rise Time | 5 |
| Insp. Trigger | Off |
| Target Volume | Off |

- 4 Start treatment on both the ventilator and the CPAP device.
- 5 Wait approximately 15 seconds before the High EPAP alarm shall be given.
- 6 Stop treatment. Test completed.

6.5.3.2 Low Pressure and Disconnection Alarms

- 1 Start treatment and disconnect the patient circuit.
- 2 Wait 15 seconds.
- 3 The Low Pressure Alarm and/or the Disconnection Alarm will be given.
- 4 Stop treatment. Test completed.

6.5.3.3 Disconnection Alarm Test

- 1 Set the disconnection alarm to On.
- 2 Disconnect the patient circuit.
⇒ The disconnection alarm shall be given.
- 3 Set the disconnection alarm to Off.

6.5.3.4 Obstruction Alarm

- 1 Start treatment; block the patient circuit completely to simulate an obstruction.
- 2 Wait approximately 10 seconds.
- 3 The Obstruction Alarm will be given.
- 4 Stop treatment. Test completed.

7 Cleaning and Maintenance

WARNING!



The Nippy 4+ should be subjected to maintenance, service and control and any applicable upgrades, in accordance with Breas service instructions.



The Nippy 4+ shall only be repaired or modified in accordance with Breas service manuals, technical bulletins, and any special service instructions, by service technicians that have been authorised after Breas Nippy 4+ service training.



Do not under any circumstances attempt to service or repair the ventilator yourself. If you do so, the manufacturer will no longer be responsible for the performance and safety of the ventilator.

Deviation from these service instructions may lead to risk of personal injury!

The patient-connected parts and the filter must be cleaned and replaced regularly to ensure correct function of the ventilator. All replaced parts must be disposed of in accordance with local environmental regulations regarding the disposal of used equipment and waste.

7.1 Cleaning the Nippy 4+

WARNING!



To avoid electrical shock, disconnect the power supply to the ventilator before cleaning. Do not immerse the ventilator into any fluids.

CAUTION!



Always be careful when cleaning to ensure that you do not damage any equipment.



Fluid must not be allowed to enter the ventilator.



Never apply any liquids directly on the ventilator by spraying, splashing or pouring. Use a moistened lint-free cloth when cleaning.



Do not use an excessive amount of liquid when cleaning the ventilator.



Do not autoclave the ventilator.

7.1.1 Main Unit

- 1 Switch off the Nippy 4+ and disconnect the power supply.
- 2 Remove the patient circuit.
- 3 Disconnect all electric cables.
- 4 Clean the outside of the Nippy 4+ using a lint-free cloth with a mild soap solution, and/or ethanol 70% for surface disinfection.
- 5 If the click-in humidifier is used, clean it as described in 5.12.14.4 *Cleaning the Water Chamber*, page 96.
- 6 Reconnect the patient circuit. Make sure all parts are dry before the ventilator is put into operation.

7.1.2 Air Pathway Disinfection

The table below lists the parts that might get contaminated by exhaled gases or bodily fluids during normal use or single fault condition.

| Condition | Parts |
|-------------------------|---|
| With bacteria filter | <ul style="list-style-type: none">• Patient circuit• EtCO₂ airway adapter (if used)• Bacteria filter |
| Without bacteria filter | <ul style="list-style-type: none">• Patient circuit• EtCO₂ airway adapter (if used)• FiO₂ sensor (if used)• Patient air outlet/Pneumatic unit• Air bypass unit/water chamber• Blower/Inlet silencer• Air inlet with filters |

In case of contamination, the internal air pathways of the Nippy 4+ may be disinfected up to 10 times by a maximum 60 minute long validated ozone gas process.

Low resistance bacteria filter, if used, should be replaced every 24 hours.

7.1.3 Patient Circuit



The patient circuit should be cleaned and replaced in accordance with the manufacturer's instructions and care provider's instructions, where applicable. For safety information, read 2.4 *Usage of Patient Circuit*, page 23.

Check the patient circuit regularly for damage. In case of damage, replace the circuit



CAUTION!

Appropriate personnel should determine the duration of use for the patient circuit based on accepted infection control procedures.

7.2 Cleaning and Replacing the Air Filters


The patient air filters are located in the filter cassette at the side of the ventilator.

There are two types of filters:

- Coarse filter (washable, grey)
- Fine filter (disposable, white)


1 Turn off the ventilator and place it on a dust free surface.

2 First Place the filters in the air inlet compartment, with the coarse filter outside the fine filter.



3 Close the side panel carefully for not displacing the filters while closing. For detailed information about closing the side panel, see 3.3.1 *Detaching and Reattaching the Side Panels*, page 35.

4 The cooling air inlet filter is located at the bottom left side of the ventilator.



7.2.1 Coarse Filter (Grey, Washable)

Replace the washable filter patient air filter and the cooling air inlet filter at least once a year. Wash the filters at least once a week.

- 1 Wash the filter using warm water and a mild soap.
- 2 Rinse thoroughly.
- 3 Dry the filter by squeezing it out in a towel. Do not wring the filter.
- 4 Make sure the filter is completely dry before inserting.

7.2.2 Fine Filter (White, Disposable)

Replace the fine filter at least every month, or more frequently when used in high pollution or pollen-rich environments.



CAUTION!

Do not wash or reuse the disposable filter.

7.3 Change of Patients

If the ventilator is used in a clinic by several patients, a low resistance bacterial filter may be used between the air outlet and the patient tube to prevent patient cross-contamination.

- 1 Follow the instructions in 7.1.1 *Main Unit*, page 140, steps 1 to 5.
- 2 Replace the patient filters according to 7.2 *Cleaning and Replacing the Air Filters*, page 141.
- 3 If a low resistance bacterial filter is used, it shall be replaced. To avoid cross-contamination when no bacterial filter has been used, a validated ozone-disinfection process may be used, see the section on disinfecting the main unit internally.
- 4 Use a new patient circuit when the ventilator is used by a new patient.

7.4 Regular Maintenance

Regular maintenance inspections and checks shall be carried out at least every 24 months, according to the ventilator Service Manual

WARNING!



Do not use the device and contact your responsible care provider for an inspection of the device in the event of:

- Unexpected patient symptoms during treatment.
- Unexplainable or sudden pressure, performance or sound changes during operation.
- Suspected damage to the device, including the occurrence of Internal Functional Failure alarms.
- Suspected damage to the click-in battery, including evidence of battery cell leakage.

7.5 Service and Repair

The service and repair of the ventilator must only be carried out by authorised service personnel in accordance with Breas service instructions. Service inspections must always be carried out following any repairs to the device.



Authorised service workshops can order the ventilator Service Manual that contains all technical documentation required for the maintenance and service of the ventilator.

7.6 Storage

Store the ventilator in a dark room, where the temperature range is within -20 to $+60^{\circ}\text{C}$ (-4 to $+140^{\circ}\text{F}$).

For instructions on how to charge the batteries after long time storage, see 5.11 *Using Batteries*, page 70.

CAUTION!



The ventilator must not be stored in a warm place, such as direct sunlight or close to a radiator. The time required for the device to cool from the maximum storage temperature of $+60^{\circ}\text{C}$ ($+140^{\circ}\text{F}$) until it is ready for use in ambient temperature of $+20^{\circ}\text{C}$ ($+68^{\circ}\text{F}$) is 30 minutes.



If stored in a cold environment, let the ventilator adapt to room temperature before using the device. The time required for the device to warm from the minimum storage temperature of -20°C (-4°F) until it is ready for use in ambient temperature of $+20^{\circ}\text{C}$ ($+68^{\circ}\text{F}$) is 30 minutes.

7.7 Disposal

The ventilator, any accessories and all replaced parts must be disposed of and recycled in accordance with the local environmental regulations regarding the disposal of used equipment and waste. Contact your service provider for information regarding the disposal procedure.



NOTE




Batteries used with the ventilator shall be recycled in accordance with the local environmental regulations.

8 Technical Specifications

8.1 System Description

Active Exhalation Valve Configuration


This diagram provides an overview of the ventilator system when used with an active exhalation valve patient circuit.



1. Nippy 4+
2. Tube
3. Active Exhalation valve
4. Patient interface connection
5. Patient


Leakage Port Configuration

This diagram provides an overview of the ventilator system when used with a leakage port patient circuit.



1. Nippy 4+
2. Tube
3. Leakage port / Patient interface connection
4. Patient

8.1.1 Pneumatic Diagram for the ventilator



| Item | Description |
|------|---|
| 1 | Air inlet with filters |
| 2 | Inlet silencer |
| 3 | Blower |
| 4 | Restriction |
| 5 | Flow sensor |
| 6 | Pressure sensors |
| 7 | Air bypass unit/Humidifier |
| 8 | Patient air outlet |
| 9 | Low pressure/bleed-in oxygen connection |

8.2 Data

8.2.1 Worst Case Accuracy

Pressure Control Modes

The worst case Nippy 4+ configuration is the 15 mm patient circuit with HCH humidifier, bacterial filter and EtCO₂ sensor.

Volume control Modes

The worst case Nippy 4+ configuration is the 15mm circuit with or without HCH humidifier, bacterial filter, FiO₂ sensor and EtCO₂ sensor.

8.2.2 Settings Specifications

This section describes the ranges and tolerances for settings that can be made on the Nippy 4+.

All stated tolerances includes measurement uncertainty. The accuracies have been tested with all allowed configurations. Stated tolerances only disclose the maximum tolerance.

Ventilation modes

- Pressure Support (PSV)
- Pressure Support with TgV (PSV+TgV)
TgV= Target Volume
- Pressure Control (PCV)
- Pressure Control with TgV (PCV+TgV)
TgV= Target Volume
- Mouthpiece - Pressure (PCV-MPV)
- SIMV-Pressure (SIMV-P)
SIMV= Synchronized Intermittent Mandatory Ventilation
- Volume Control (VCV)
- Mouthpiece - Volume (VCV-MPV)
- SIMV-Volume (SIMV-V)
SIMV= Synchronized Intermittent Mandatory Ventilation
- CPAP

Device modes

- Clinical
- Home

IPAP

Range/Performance: 4 to 50 cmH₂O.

Tolerance: ± 0.5 cmH₂O or $\pm 5\%$, whichever is greatest.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

EPAP

Range/Performance: 2 cmH₂O to 20 cmH₂O, IPAP -2 cmH₂O or Min IPAP -2 cmH₂O.

Tolerance: ± 0.5 cmH₂O or $\pm 5\%$, whichever is greatest.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

Breath Rate

Range/Performance:

Tolerance: $\pm 2\%$

Resolution: 1 bpm

SIMV Rate

Range/Performance: 4 to 40 bpm (Adult) ,6 to 60 bpm (Paediatric).

Tolerance: $\pm 2\%$

Resolution: 1 bpm

Inspiratory Time

Range/Performance:

Tolerance: $\pm (20 \text{ ms} + 5\% \text{ of setting})$ or $\pm 0.1 \text{ s}$

Resolution: 0.1 s

Backup Inspiratory Time

Range/Performance:

Resolution: 0.1 s

Sigh

Range/Performance rate (frequency): Sigh rate (frequency): Off, every 10 to 250 breaths.

Range/Performance Sigh % (pressure or volume) 25% to 200% of actual set pressure or volume. Limited to 50 cmH₂O or 2500 ml.

Resolution rate: 10 breaths

Resolution %: 25% .

Rise Time

Range/Performance:

1 to 9 (PSV, PCV)

Resolution:

1

Inspiratory Trigger

Range/Performance: 1 to 9 (PSV), 1 to 9, Off (PCV).

Resolution: 1

Support Pressure (in SIMV modes)

Range/Performance: 4 to 60 cmH₂O, 4 to 50 cmH₂O.

Tolerance: ± 0.5 cmH₂O or $\pm 5\%$, whichever is greatest.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

Expiratory Trigger

Range/Performance: 1 to 9.

Resolution: 1

Min Inspiration Time

Range/Performance:

Resolution: 0.1 s

Max Inspiration Time

Range/Performance:

Resolution: 0.1 s

Backup Rate

Range/Performance:

Resolution: 1 bpm

Target Volume

Range/Performance:

Off, 100 to 2000 ml

Tolerance: ± 12 ml or $\pm 10\%$, whichever is greatest.

Resolution:

10 ml below 500 ml

50 ml above 500 ml

Max Pressure

Range/Performance: Min IPAP to 50 cmH₂O.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

Min IPAP

Range/Performance: 4 cmH₂O to Max Pressure.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

Tidal Volume

Range/Performance:

100 to 2000 ml

Tolerance: ± 12 ml or $\pm 10\%$, whichever is greatest.

Resolution: 10 ml below 500 ml, 50 ml above 500 ml

Flow Pattern

Range/Performance: Square, Decelerating

CPAP

Range/Performance:

4 to 20 cmH₂O.

Tolerance: ± 0.5 cmH₂O or $\pm 5\%$, whichever is greatest.

Resolution: 0.5 below 10 cmH₂O, 1.0 above 10 cmH₂O

Audible alarm level

Range/Performance: 1 to 5, where 1 is the lowest volume setting and 5 is the highest volume setting.

Resolution: 1

8.2.3 Monitored Values Specifications

This section describes the ranges and tolerances for monitored values on the Nippy 4+.

All stated tolerances includes measurement uncertainty. The accuracies have been tested with all allowed configurations. Stated tolerances only disclose the maximum tolerance.

P_{peak}

Range/Performance: 4 to 99 cmH₂O.

Resolution: ± 0.5 cmH₂O or $\pm 10\%$, whichever is greatest

EPAP

Range/Performance: 0 to 99 cmH₂O.

Resolution: ± 0.5 cmH₂O or $\pm 10\%$, whichever is greatest

P_{mean}

Range/Performance: 0 to 99 cmH₂O.

Resolution: ± 0.5 cmH₂O or $\pm 10\%$, whichever is greatest

Leakage

Range/Performance: 0 to 99.9 l/min (BTPS*).

Resolution: $\pm 10\%$

V_te

Range/Performance: 0 to 9999 ml (BTPS*).

Resolution: ± 15 ml or 15%, whichever is greatest

FiO₂

Range/Performance: 0 to 100%.

Resolution: $\pm 2\%$

% in TgV

Range/Performance: 0 to 100%.

Resolution: $\pm 1\%$

Total Rate

Range/Performance: 0 to 99 bpm.

Resolution: ± 1 bpm

Spont Rate

Range/Performance: 0 to 99 bpm.

Resolution: ± 1 bpm

% Spont

Range/Performance: 0 to 100%.

SpO₂

Range/Performance: 70 to 100%.

Resolution: ± 3 digits. No motion and flex sensor.

Pulse Rate

Range/Performance: 25 to 240 bpm.

Resolution: ± 3 digits. No motion and flex sensor.

I:E

Range/Performance: 1:10 to 10:1.

Resolution: ± 0.1 unit for E < 9.9, ± 1 unit otherwise.

Insp. Time

Range/Performance: 0.3 to 5 s.

Resolution: ± 0.1 s

Rise Time

Range/Performance: 0.1 to 5 s.

Resolution: $\pm 10\%$ or ± 0.1 s, whichever is greatest

EtCO₂

Range/Performance: 0 to 25%.

Resolution: 0 to 15%: $\pm(0.3 \text{ vol}\% + 4\% \text{ of reading})$. 15 to 25%: unspecified

InspCO₂

Range/Performance: 0 to 25%.

Resolution: 0 to 15%: $\pm(0.3 \text{ vol}\% + 4\% \text{ of reading})$. 15 to 25%: unspecified

8.2.4 Power Supply

AC supply: 100 to 240 V AC, tolerance: $+10\%/-20\%$, 50 to 60 Hz, 1.0 - 2.0 A.

External DC: 19 V DC, tolerance: $19 \text{ V} \pm 6 \text{ V}$. Max 90 W.

Click-in battery: Capacity: 65Wh. Li-ion.

Internal battery: Capacity: 25Wh. Li-ion. Expected service life: 500 full charging cycles.

8.2.5 Environmental Conditions


Operating temperature range: 5 to 40°C (41 to 104°F)

Storage and transport temperature: -20 to +60°C (-4 to +140°F)

Ambient pressure range:

700 to 1100 mbar, corresponding to ~4200 metres (13800 feet) above sea level to ~700 metres (2300 feet) below sea level, at normal atmospheric pressure.

As seen in the graph above, the ventilator is unable to deliver set max pressure at a very low ambient pressure.



Ingress Protection:

IP22

Solid particle protection: Hazardous parts are protected from touch by fingers and by objects greater than 12 mm.

Liquid ingress protection: The protection withstands dripping water less than 15 degrees from vertical.

The ingress protection has been tested by water drips equivalent to 3mm rain/minute for 10 minutes (2.5 minutes for each tilting direction).

8.2.6 Other**Patient Circuit Leakage**

Recommended leakage: 20 to 50 l/min at 10 cmH₂O (leakage circuit)

Minimum leakage: 12 l/min at 4 cmH₂O (leakage circuit)

Oxygen Inlet

Oxygen inlet port: Maximum flow: 15 l/min (medical oxygen). Oxygen coupling is type CPC PMCD181032.

Startup Time

Startup from unpowered state: about 30 seconds.

Sound Power Level

Sound level at 10 cmH₂O in CPAP mode: Less than 30 dB(A). Measured at 1 m.

Alarm sound level : Adjustable 50–80 dB(A), Measured at 1m. Tolerance: ± 5 dB(A).

Miscellaneous

Maximum flow: > 300 l/min

Maximum flow at 20 mbar: > 150 l/min

Maximum limited pressure during single fault condition: 80 cmH₂O (PCV, PSV & VCV) 30 cmH₂O (CPAP)

Breathing resistance under single-fault: <6 cmH₂O at 30 l/min, <6 cmH₂O at 60 l/min

Bias-flow when using active exhalation valve: 8 l/min

Nippy 4+ Dimensions

W × H × D: 216 × 159 × 152 mm

Weight: 2.4 kg

Patient air outlet: 22 mm male, conical standard connector

EtCO₂ Sensor

W × H × D: 8 × 37 × 34 mm

Cable length: 2.4 m

Weight: 75 g

Warm-up time: 10 s

Total system response time: 30 s

Interference from medical gases: O₂: <-0.1% relative CO₂ per % O₂
(calibrated at 21% O₂)

FiO₂ Sensor

Total system response time : 20 s

Filtering/Smoothing Techniques

Pressure: Low pass average time constant 16 ms

Inspiration trigger: Differential mass flow resolution 4 ms

Expiration trigger: Flow low pass filtering with level sensing

SpO₂: No data post-processing done by the ventilator

Effort Belt: Low pass filter: 5Hz, High pass filter: 0.1Hz

8.3 Emission and Immunity Declaration

According to IEC 60601-1-2:2014.

The performance of all functions of the ventilator is considered as essential performance for the purpose of immunity testing.

8.3.1 Nippy 4+ Essential Performance

The ventilator will deliver ventilation at the patient-connection port within its published accuracy specifications and within the alarm limits set by the operator, or generate an alarm condition for high pressure, low pressure, high EPAP, low tidal volume, low minute volume, low breath rate, high and low FiO₂, obstruction, low last power source, or power failure.

The ventilator will provide SpO₂ and pulse rate values within its published accuracy specifications and generate an alarm upon a low SpO₂ condition. The ventilator will

provide indication when the SpO₂ value or pulse rate is potentially incorrect, and generate an alarm condition to indicate when the SpO₂ value update period has exceeded 30 seconds.

The ventilator will provide EtCO₂ and FiO₂ values within its published accuracy specifications and generate an alarm condition upon high and low EtCO₂ and FiO₂ conditions.

Under the immunity test conditions, the following allowances are acceptable:

- Error of delivered volume and EPAP of individual breaths up to 35% and error of the delivered volume and EPAP averaged over a one-minute interval up to 25%.
- Any temporary degradation of SpO₂, EtCO₂ or FiO₂ performance following transient immunity test exposure shall recover from any disruption within 30 seconds.

Additionally, the following shall not be allowed:

- permanent damage or unrecoverable loss of function,
- changes in programmable parameters or settings,
- reset to default settings,
- change of operating mode,
- initiation of unintended operation.

8.3.2 Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The ventilator is intended for use in the electromagnetic environment specified below. The customer or the user of the ventilator should assure that it is used in such an environment.

| Immunity Test | Compliance Level | Electromagnetic Environment - Guidance |
|--|--|---|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±15 kV air | The relative humidity should be at least 5 %. |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | Mains power quality should be that of a typical commercial, hospital and residential environment. |
| Surge IEC 61000-4-5 | ±1 kV line to line | Mains power quality should be that of a typical commercial, hospital and residential environment. |

| Immunity Test | Compliance Level | Electromagnetic Environment - Guidance |
|--|---|---|
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial, hospital and residential environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% UT, 0.5 cycle (multiple phase analysis); 0% UT, 1 cycle; 70% UT, 25/30 cycles (50/60 Hz); 0% UT, 250/300 cycles (50/60 Hz); | Nippy 4+ runs on internal battery during voltage dips, short interruptions and voltage variations on power supply input lines. |




UT is the mains voltage prior to application of the test level.



WARNING!

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ventilator, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

| Immunity Test | Compliance Level | Electromagnetic Environment - Guidance |
|----------------------------|---------------------------------------|---|
| Conducted RF IEC 61000-4-6 | 10 V _{rms} 150 kHz to 80 MHz | $d=0.35*\sqrt{P}$ m at 150 kHz to 80 MHz |
| Radiated RF IEC 61000-4-3 | 20 V/m 80 MHz to 2.5 GHz | $d= 0.6*\sqrt{P}$ m at 80 MHz to 800 MHz $d= 1.2*\sqrt{P}$ m at 800 MHz to 2.5 GHz Equation description: P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with this symbol:  |

NOTE



At 80 MHz and 800 MHz, the higher frequency range applies.



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ventilator is used exceeds the applicable RF compliance level above, the ventilator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ventilator.



b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

8.3.3 Guidance and Manufacturer's Declaration – Electromagnetic Emission

The ventilator are intended for use in the electromagnetic environment specified below. The customer or the user of the ventilator should assure that it is used in such an environment

| Emissions test | Compliance Level | Electromagnetic Environment - Guidance |
|--|------------------|--|
| RF emissions CISPR 11 | Group 1 | The ventilator use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The ventilator are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/flicker emission IEC 61000-3-3 | Complies | |

8.3.4 Recommended separation distances between portable and mobile RF communications equipment and the ventilator

The ventilator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ventilator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ventilator as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to the frequency of transmitter (m) | | |
|---|---|--|---|
| | 150 kHz to 80 MHz $d = 0.35 \cdot \sqrt{P}$ m | 80 MHz to 800 MHz $d = 0.6 \cdot \sqrt{P}$ m | 800 MHz to 2.5 GHz $d = 1.2 \cdot \sqrt{P}$ m |
| 0.01 | 0.035 | 0.06 | 0.12 |
| 0.1 | 0.11 | 0.19 | 0.36 |
| 1 | 0.35 | 0.60 | 1.2 |
| 10 | 1.1 | 1.9 | 3.6 |
| 100 | 3.5 | 6.0 | 12 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE



At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8.3.5 Recommended separation distances between external power conductors and the ventilator

| Rated maximum current in conductor (A) | Separation distance (m) |
|--|---------------------------------|
| | 50-60 Hz $d = I/2\pi H = I/188$ |
| 1 | 0.005 |
| 10 | 0.05 |
| 30 | 0.16 |

For conductors rated at a maximum current not listed above, the recommended separation distance d in metres (m) can be estimated using the equation $d = I/2\pi H$, where I is the maximum current rating of the conductor in amperes (A) according to the transmitter manufacturer; H is the ventilator immunity compliance level to electromagnetic fields in the 50-60 Hz frequency span (30 A/m).

8.3.6 Recommended separation distances between external power conductors and the ventilator

| Rated maximum current in conductor (A) | Separation distance (m) |
|--|---------------------------------|
| | 50-60 Hz $d = I/2\pi H = I/188$ |
| 1 | 0.005 |
| 10 | 0.05 |
| 30 | 0.16 |

For conductors rated at a maximum current not listed above, the recommended separation distance d in metres (m) can be estimated using the equation $d = I/2\pi H$, where I is the maximum current rating of the conductor in amperes (A) according to the transmitter manufacturer; H is the ventilator immunity compliance level to electromagnetic fields in the 50-60 Hz frequency span (30 A/m).

8.4 Delivery Settings

Delivery settings: modes and functions

Ventilation Mode: Pressure Support

Device Mode: Clinical

Profile 1: Active

Profile 2: Off

Profile 3: Off

Delivery Settings, Parameters

IPAP: 15 cmH₂O

EPAP: 5 cmH₂O

Breath Rate: 12 bpm

SIMV Rate: 12 bpm

Inspiration Time: 1.5 s

Rise Time: 3

Inspiratory Trigger: 3

Expiratory Trigger: 3

Maximum Inspiratory Time: Off

Minimum Inspiratory Time: Off

Backup Rate: 12 bpm

Backup Inspiration Time: 1.5 s

Sigh: Off

Sigh Rate: 50 bpm

Sigh %: 125%

Target Volume: Off

Max Pressure: 15 cmH₂O

Min Pressure: 15 cmH₂O

CPAP: 10 cmH₂O

Delivery Settings, Alarms

High Pressure Alarm:

25 cmH₂O

Low Pressure Alarm: 10 cmH₂O

High EPAP Alarm: Off

Low EPAP Alarm: Off

High Flow Alarm: 100 l/min

Low Flow Alarm: 20 l/min

High V_t_i Alarm: 500 ml

High V_t_e Alarm: 500 ml

Low V_t_i Alarm: 300 ml

Low V_t_e Alarm: 300 ml

High MV_i Alarm: Off

High MV_e Alarm: Off

Low MV_i Alarm: Off

Low MV_e Alarm: Off

High Breath Rate Alarm: Off

Low Breath Rate Alarm: Off

Apnoea Alarm: Off

Disconnection Alarm: On

Rebreathing Alarm: On

Obstruction Alarm: Off

High FiO₂ Alarm: Off

Low FiO₂ Alarm: Off

High SpO₂ Alarm: Off

Low SpO₂ Alarm: Off

High EtCO₂ Alarm: 51 mmHg

Low EtCO₂ Alarm: Off

High InspCO₂ Alarm: Off

Low Pulse Rate: Off

High Pulse Rate: Off

Other

Patient operating time: 0 h

Display light: On

Light Intensity: 9

Alarm sound level: 5

CO₂ Unit: mmHg

Auto keypad lock: Off

Pre-use Test:

9 Accessories and Parts

9.1 Breas Accessories List

WARNING!



Only use accessories recommended by Breas Medical. Breas Medical cannot guarantee the performance and safety for the use of other accessories with the ventilator.

NOTE




Accessory equipment connected to the analogue and digital interfaces must be certified according to the respective IEC standards (e.g. IEC 60950 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all configurations must comply with the valid version of the system standard IEC 60601-1-1. Anybody who connects additional equipment to the signal input part or signal output part is configuring a medical system, and is therefore responsible for ensuring the system complies with the requirements of the valid version of the system standard IEC 60601-1-1. If in doubt, consult the technical service department or your local representative.

The following Breas accessories are approved for the Nippy 4+:

Circuit: Single limb 22 mm, disposable

Function: Delivers air to the patient, applied part


Part No: 005060



Circuit: Single limb heated wire 15 mm, disposable

Function: Deliver heated air to the patient, non-invasively


Part No: 006193



Circuit: Single limb with active exhalation valve, disposable

Function: Deliver air to the patient
(applied part)


Part No: 005050



Circuit: Single limb 22 mm with exhalation valve, disposable

Function: Delivers air to the patient
(applied part)


Part No: 007387



Circuit: Single limb 22 mm with exhalation valve and pilot line connector, disposable

Function: Delivers air to the patient
(applied part)


Part No: 007474



Circuit: Single limb 22 mm with leakage port and bacterial filter, disposable

Function: Delivers air to the patient
(applied part)


Part No: 007615



Circuit: Single limb 15 mm

Function: Deliver air to the patient


Part No: 006712



Click-in water chamber

Function: Humidify the patient air


Part No: 006490



Click-in battery

Function: Power source for transportation

Part No: 006265




Click-in Battery Charger

Function: External charger for click-in batteries, available with bank for 2 or 4 batteries)

Part no:

07728 (2 batteries charger)


007730 (4 batteries charger)



Leakage Port

Function: Providing a leakage for clearing exhaled gases.


Part No: 004426



Memory card

Function: Storage and transfer of settings, patient data and usage data


Part No: 006705



Memory card reader/writer

Function: Read/write memory card


Part No: 002185



Remote alarm with cable

Function: Monitor Nippy 4+ alarms remotely


Part No: 10 m: 006348, 25 m: 006349



Remote alarm cable

Function:

Part No: 10 m: 006359, 25 m: 006360, 50 m: 006361



Nurse call cable

Function: Connect the ventilator to a hospital nurse call system


Part No:

NO: 006365

NC: 006364

10 k Ω , NO: 006363


10 k Ω , NC: 006362



Remote start/stop

Function: Start and stop the ventilator remotely. Also, pause audio remotely.


Part No: 006649



Trolley

Function: Mobile use, transportation


Part No: 007384



Universal rail clamp

Function: Attach a humidifier to a trolley.



Part No: 007858



E-cylinder holder

Function: Attach an E-cylinder to a trolley.


Part No: 005128



IV-pole

Function: Pole with hooks to hang IV fluid bags.


Part No: 007859



Mounting bracket

Function: Mount the ventilator to a stand / trolley / rail system.

Part No: 006761




Low resistance bacterial filter (303 Respirgard-II Filter)

Function: Filter air at ventilator outlet

Characteristics

- Resistance: 1.8 cmH₂O @ 60 l/m
- Deadspace: 30 ml
- BFE (Bacterial Filtration Efficiency): 99.9%
- VFE (Viral Filtration Efficiency): 99.8 %




Part No: 004185

FiO₂ sensor

Function: Measure FiO₂ to the patient.


Part No: 006172



SpO₂ module

Function: Connection interface

Part No: 006369




SpO₂ sensor

Function: Finger Clip SpO₂ sensor

Part No:

Adult: 006589


Paediatric: 006590



SpO₂ sensor

Function: Multisite SpO₂ sensor


Part No: 006591



Low pressure oxygen adapter

Function: Oxygen tube adapter with connector for the Nippy 4+.


Part No: 005032



Cable, external DC

Function: External DC cable.


Part No:006709



Cable, external DC to Ventilator Adapter


Function: Connect the ventilator to external DC

Part No: 006710



Cable, Y-adapter, Mains AC and External DC to Ventilator

Function: Connects the ventilator to both mains and external DC at the same time. If the Mains power source is available, it will have precedence over the DC power source.



Part No: 006711

Protective cover

Function: Shock protection


Part No: 006067



Lightweight Mobility Bag

Function: Mobile use


Part No: 007555



EtCO₂ sensor

Function: Measure CO₂ in the airflow


Part No: 006346



Airway adapter

Function: Connects the EtCO₂ sensor to the patient circuit


Part No: 005263 (25 pcs)



PtcCO₂ Cable, Sentec

Function: Connects the ventilator to a Sentec PtcCO₂ monitor.


Part No: 006179



PtcCO₂ Cable, Radiometer

Function: Connects the ventilator to a Radiometer PtcCO₂ monitor.


Part No: 006180



Effort belt communication box

Function: Connects the ventilator to one or two effort belts.


Part No: 006182



Effort belt wireset

Function: Connects an effort belt to the communication box

Part No: 007083



Effort belt

Function: Measures respiratory effort

Part No:


Adult:

24"–74" (107–188 cm): 007085

45"–123" (114–312 cm): 007091

Paediatric:


16"–42" (41–107 cm): 007084



Circuit: Single limb for Mouthpiece ventilation (MPV)

Function: Deliver air to the patient


Part No: 006093



Mouthpiece

Function: Patient interface for Mouthpiece ventilation (MPV)


Part No: 006094



MPV arm

Function: Hold an MPV circuit so Mouthpiece can be mounted close to the patient

Part No: 006095



10 Patient Settings

This section can be copied and used for noting the patient's settings.

Patient Settings - Nippy 4+

Patient

Date

Clinic

Set by

Ventilation mode:.....

| | |
|------------------|-------------------------|
| Patient Circuit | |
| IPAP | Inspiratory Trigger |
| EPAP | Expiratory Trigger |
| Breath Rate | Min Inspiratory Time |
| Inspiratory Time | Max Inspiratory Time |
| Backup Rate | Backup Inspiratory Time |
| Target Volume | Min Pressure |
| Max Pressure | CPAP |
| SIMV Rate | SIMV Support Pressure |

Notes

.....

.....

.....

.....

11 FAA Compliance

To whom it may concern:

The US Department of Transportation (DOT) Final Rule, “Nondiscrimination on the Basis of Disability in Air Travel” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009 provides important requirements for the accommodation of passengers with respiratory assistive devices (Ventilators, Respirators and CPAP machines).

In line with these requirements, respiratory assistive devices may be used onboard an aircraft, without further testing by the carrier, provided they have been tested for Electromagnetic Compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M.

Breas Medical has successfully completed testing for the ventilator System. The ventilator System complies with RTCA/DO-160, Section 21, Category M and can be considered FAA compliant.

Some airlines may require advance notification before travel, and devices may need to be operated by battery. Breas Medical recommends that customers check with their airline.

FAA Compliance (English text)

To whom it may concern:

The US Department of Transportation (DOT) Final Rule, “Nondiscrimination on the Basis of Disability in Air Travel” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009 provides important requirements for the accommodation of passengers with respiratory assistive devices (Ventilators, Respirators and CPAP machines).

In line with these requirements, respiratory assistive devices may be used onboard an aircraft, without further testing by the carrier, provided they have been tested for Electromagnetic Compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M.

Breas Medical has successfully completed testing for the ventilator System. The ventilator System complies with RTCA/DO-160, Section 21, Category M and can be considered FAA compliant.

Some airlines may require advance notification before travel, and devices may need to be operated by battery. Breas Medical recommends that customers check with their airline.

Index

+/- buttons..... 33

A

AC power
connect 43

Accessories
usage 76

Active exhalation valve patient circuit
connect 44

Adjust
patient settings 46

Air inlet, position 34

Air pathway
disinfection 140

Air Temp. Sensor Fail alarm 131

Alarm

Air Temp. Sensor Fail..... 131

Alarm Battery Error..... 128

Ambient Pressure Compensation Lost... 126

Ambient Temperature Compensation
Lost..... 127

Apnoea 112

Check CO₂ Adapter 125

Clock failure..... 133

CO₂ Disconnected (CO₂ Sensor Failure/
Disconnection Alarm)..... 125

CO₂ Sensor Error 126

CO₂ Sensor Error (CO₂ Accuracy Error
Alarm)..... 125

Cooling Fan Error..... 132

Crit. Low Last Power Source..... 122

Database Integrity Fail 132

Disconnection..... 112

Exhalation Valve Control Error 123

FiO₂ Disconnected (FiO₂ Sensor Failure/
Disconnection Alarm)..... 126

Heated Circuit Fault..... 131

Heated Circuit Temp. 129

High breath rate 110

High EPAP 105

High EtCO₂..... 117

High FiO₂ 115

High Flow..... 101

High Humidifier Temp..... 130

High InspCO₂ 118

High MVe 109

High Mvi..... 107

High Patient Air Temp (High Patient Air
Temperature)..... 121

High pressure..... 103

High pulse rate..... 119

High SpO₂ 116

High Vte 108

High Vti..... 106

Humidifier Fault..... 130

Humidifier/bypass loose..... 133

Humidity Compensation Lost..... 127

Internal Error..... 132

Internal Function Failure..... 131

Internal temp high..... 133

Internal/Click-In Battery Hot 129

LED Failure..... 127

Lost Mains Power 123

Low Alarm Battery..... 128

Low Battery 122

Low breath rate..... 111

Low EPAP 105

Low EtCO₂ 118

Low FiO₂..... 116

Low Flow..... 102

Low MVe 109

Low MVi..... 107

Low Patient Air Temp. (Low Patient Air
Temperature Alarm)..... 122

Low pressure 104

Low pulse rate..... 119

Low SpO₂ 117

Low Vte 108

Low Vti..... 106

Obstruction 115

Poor SpO₂ Signal 124

Power Failure Alarm 121

PtcCO₂ Alarm 120

Rebreathing..... 114

SpO₂ Disconnected (SpO₂ Sensor Failure/
Disconnection Alarm)..... 124

SpO₂ Signal Lost..... 124

Alarm Battery Error alarm 128

Alarm sound level
set 66

Alarms 98

technical..... 121

Ambient Pressure Compensation Lost
alarm..... 126

Ambient Temperature Compensation
Lost alarm..... 127

Apnoea alarm..... 112

Artificial nose..... 26

Audible range
operator's position 100

Audio pause and reset 100

B

| | |
|------------------------|----|
| Backup insp. time..... | 56 |
| Backup rate | 54 |
| Bacterial filter..... | 26 |
| Breath rate | 54 |

C

| | |
|--|-----|
| Cables | |
| inspect..... | 45 |
| Calibration | |
| CO ₂ | 79 |
| Caution | |
| Icon | 14 |
| Change | |
| patients | 142 |
| Check | |
| before first use | 41 |
| Check CO ₂ Adapter alarm..... | 125 |
| Circuit | |
| connect | 44 |
| Cleaning..... | 139 |
| CO ₂ sensor..... | 79 |
| main unit..... | 139 |
| patient..... | 140 |
| patient air filters | 141 |
| Cleaning and maintenance | |
| safety information | 28 |
| Click-in humidifier | |
| using | 90 |
| Clock failure alarm | 133 |
| CO ₂ Accuracy Error Alarm (CO ₂ Sensor Error) alarm..... | 125 |
| CO ₂ sensor | |
| clean..... | 79 |
| connect | 79 |
| safety information | 79 |
| specifications..... | 150 |
| usage..... | 79 |
| CO ₂ Sensor Error alarm..... | 126 |
| CO ₂ Sensor Failure/Disconnection | |
| Alarm (CO ₂ Disconnected) alarm | 125 |
| CO ₂ zeroing | 81 |
| Compliance data | |
| show/hide..... | 66 |
| view | 67 |
| Connect | |
| CO ₂ sensor..... | 79 |
| nurse call..... | 77 |
| to mains | 43 |
| Cooling Fan alarm..... | 132 |
| CPAP | 62 |
| Crit. Low Last Power Source alarm | 122 |

D

| | |
|---------------------------|----------|
| Date | |
| set | 66 |
| DC power | |
| click-in battery LED..... | 33 |
| external DC..... | 33 |
| external DC LED..... | 33 |
| internal battery LED..... | 33 |
| Delivery settings..... | 160 |
| Dimensions..... | 146, 150 |
| Disconnection alarm | 112 |
| Disinfection | |
| air pathway..... | 140 |
| Display light | |
| set | 66 |

E

| | |
|---|-----|
| Effort belts | |
| using | 85 |
| Electrical safety | |
| user precautions | 21 |
| Environmental conditions | |
| safety information | 22 |
| specifications..... | 146 |
| EPAP..... | 53 |
| Equipment designation and safety label..... | 38 |
| Exhalation Valve Control Error alarm..... | 123 |
| Exp. Trigger..... | 58 |

F

| | |
|--|-----|
| Filter | |
| bacterial..... | 26 |
| disposable | 31 |
| safety information | 26 |
| Filtering/smoothing techniques | |
| specifications..... | 146 |
| Filters | |
| cleand and replace..... | 141 |
| FiO ₂ Sensor Failure/Disconnection | |
| Alarm (FiO ₂ Disconnected) alarm..... | 126 |
| First use check | 41 |
| Font panel | |
| main unit..... | 33 |
| Function | |
| navigation buttons..... | 33 |

G

| | |
|-------------------------------|----|
| General user precautions..... | 17 |
|-------------------------------|----|

H

| | |
|-----------------------------------|----|
| HCH..... | 26 |
| Heat and moisture exchanger | 26 |

| | | | |
|---------------------------------------|-----|--|-----|
| Heated circuit (on/off. setting)..... | 63 | Internal temp high alarm..... | 133 |
| Heated circuit (temp. setting)..... | 63 | Internal/Click-In Battery Hot alarm..... | 129 |
| Heated Circuit Fault alarm..... | 131 | IPAP..... | 53 |
| Heated Circuit Temp alarm..... | 129 | L | |
| Heated patient circuit | | Label..... | 38 |
| connect..... | 44 | Leakage circuit (single limb) | |
| using..... | 96 | pneumatic diagram..... | 146 |
| High breath rate alarm..... | 110 | LED | |
| High EPAP alarm..... | 105 | click-in battery..... | 33 |
| High EtCO ₂ alarm..... | 117 | external DC..... | 33 |
| High FiO ₂ alarm..... | 115 | front panel..... | 33 |
| High Flow..... | 101 | internal battery..... | 33 |
| High Humidifier Temp. alarm..... | 130 | LED Failure alarm..... | 127 |
| High InspCO ₂ alarm..... | 118 | Lost Mains Power alarm..... | 123 |
| High MVe alarm..... | 109 | Low Alarm Battery alarm..... | 128 |
| High Mvi alarm..... | 107 | Low Battery alarm..... | 122 |
| High Patient Air Temperature Alarm | | Low breath rate alarm..... | 111 |
| (High Patient Air Temp) alarm..... | 121 | Low EPAP alarm..... | 105 |
| High pressure alarm..... | 103 | Low EtCO ₂ alarm..... | 118 |
| High pulse rate alarm..... | 119 | Low FiO ₂ alarm..... | 116 |
| High SpO ₂ alarm..... | 116 | Low Flow..... | 102 |
| High Vte alarm..... | 108 | Low MVe alarm..... | 109 |
| High Vti alarm..... | 106 | Low MV _i alarm..... | 107 |
| HME..... | 26 | Low Patient Air Temperature Alarm | |
| Humidification | | (Low Patient Air Temp) alarm..... | 122 |
| safety information..... | 26 | Low pressure alarm..... | 104 |
| Humidifier | | Low pulse rate..... | 119 |
| detaching..... | 93 | Low SpO ₂ alarm..... | 117 |
| filling..... | 94 | Low Vte alarm..... | 108 |
| installing..... | 91 | Low Vti alarm..... | 106 |
| Humidifier (on/off setting)..... | 62 | M | |
| Humidifier (value setting)..... | 62 | Main components..... | 31 |
| Humidifier Fault alarm..... | 130 | Main unit | |
| Humidifier/bypass loose alarm..... | 133 | clean..... | 139 |
| Humidity Compensation Lost alarm..... | 127 | Mains | |
| Hygroscopic condenser humidifier..... | 26 | connect..... | 43 |
| Hypoventilation..... | 111 | Maintenance | |
| I | | service information..... | 142 |
| Icon | | Max Pressure..... | 61 |
| Caution..... | 14 | Maximum inspiratory trigger..... | 59 |
| Note..... | 14 | Memory card slot | |
| Reference..... | 14 | position..... | 34 |
| Warning..... | 14 | Min Pressure..... | 61 |
| Insp. Time..... | 55 | Minimum inspiratory trigger..... | 59 |
| Insp. trigger..... | 58 | N | |
| Inspect | | Note | |
| before use..... | 45 | Icon..... | 14 |
| cables..... | 45 | Nurse call | |
| placement..... | 45 | connect..... | 77 |
| Internal Error alarm..... | 132 | | |
| Internal Function Failure alarm..... | 131 | | |

| | | | |
|---|-----|--|-----|
| electrical safety | 21 | Regular maintenance | 142 |
| O | | Remote alarm | |
| Obstruction alarm..... | 115 | using | 79 |
| Operating conditions | | Remote Start/Stop | 87 |
| specifications..... | 146 | Repair | 143 |
| Operator's position | | Replace patient air filters | 141 |
| audible range..... | 100 | Reset alarm | 100 |
| Oxygen | | Rise time | 57 |
| safety information | 29 | S | |
| usage..... | 29 | Safety information | |
| Oxygen inlet | | cleaning and maintenance..... | 28 |
| position..... | 38 | CO2 sensor..... | 79 |
| specification | 146 | environmental conditions..... | 22 |
| P | | filters..... | 22 |
| Part numbers | | humidification | 26 |
| main components..... | 31 | oxygen | 29 |
| Pathway | | patient circuit | 23 |
| internal disinfection..... | 140 | Safety label..... | 38 |
| Patient air outlet | | Service | 143 |
| position..... | 34 | Settings | |
| Patient circuit | | at delivery..... | 160 |
| connect | 44 | specification | 146 |
| safety information | 23 | Side panel..... | 34 |
| Patient settings | | Sigh..... | 56 |
| adjust | 46 | SIMV rate | 55 |
| Patients | | SIMV support pressure | 58 |
| change..... | 142 | Single limb circuit | |
| Perform pre-use test..... | 47 | pneumatic diagram..... | 146 |
| Placement | 41 | Sleep mode | 49 |
| inspect..... | 45 | Sound level | |
| Pneumatic diagram | | specification | 146 |
| leakage circuit (single limb)..... | 146 | Specifications | |
| single limb circuit | 146 | CO2 sensor..... | 146 |
| Poor SpO ₂ Signal alarm..... | 124 | dimension | 146 |
| Position | | environmental conditions..... | 146 |
| operator, audible range..... | 100 | Filtering/smoothing techniques | 146 |
| Power Failure Alarm alarm..... | 121 | miscellaneous | 146 |
| Power supply specification | 146 | monitored values..... | 146 |
| Product label..... | 38 | operating conditions..... | 146 |
| Profiles | | oxygen inlet..... | 146 |
| about..... | 52 | power supply..... | 146 |
| delivery settings..... | 160 | settings..... | 146 |
| select..... | 52 | sound level..... | 146 |
| Protective cover | 88 | SpO ₂ | |
| PtcCO ₂ Alarm | 120 | usage..... | 85 |
| R | | SpO ₂ Sensor Failure/Disconnection | |
| Rebreathing alarm | 114 | Alarm (SpO ₂ Disconnected) alarm..... | 124 |
| Reference | | Start-Stop | |
| Icon | 14 | using remote | 87 |
| | | Startup time | |
| | | specification | 146 |

| | |
|---|-----|
| Storage..... | 143 |
| Sup. pressure..... | 58 |
| Switch off..... | 49 |
| Symbols | |
| equipment designation and safety label..... | 38 |
| T | |
| Target volume..... | 59 |
| Technical alarms..... | 121 |
| Time | |
| set..... | 66 |
| Turn off..... | 49 |
| U | |
| Units | |
| set..... | 66 |
| Up/Down buttons..... | 33 |
| User precautions | |
| electrical safety..... | 21 |
| environmental conditions..... | 22 |
| general..... | 17 |
| patient circuit..... | 23 |
| W | |
| Warning | |
| Icon..... | 14 |
| Water chamber | |
| detaching..... | 93 |
| filling..... | 94 |
| installing..... | 91 |
| open lid..... | 95 |
| Water trap..... | 26 |
| Weight..... | 146 |
| Wheelchair | |
| protective cover..... | 88 |